
**English**

## IB10 sphingotest® penKid®

For the Quantitative Determination of  
Proenkephalin in EDTA Whole Blood and Plasma

Nexus Dx, Inc.,  
6759 Mesa Ridge Road  
San Diego, California 92121 USA  
Telephone: 1 (858) 410 4600  
Fax: 1 (858) 410 4700


TheraGenesis GmbH  
Bahnhofstrasse 5  
55276 Oppenheim, Germany  
Tel: +49 (0) 151 506 403 14

022-00075 REV D

### Explanation of Symbols

	CE Mark of Conformity
	Manufacturer
	Catalog Number
	Expiry date/Use by
	Lot number
	<i>In Vitro</i> Diagnostic medical device
	Consult instructions for use
	Store between 2 °C and 8 °C
	European Authorized Representative
	Contains sufficient for <n> test
	Do not reuse
	Serial Number
	Insert disc label side up

# IB10 sphingotest® penKid®

## For *In Vitro* Diagnostic Use

### INTENDED USE

The IB10 sphingotest® penKid® is a rapid point-of-care (POC) immunoassay for the in vitro quantitative determination of Proenkephalin A 119-159 in human EDTA whole blood and plasma. This Test is intended for use in conjunction with the Nexus IB10 Analyzer and provides quantitative results within 20 minutes.

The IB10 sphingotest® penKid® is designed for professional use only and may be used in hospital central laboratories and in alternate care settings such as emergency departments, critical care units, and other sites where near patient testing is practiced.

### SUMMARY AND EXPLANATION OF THE TEST

The IB10 sphingotest® penKid® quantitatively measures the level of Proenkephalin A 119-159 (Proenkephalin) circulating in the bloodstream. Proenkephalin is a stable mid-regional fragment of Proenkephalin A and has been established as a surrogate blood marker for the instable Enkephalin peptides that are derived from the same precursor.<sup>1,2</sup> The blood concentration of Proenkephalin is a reliable surrogate marker for the true glomerular filtration rate (GFR) and is highly dynamic, reflecting worsening and improvement of renal function up to 2 days earlier than standard of care.<sup>2,7</sup> In contrast to most other kidney markers, Proenkephalin is unaffected by inflammation and other comorbidities which makes it suitable for timely monitoring of the actual kidney status in several acute diseases.<sup>3,8-10</sup> It proved to be applicable in a range of populations including critically ill patients, sepsis, heart failure, cardiac surgery, acute myocardial infarction, renal transplantation, and patients with chronic kidney disease.<sup>5,8-15</sup>

The measurement of Proenkephalin levels by IB10 sphingotest® penKid® may hence enable physicians to predict and assess worsening of kidney function in real-time and to improve risk stratification and patient management in acute care settings. Clinicians should use the IB10 sphingotest® penKid® results in conjunction with the patient's other laboratory findings and clinical signs, and interpret the concrete patient values in the context of the patient's clinical situation.

### PROENKEPHALIN

Enkephalins are short-lived pentapeptide hormones that are proteolytically cleaved from Preproenkephalin A and co-released together with stable prohormone fragments, one of which is 4.5 kDa Proenkephalin A 119-159 (also described as Proenkephalin, pro-ENK, and PENK).<sup>1</sup> The mature Enkephalins are endogenous ligands of opioid receptors, and are involved in various physiological processes such as pain signaling, stress response, inflammation and cardiovascular regulation with natriuretic and diuretic effects in the kidneys.<sup>3,16,17</sup> Enkephalins and their receptors are expressed in the central nervous system and peripheral organs, and are highly abundant in the kidneys making direct stimulatory effects on the kidneys very likely.<sup>3,18</sup> Proenkephalin does not interact with plasma proteins and is freely filtrated through the glomerulus, thus blood levels demonstrated to be a functional biomarker of GFR.<sup>3</sup>

### PRINCIPLE

The Nexus IB10 immunochemistry system combines chemistry with microfluidics and centrifugal flow to rapidly prepare a cell free plasma from whole blood that can then be moved through a channel to rehydrate, solubilize and mix with freeze dried immunoconjugates. Using a combination of active flow and capillary action, the test sample is quantitatively measured within 20 minutes with an optical signal level proportional to the analyte(s) concentration.

After addition of the patient sample, the entire test is performed within the Nexus IB10 Analyzer, which provides control of the temperature of the disc, as well as the reaction sequence, centrifugal flow, mixing, incubation time, final signal measurement, quantitation and reporting of results. The Test disc includes a positive internal control to ensure that the Test has operated properly. Each lot is calibrated to provide accurate analyte concentration. Lot specific calibration along with additional information such as the lot expiration date is contained on a QR code label affixed to each disc. It is recommended that external controls be tested at appropriate time intervals to confirm that the system and test lot are performing within acceptable limits.

### REAGENTS

The IB10 sphingotest® penKid® contains all required reagents to evaluate the level of Proenkephalin including dye-conjugated monoclonal anti-Proenkephalin antibody, biotin-conjugated monoclonal anti-Proenkephalin antibodies and streptavidin immobilized at the detection area on the disc.

### MATERIALS PROVIDED

Each box contains the following:

- 10 IB10 sphingotest® penKid® discs, each individually sealed in a foil pouch with a desiccant.
- Instructions for Use (IFU).

## MATERIALS/EQUIPMENT REQUIRED BUT NOT PROVIDED

1. Nexus IB10 Analyzer - Model #BCA-IB10.
2. Commercially available Proenkephalin Controls for external Quality Control (QC). **Contact the Distributor in your area** for recommended external QC materials or related technical assistance.
3. Calibrated reusable fixed or variable volume pipette gun with high precision and accuracy capable of delivering 500 µL of whole blood or plasma.
4. Disposable pipette tips capable of accepting and delivering 500 µL of whole blood or plasma.

## PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Carefully follow the Instructions for Use.
- Prior to testing controls or patient samples, ensure that the Analyzer software is updated to the latest version. (Refer to the Nexus IB10 manual for specific instructions).
- Wear disposable gloves while handling samples.
- Handle samples with care. Samples and used test discs should be treated as potentially infectious and should be discarded as biohazardous material according to local regulations.
- Thoroughly wash hands following handling.
- The result obtained from the IB10 sphingotest® penKid® does not provide a definitive diagnosis and should be interpreted by a physician in conjunction with other laboratory test results and according to current medical guidelines and patient clinical findings.
- Keep the test disc in the sealed pouch until ready for use.
- Do not use the test disc if the pouch is damaged or the seal is broken.
- Do not use the test disc after the expiration date printed on the pouch.
- Prior to use, place the unopened pouch at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
- Always pay attention to cleanliness when handling the test disc. Avoid any contamination from fingerprints or foreign substances. Do not contaminate the sample channel inlet.
- Do not drop or damage the test disc.
- The test disc should be inserted with the label side up, into the Nexus IB10 Analyzer tray immediately after injecting the sample into the disc.
- Do not turn the disc upside down.
- This is a quantitative test; therefore no visual interpretation of the results should be made.

## STORAGE AND STABILITY

- Store the IB10 sphingotest® penKid® disc between 2 and 8 °C (35 and 46 °F) until the expiration date printed on the pouch is reached.
- The IB10 sphingotest® penKid® in its sealed pouch is stable between 18 and 30 °C (64 and 86 °F) for 30 days, provided the expiration date printed on the pouch is not exceeded.

## SAMPLE COLLECTION AND PREPARATION

- The IB10 sphingotest® penKid® is to be run using EDTA whole blood or plasma samples.
- It is recommended that samples be tested as soon as possible after collection.
- Plasma samples should be kept frozen at -20 °C (-4 °F) or lower if longer storage is required.
- Allow samples to equilibrate to room temperature (19 to 25 °C/66 to 77 °F) prior to testing.

## PROCEDURE

### Nexus IB10 Analyzer

 **Consult the Nexus IB10 Analyzer User Manual**

For Analyzer installation, start up and complete instructions for use, refer to the **Nexus IB10 User Manual**. Operator must consult the User Manual prior to use to become familiar with the proper operation and quality control procedures.

## PERFORMING SYSTEM CHECK AND DISC CALIBRATION

Each time the Nexus IB10 Analyzer is turned on, a Self Check is automatically performed. The QR code on each test disc contains information for disc calibration which the Analyzer automatically reads when running a test.

## RUNNING QC WITH EXTERNAL CONTROLS

The manufacturer recommends the use of commercially available Proenkephalin Controls (please refer to section **Materials/Equipment Required But Not Provided**). Ensure that the Proenkephalin controls are handled and prepared according to the corresponding Instructions for Use (IFU).

1. Remove an unopened test pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes prior to testing.
2. Open the pouch and remove the test disc.
3. Place the test disc on a level surface.
4. On the Nexus IB10 Analyzer press **New Analysis**.
5. The Analyzer will perform a general system check.
6. Enter the external control ID manually (up to 20 characters can be used for the ID) or enter the external control ID by using the barcode scanner.
7. Mix the external quality control vial by gently inverting the vial several times before sampling.
8. **Testing External Control Samples on the IB10 sphingotest® penKid®**
  - Using a precision pipette (fixed or adjusted to 500 µL), slowly draw well-mixed external quality control sample into the pipette tip.
  - Positioning the tapered pipette tip at a 45° angle, pierce the **X** on the red dot to expose the sample channel inlet.
  - Slowly express the external control sample into the inlet applying minimum, but continuous force, on the pipette plunger.
  - Express sample to the first stop on the pipette, at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
  - Press **QC** on the Nexus IB10 Analyzer display.
  - When the tray opens, insert the filled Test disc into the tray and press **Run**.
  - The tray will close and perform a disc validity check.
  - A screen appears for selecting quality control materials (please refer to the **Quality Control Settings** section of the **Nexus IB10 Analyzer User Manual** for how to update quality control material [external controls]).
  - Select the quality control material to be tested.
  - Press **OK** on the Nexus IB10 Analyzer display.
  - Within 20 minutes, the Nexus IB10 Analyzer will display the results on the screen.
  - Results will print automatically (if selected during Set Up) or press **Print**.
  - When the test is complete, analyze and compare the result with the Expected Value reported in the external control IFU for the external control level as measured using the IB10 sphingotest® penKid®.
  - Remove the test disc and discard in appropriate receptacle.
  - If the external control result is outside the expected range, refer to the Quality Control section below.

*Note: If the test run is cancelled before a test result is displayed, the test disc cannot be reused and should be disposed of appropriately.*

## TESTING PATIENT SAMPLES ON THE NEXUS IB10 ANALYZER

1. Remove an unopened Test pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
2. Open the pouch and remove the Test disc.
3. Place the Test disc on a level surface.
4. On the Nexus IB10 Analyzer press **New Analysis**.
5. The Analyzer will perform a general system check.
6. Enter the Patient ID manually (up to 20 characters can be used for the ID) or enter the Patient ID by using the barcode scanner.
7. Mix the whole blood patient sample tube by gently inverting the tube several times before testing.

## 8. Testing Patient Samples on the IB10 sphingotest® penKid®

- Using a precision pipette (fixed or adjusted to 500 µL), slowly draw well-mixed patient sample into the pipette tip.
- Positioning the tapered pipette tip at a 45° angle, pierce the **X** on the red dot to expose the sample channel inlet.
- Slowly express the patient sample into the inlet applying minimum, but continuous force, on the pipette plunger.
- Express sample to the first stop on the pipette, at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
- Press **OK** on the Nexus IB10 Analyzer display.
- When the tray opens, insert the filled Test disc into the tray and press **Run**.
- Within 20 minutes, the Nexus IB10 Analyzer will display the results on the screen.
- Results will print automatically (if selected during Set Up) or press **Print**.
- Remove the test disc and discard in appropriate receptacle.

*Note: If the test run is cancelled before a test result is displayed, the test disc cannot be reused and should be disposed of appropriately.*

## INTERPRETATION OF RESULTS

The range of Proenkephalin concentrations reported by the Nexus IB10 is 50 pmol/L to 500 pmol/L. Results below or above this range will be shown as “<50 pmol/L” or “>500 pmol/L”, respectively.

## Quality Control

### EXTERNAL CONTROLS

Good laboratory practice includes the use of external controls to ensure proper test performance. It is recommended that prior to using a new lot or shipment of IB10 sphingotest® penKid®, the performance of the lot should be confirmed by testing with external controls (see section **Materials/Equipment Required But Not Provided**) to ensure the test will deliver the correct result. The frequency of quality control testing should be determined according to individual laboratory standard quality control procedures. Upon confirmation of the expected results, the test discs are ready for use with patient samples. Controls should also be used any time the validity of test results is questionable. If external controls do not perform as expected, do not use the IB10 sphingotest® penKid® and contact the Distributor in your area for Technical Assistance.

### INTERNAL CONTROL

The IB10 sphingotest® penKid® has a built in positive procedural control. The Nexus IB10 Analyzer automatically determines the presence of this control thereby confirming that the test run has delivered a valid result. If the control does not form or if it is not recognized by the Analyzer, the test result is considered “invalid” and the test must be repeated.

### LIMITATIONS

The test result should be evaluated and interpreted in conjunction with the patient’s medical history, symptoms and other clinical information.

As with any assay employing mouse antibodies, the possibility exists for interference by heterophilic antibodies in the sample.<sup>19,20</sup> Specimens from individuals who have been regularly exposed to animals or who have been treated with animal products may contain these antibodies.

Other substances and/or factors not listed, e.g. technical or procedural error, may interfere with the test and cause inaccurate results.

*Nexus Dx, Inc. offers products for their intended use. Refer to the specific product literature for the intended use statements for each product. Product claims are subject to change. Nexus Dx, Inc.’s expressed and implied warranties (inclusive of implied warranties of merchantability and fitness) are conditional upon adherence to, or observance of Nexus Dx, Inc. published directions with regard to the use of Nexus Dx, Inc. products. Under no circumstances will Nexus Dx, Inc. be liable for any indirect or consequential damages.*

For Technical Assistance please contact your local Distributor.

## Performance Characteristics

### MEASURING RANGE

The IB10 sphingotest® penKid® has been demonstrated to provide measurable results at Proenkephalin levels from 50 pmol/L to 500 pmol/L.

### ANALYTICAL SENSITIVITY

The LoD (Limit of Detection) of the IB10 sphingotest® penKid® is 50.0 pmol/L, determined according to Clinical and Laboratory Standards Institute (CLSI) Approved Guideline EP17-A2.<sup>21</sup> The proportion of false positives ( $\alpha$ ) and false negatives ( $\beta$ ) is less than 5% (limit of blank = 21.5 pmol/L).

LoQ (Limit of Quantitation) is the lowest Proenkephalin concentration that can be reproducibly measured with a total coefficient of variation of at most 20%. It was determined to be 50.0 pmol/L.

## INTERFERING SUBSTANCES

The following substances were tested for potential interference with Proenkephalin measurement in the IB10 sphingotest® penKid® based on CLSI Approved Guideline EP7-A2 and EP37 (Table 1).<sup>22-23</sup> No significant interference was observed at the maximum concentration for substances listed in the table below.

Table 1.

Potential Interferents	Maximum Concentration
Bilirubin (Conjugated)	40 mg/dL
Bilirubin (Free)	60 mg/dL
Biotin	600 mg/mL
Cholesterol	400 mg/dL
Citrate	16.3 mmol/L
Creatinine	15 mg/dL
Glucose	250 mg/dL
Hemoglobin	5 mg/dL
Unfractionated Heparin	300 U/dL
Fractionated Heparin	3000 U/L
Human anti-mouse antibody (HAMA)	645 ng/mL
Rheumatoid factor (RF)	1500 IU/mL
Total Protein	3000 mg/dL
Triglycerides	1500 mg/dL
Urea	120 mg/dL

## HOOK EFFECT

No high dose hook effect is observed for Proenkephalin concentrations up to 250,000 pmol/L.

## PRECISION

The precision of Proenkephalin measurement of the IB10 sphingotest® penKid® was determined using samples where Proenkephalin was added to normal human plasma at three concentrations (Table 2). The within-run and total precision was performed in two runs per day, in replicates of 4 per run at each concentration level for a 20 day period for a total number of repetitions of 160 at each concentration level. The within-run and total precision were computed according to the CLSI Approved Guideline EP05-A3.<sup>24</sup>

Table 2.

Sample	Test lot	Mean (pmol/l)	Within-run Precision		Total Precision	
			Std. dev. (pmol/L)	CV (%)	Std. dev. (pmol/L)	CV (%)
1	A	82.9	8.4	10.1%	9.4	11.4%
	B	80.1	8.0	10.0%	8.5	10.6%
2	A	203.0	19.0	9.4%	20.7	10.2%
	B	204.5	24.5	12.0%	25.0	12.2%

## WHOLE BLOOD Vs. PLASMA CORRELATION

A comparison study was performed using matching EDTA whole blood and EDTA plasma samples. When performing a Passing-Bablok regression analysis comparing the whole blood concentrations versus the corresponding plasma concentrations from the same subject samples (N=30), a slope of 1.06 (95% C.I. = [0.83-1.17]), an intercept of -9.3 pmol/L, and a correlation coefficient of  $r = 0.96$  were obtained.

## METHOD COMPARISON

An equivalence study was performed between the IB10 sphingotest® penKid® using the Nexus IB10 Analyzer, and the sphingotest® penKid® immuoluminometric (ILMA) assay. A total of 197 samples were tested within a concentration range of Proenkephalin of 50.0 pmol/L to 497.9 pmol/L. The Passing-Bablok regression was:

**IB10 sphingotest® penKid® = 0.93 x (sphingotest® penKid® ILMA) – 0.3 pmol/L**

**Correlation coefficient,  $r = 0.92$**

## REFERENCES

1. Ernst A, Köhrle J, Bergmann A. Proenkephalin A 119-159, a stable proenkephalin A precursor fragment identified in human circulation. *Peptides*. 2006;27: 1835-40.
2. Donato LJ, Meeusen JW, Lieske JC, Bergmann D, Sparwasser A, Jaffe AS. Analytical performance of an immunoassay to measure proenkephalin. *Clin Biochem*. 2018;58: 72-77.
3. Beunders R, Struck J, Wu AH, Zarbock A, Di Somma S, Mehta RL, Koyner JL, Nadim MK, Maisel AS, Murray PT, Neath SX, Jaffe AS, Pickkers P. Proenkephalin (PENK) as a novel biomarker for kidney function. *J Appl Lab Med*. 2017; DOI: 10.1373/jalm.2017.023598 400-12.
4. Matsue Y, Ter Maaten JM, Struck J, Metra M, O'Connor CM, Ponikowski P, Teerlink JR, Cotter G, Davison B, Cleland JG, Givertz MM, Bloomfield DM, Dittrich HC, van Veldhuisen DJ, van der Meer P, Damman K, Voors AA. Clinical correlates and prognostic value of proenkephalin in acute and chronic heart failure. *J Card Fail*. 2017;23: 231-39.
5. Kieneker LM, Hartmann O, Struck J, Bergmann A, Gansevoort RT, Joosten MM, van den Berg E, de Boer RA, Bakker SJL. Plasma proenkephalin and poor long-term outcome in renal transplant recipients. *Transplant Direct*. 2017;3: e190.
6. Leitje G, Beunders R, van Groenendal R, Peters E, Kox M, Pickkers P. Proenkephalin, the new marker for kidney function on the Intensive care unit? *Intensive Care Med Exp*. 2017;5 (Suppl 2): 51-52.
7. Beunders R, Meekes M, Struck J, Pickkers P. Proenkephalin, a novel biomarker for kidney function, is earlier in detecting acute kidney injury compared to creatinine. *Intensive Care Med Exp*. 2019;7 (Suppl 3): 000960.
8. Marino R, Struck J, Hartmann O, Maisel AS, Rehfeldt M, Magrini L, Melander O, Bergmann A, Di Somma S. Diagnostic and short-term prognostic utility of plasma pro-enkephalin (pro-ENK) for acute kidney injury in patients admitted with sepsis in the emergency department. *J Nephrol*. 2015;28: 717-24.
9. Hollinger A, Wittebole X, Francois B, Pickkers P, Antonelli M, Gayat E, Chousterman BG, Lascarrou JB, Dugernier T, Di Somma S, Struck J, Bergmann A, Beishuizen A, Constatin JM, Damoiseil C, Deye N, Gaudry S, Huberlant V, Marx G, Mercier E, Oueslati H, Hartmann O, Sonnevile R, Laterre PF, Mebazaa A, Legrand M. Proenkephalin A 119-159 (Penkid) is an early biomarker of septic acute kidney injury: The Kidney in Sepsis and Septic Shock (Kid-SSS) study. *Kidney Int Rep*. 2018;3: 1424-33.
10. Gayat E, Cariou A, Deye N, Vieillard-Baron A, Jaber S, Damoiseil C, Lu Q, Monnet X, Rennuit I, Azoulay E, Leone M, Oueslati H, Guidet B, Friedman D, Tesniere A, Sonnevile R, Montravers P, Pili-Floury S, Lefrant JY, Duranteau J, Laterre PF, Brechot N, Chevrel K, Michel M, Cholley B, Legrand M, Launay JM, Vicaut E, Singer M, Resche-Rigon M, Mebazaa A. Determinants of long-term outcome in ICU survivors: results from the FROG-ICU study. *Crit Care*. 2018;22: 8.
11. Ng LL, Sandhu JK, Narayan H, Quinn PA, Squire IB, Davies JE, Bergmann A, Maisel A, Jones DJ. Proenkephalin and prognosis after acute myocardial infarction. *J Am Coll Cardiol*. 2014;63: 280-9.
12. Shah KS, Taub P, Patel M, Rehfeldt M, Struck J, Clopton P, Mehta RL, Maisel AS. Proenkephalin predicts acute kidney injury in cardiac surgery patients. *Clin Nephrol*. 2015;83: 29-35.
13. Ng LL, Squire IB, Jones DJL, Cao TH, Chan DCS, Sandhu JK, Quinn PA, Davies JE, Struck J, Hartmann O, Bergmann A, Mebazaa A, Gayat E, Arrigo M, Akiyama E, Sabti Z, Lohrmann J, Twerenbold R, Herrmann T, Schumacher C, Kozhuharov N, Mueller C, and Great Network. Proenkephalin, renal dysfunction, and prognosis in patients with acute heart failure: A GREAT Network study. *J Am Coll Cardiol*. 2017;69: 56-69.
14. Breidthardt T, Jaeger C, Christ A, Klima T, Mosimann T, Twerenbold R, Boeddinghaus J, Nestelberger T, Badertscher P, Struck J, Bergmann A, Hartmann O, Kalbermatter S, Marenzi G, Mueller C. Proenkephalin for the early detection of acute kidney injury in hospitalized patients

with chronic kidney disease. *Eur J Clin Invest*. 2018: e12999.

15. Emmens JE, Ter Maaten JM, Damman K, van Veldhuisen DJ, de Boer RA, Struck J, Bergmann A, Sama IE, Streng KW, Anker SD, Dickstein K, Lang CC, Metra M, Samani NJ, Ng LL, Voors AA. Proenkephalin, an opioid system surrogate, as a novel comprehensive renal marker in heart failure. *Circ Heart Fail*. 2019;12: e005544.
16. Holaday, JW. Cardiovascular effects of endogenous opiate systems. *Annu Rev Pharmacol Toxicol*. 1983;23: 541-94.
17. Sezen, SF, Kenigs VA, Kapusta DR. Renal excretory responses produced by the delta opioid agonist, BW373U86, in conscious rats. *J Pharmacol Exp Ther*. 1998;287: 238-45.
18. Denning, GM, Ackermann LW, Barna TJ, Armstrong JG, Stoll LL, Weintraub NL, Dickson EW. Proenkephalin expression and enkephalin release are widely observed in non-neuronal tissues. *Peptides*. 2008;29: 83-92.
19. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45: 942-6.
20. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15: 108-15.
21. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
22. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
23. McEnroe RJ, Dimeski G, Miller WG, et al. Supplemental Tables for Interference Testing in Clinical Chemistry. CLSI Document EP37, First Edition, 2018.
24. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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- DM, Dittrich HC, van Veldhuisen DJ, van der Meer P, Damman K, Voors AA. Clinical correlates and prognostic value of proenkephalin in acute and chronic heart failure. *J Card Fail.* 2017;23: 231-39.
5. Kieneker LM, Hartmann O, Struck J, Bergmann A, Gansevoort RT, Joosten MM, van den Berg E, de Boer RA, Bakker SJL. Plasma proenkephalin and poor long-term outcome in renal transplant recipients. *Transplant Direct.* 2017;3: e190.
  6. Leijte G, Beunders R, van Groenendal R, Peters E, Kox M, Pickkers P. Proenkephalin, the new marker for kidney function on the Intensive care unit? *Intensive Care Med Exp.* 2017;5 (Suppl 2): 51-52.
  7. Beunders R, Meekes M, Struck J, Pickkers P. Proenkephalin, a novel biomarker for kidney function, is earlier in detecting acute kidney injury compared to creatinine. *Intensive Care Med Exp.* 2019;7 (Suppl 3): 000960.
  8. Marino R, Struck J, Hartmann O, Maisel AS, Rehfeldt M, Magrini L, Melander O, Bergmann A, Di Somma S. Diagnostic and short-term prognostic utility of plasma pro-enkephalin (pro-ENK) for acute kidney injury in patients admitted with sepsis in the emergency department. *J Nephrol.* 2015;28: 717-24.
  9. Hollinger A, Wittebole X, Francois B, Pickkers P, Antonelli M, Gayat E, Chousterman BG, Lascarrou JB, Dugernier T, Di Somma S, Struck J, Bergmann A, Beishuizen A, Constatin JM, Damoiseil C, Deye N, Gaudry S, Huberlant V, Marx G, Mercier E, Oueslati H, Hartmann O, Sonnevile R, Laterre PF, Mebazaa A, Legrand M. Proenkephalin A 119-159 (Penkid) is an early biomarker of septic acute kidney injury: The Kidney in Sepsis and Septic Shock (Kid-SSS) study. *Kidney Int Rep.* 2018;3: 1424-33.
  10. Gayat E, Cariou A, Deye N, Vieillard-Baron A, Jaber S, Damoiseil C, Lu Q, Monnet X, Rennuit I, Azoulay E, Leone M, Oueslati H, Guidet B, Friedman D, Tesniere A, Sonnevile R, Montravers P, Pili-Floury S, Lefrant JY, Duranteau J, Laterre PF, Brechot N, Chevreur K, Michel M, Cholley B, Legrand M, Launay JM, Vicaut E, Singer M, Resche-Rigon M, Mebazaa A. Determinants of long-term outcome in ICU survivors: results from the FROG-ICU study. *Crit Care.* 2018;22: 8.
  11. Ng LL, Sandhu JK, Narayan H, Quinn PA, Squire IB, Davies JE, Bergmann A, Maisel A, Jones DJ. Proenkephalin and prognosis after acute myocardial infarction. *J Am Coll Cardiol.* 2014;63: 280-9.
  12. Shah KS, Taub P, Patel M, Rehfeldt M, Struck J, Clopton P, Mehta RL, Maisel AS. Proenkephalin predicts acute kidney injury in cardiac surgery patients. *Clin Nephrol.* 2015;83: 29-35.
  13. Ng LL, Squire IB, Jones DJL, Cao TH, Chan DCS, Sandhu JK, Quinn PA, Davies JE, Struck J, Hartmann O, Bergmann A, Mebazaa A, Gayat E, Arrigo M, Akiyama E, Sabti Z, Lohrmann J, Twerenbold R, Herrmann T, Schumacher C, Kozhuharov N, Mueller C, and Great Network. Proenkephalin, renal dysfunction, and prognosis in patients with acute heart failure: A GREAT Network study. *J Am Coll Cardiol.* 2017;69: 56-69.
  14. Breidhardt T, Jaeger C, Christ A, Klima T, Mosimann T, Twerenbold R, Boeddinghaus J, Nestelberger T, Badertscher P, Struck J, Bergmann A, Hartmann O, Kalbermatter S, Marenzi G, Mueller C. Proenkephalin for the early detection of acute kidney injury in hospitalized patients with chronic kidney disease. *Eur J Clin Invest.* 2018: e12999.
  15. Emmens JE, Ter Maaten JM, Damman K, van Veldhuisen DJ, de Boer RA, Struck J, Bergmann A, Sama IE, Streng KW, Anker SD, Dickstein K, Lang CC, Metra M, Samani NJ, Ng LL, Voors AA. Proenkephalin, an opioid system surrogate, as a novel comprehensive renal marker in heart failure. *Circ Heart Fail.* 2019;12: e005544.
  16. Holaday, JW. Cardiovascular effects of endogenous opiate systems. *Annu Rev Pharmacol Toxicol.* 1983;23: 541-94.
  17. Sezen, SF, Kenigs VA, Kapusta DR. Renal excretory responses produced by the delta opioid agonist, BW373U86, in conscious rats. *J Pharmacol Exp Ther.* 1998;287: 238-45.
  18. Denning, GM, Ackermann LW, Barna TJ, Armstrong JG, Stoll LL, Weintraub NL, Dickson EW. Proenkephalin expression and enkephalin release are widely observed in non-neuronal tissues. *Peptides.* 2008;29: 83-92.
  19. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45: 942-6.
  20. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15: 108-15.
  21. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
  22. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
  23. McEnroe RJ, Dimeski G, Miller WG, et al. Supplemental Tables for Interference Testing in Clinical Chemistry. CLSI Document EP37, First Edition, 2018.
  24. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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proenkephalin and poor long-term outcome in renal transplant recipients. *Transplant Direct*. 2017;3: e190.

6. Leitte G, Beunders R, van Groenendal R, Peters E, Kox M, Pickkers P. Proenkephalin, the new marker for kidney function on the Intensive care unit? *Intensive Care Med Exp*. 2017;5 (Suppl 2): 51-52.
7. Beunders R, Meekes M, Struck J, Pickkers P. Proenkephalin, a novel biomarker for kidney function, is earlier in detecting acute kidney injury compared to creatinine. *Intensive Care Med Exp*. 2019;7 (Suppl 3): 000960.
8. Marino R, Struck J, Hartmann O, Maisel AS, Rehfeldt M, Magrini L, Melander O, Bergmann A, Di Somma S. Diagnostic and short-term prognostic utility of plasma pro-enkephalin (pro-ENK) for acute kidney injury in patients admitted with sepsis in the emergency department. *J Nephrol*. 2015;28: 717-24.
9. Hollinger A, Wittebole X, Francois B, Pickkers P, Antonelli M, Gayat E, Chousterman BG, Lascarrou JB, Dugernier T, Di Somma S, Struck J, Bergmann A, Beishuizen A, Constatin JM, Damoiseil C, Deye N, Gaudry S, Huberlant V, Marx G, Mercier E, Oueslati H, Hartmann O, Sonnevile R, Laterre PF, Mebazaa A, Legrand M. Proenkephalin A 119-159 (Penkid) is an early biomarker of septic acute kidney injury: The Kidney in Sepsis and Septic Shock (Kid-SSS) study. *Kidney Int Rep*. 2018;3: 1424-33.
10. Gayat E, Cariou A, Deye N, Vieillard-Baron A, Jaber S, Damoiseil C, Lu Q, Monnet X, Rennuit I, Azoulay E, Leone M, Oueslati H, Guidet B, Friedman D, Tesniere A, Sonnevile R, Montravers P, Pili-Floury S, Lefrant JY, Duranteau J, Laterre PF, Brechot N, Chevreur K, Michel M, Chollet B, Legrand M, Launay JM, Vicaut E, Singer M, Resche-Rigon M, Mebazaa A. Determinants of long-term outcome in ICU survivors: results from the FROG-ICU study. *Crit Care*. 2018;22: 8.
11. Ng LL, Sandhu JK, Narayan H, Quinn PA, Squire IB, Davies JE, Bergmann A, Maisel A, Jones DJ. Proenkephalin and prognosis after acute myocardial infarction. *J Am Coll Cardiol*. 2014;63: 280-9.
12. Shah KS, Taub P, Patel M, Rehfeldt M, Struck J, Clopton P, Mehta RL, Maisel AS. Proenkephalin predicts acute kidney injury in cardiac surgery patients. *Clin Nephrol*. 2015;83: 29-35.
13. Ng LL, Squire IB, Jones DJL, Cao TH, Chan DCS, Sandhu JK, Quinn PA, Davies JE, Struck J, Hartmann O, Bergmann A, Mebazaa A, Gayat E, Arrigo M, Akiyama E, Sabti Z, Lohrmann J, Twerenbold R, Herrmann T, Schumacher C, Kozhuharov N, Mueller C, and Great Network. Proenkephalin, renal dysfunction, and prognosis in patients with acute heart failure: A GREAT Network study. *J Am Coll Cardiol*. 2017;69: 56-69.
14. Breidhardt T, Jaeger C, Christ A, Klima T, Mosimann T, Twerenbold R, Boeddinghaus J, Nestelberger T, Badertscher P, Struck J, Bergmann A, Hartmann O, Kalbermatter S, Marenzi G, Mueller C. Proenkephalin for the early detection of acute kidney injury in hospitalized patients with chronic kidney disease. *Eur J Clin Invest*. 2018: e12999.
15. Emmens JE, Ter Maaten JM, Damman K, van Veldhuisen DJ, de Boer RA, Struck J, Bergmann A, Sama IE, Streng KW, Anker SD, Dickstein K, Lang CC, Metra M, Samani NJ, Ng LL, Voors AA. Proenkephalin, an opioid system surrogate, as a novel comprehensive renal marker in heart failure. *Circ Heart Fail*. 2019;12: e005544.
16. Holaday, JW. Cardiovascular effects of endogenous opiate systems. *Annu Rev Pharmacol Toxicol*. 1983;23: 541-94.
17. Sezen, SF, Kenigs VA, Kapusta DR. Renal excretory responses produced by the delta opioid agonist, BW373U86, in conscious rats. *J Pharmacol Exp Ther*. 1998;287: 238-45.
18. Denning, GM, Ackermann LW, Barna TJ, Armstrong JG, Stoll LL, Weintraub NL, Dickson EW. Proenkephalin expression and enkephalin release are widely observed in non-neuronal tissues. *Peptides*. 2008;29: 83-92.
19. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45: 942-6.
20. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15: 108-15.
21. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
22. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
23. McEnroe RJ, Dimeski G, Miller WG, et al. Supplemental Tables for Interference Testing in Clinical Chemistry. CLSI Document EP37, First Edition, 2018.
24. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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unit? *Intensive Care Med Exp.* 2017;5 (Suppl 2): 51-52.

7. Beunders R, Meekes M, Struck J, Pickkers P. Proenkephalin, a novel biomarker for kidney function, is earlier in detecting acute kidney injury compared to creatinine. *Intensive Care Med Exp.* 2019;7 (Suppl 3): 000960.
8. Marino R, Struck J, Hartmann O, Maisel AS, Rehfeldt M, Magrini L, Melander O, Bergmann A, Di Somma S. Diagnostic and short-term prognostic utility of plasma pro-enkephalin (pro-ENK) for acute kidney injury in patients admitted with sepsis in the emergency department. *J Nephrol.* 2015;28: 717-24.
9. Hollinger A, Wittebole X, Francois B, Pickkers P, Antonelli M, Gayat E, Chousterman BG, Lascarrou JB, Dugernier T, Di Somma S, Struck J, Bergmann A, Beishuizen A, Constatin JM, Damoisel C, Deye N, Gaudry S, Huberlant V, Marx G, Mercier E, Oueslati H, Hartmann O, Sonnevile R, Laterre PF, Mebazaa A, Legrand M. Proenkephalin A 119-159 (Penkid) is an early biomarker of septic acute kidney injury: The Kidney in Sepsis and Septic Shock (Kid-SSS) study. *Kidney Int Rep.* 2018;3: 1424-33.
10. Gayat E, Cariou A, Deye N, Vieillard-Baron A, Jaber S, Damoisel C, Lu Q, Monnet X, Rennuit I, Azoulay E, Leone M, Oueslati H, Guidet B, Friedman D, Tesniere A, Sonnevile R, Montravers P, Pili-Floury S, Lefrant JY, Duranteau J, Laterre PF, Brechot N, Chevrel K, Michel M, Cholley B, Legrand M, Launay JM, Vicaut E, Singer M, Resche-Rigon M, Mebazaa A. Determinants of long-term outcome in ICU survivors: results from the FROG-ICU study. *Crit Care.* 2018;22: 8.
11. Ng LL, Sandhu JK, Narayan H, Quinn PA, Squire IB, Davies JE, Bergmann A, Maisel A, Jones DJ. Proenkephalin and prognosis after acute myocardial infarction. *J Am Coll Cardiol.* 2014;63: 280-9.
12. Shah KS, Taub P, Patel M, Rehfeldt M, Struck J, Clopton P, Mehta RL, Maisel AS. Proenkephalin predicts acute kidney injury in cardiac surgery patients. *Clin Nephrol.* 2015;83: 29-35.
13. Ng LL, Squire IB, Jones DJL, Cao TH, Chan DCS, Sandhu JK, Quinn PA, Davies JE, Struck J, Hartmann O, Bergmann A, Mebazaa A, Gayat E, Arrigo M, Akiyama E, Sabti Z, Lohrmann J, Twerenbold R, Herrmann T, Schumacher C, Kozhuharov N, Mueller C, and Great Network. Proenkephalin, renal dysfunction, and prognosis in patients with acute heart failure: A GREAT Network study. *J Am Coll Cardiol.* 2017;69: 56-69.
14. Breidthardt T, Jaeger C, Christ A, Klima T, Mosimann T, Twerenbold R, Boeddinghaus J, Nestelberger T, Badertscher P, Struck J, Bergmann A, Hartmann O, Kalbermatter S, Marenzi G, Mueller C. Proenkephalin for the early detection of acute kidney injury in hospitalized patients with chronic kidney disease. *Eur J Clin Invest.* 2018: e12999.
15. Emmens JE, Ter Maaten JM, Damman K, van Veldhuisen DJ, de Boer RA, Struck J, Bergmann A, Sama IE, Streng KW, Anker SD, Dickstein K, Lang CC, Metra M, Samani NJ, Ng LL, Voors AA. Proenkephalin, an opioid system surrogate, as a novel comprehensive renal marker in heart failure. *Circ Heart Fail.* 2019;12: e005544.
16. Holaday, JW. Cardiovascular effects of endogenous opiate systems. *Annu Rev Pharmacol Toxicol.* 1983;23: 541-94.
17. Sezen, SF, Kenigs VA, Kapusta DR. Renal excretory responses produced by the delta opioid agonist, BW373U86, in conscious rats. *J Pharmacol Exp Ther.* 1998;287: 238-45.
18. Denning, GM, Ackermann LW, Barna TJ, Armstrong JG, Stoll LL, Weintraub NL, Dickson EW. Proenkephalin expression and enkephalin release are widely observed in non-neuronal tissues. *Peptides.* 2008;29: 83-92.
19. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45: 942-6.
20. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15: 108-15.
21. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
22. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
23. McEnroe RJ, Dimeski G, Miller WG, et al. Supplemental Tables for Interference Testing in Clinical Chemistry. CLSI Document EP37, First Edition, 2018.
24. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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**IB10 sphingotest® penKid® = 0,93 x (sphingotest® penKid® ILMA) – 0,3 pmol/L**

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## REFERÊNCIAS

- Ernst A, Köhrlé J, Bergmann A. Proenkephalin A 119-159, a stable proenkephalin A precursor fragment identified in human circulation. *Peptides*. 2006;27: 1835-40.
- Donato LJ, Meeusen JW, Lieske JC, Bergmann D, Sparwasser A, Jaffe AS. Analytical performance of an immunoassay to measure proenkephalin. *Clin Biochem*. 2018;58: 72-77.
- Beunders R, Struck J, Wu AH, Zarbock A, Di Somma S, Mehta RL, Koyner JL, Nadim MK, Maisel AS, Murray PT, Neath SX, Jaffe AS, Pickkers P. Proenkephalin (PENK) as a novel biomarker for kidney function. *J Appl Lab Med*. 2017; DOI: 10.1373/jalm.2017.023598 400-12.
- Matsue Y, Ter Maaten JM, Struck J, Metra M, O'Connor CM, Ponikowski P, Teerlink JR, Cotter G, Davison B, Cleland JG, Givertz MM, Bloomfield DM, Ditttrich HC, van Veldhuisen DJ, van der Meer P, Damman K, Voors AA. Clinical correlates and prognostic value of proenkephalin in acute and chronic heart failure. *J Card Fail*. 2017;23: 231-39.
- Kieneker LM, Hartmann O, Struck J, Bergmann A, Gansevoort RT, Joosten MM, van den Berg E, de Boer RA, Bakker SJL. Plasma proenkephalin and poor long-term outcome in renal transplant recipients. *Transplant Direct*. 2017;3: e190.
- Leitje G, Beunders R, van Groenendal R, Peters E, Kox M, Pickkers P. Proenkephalin, the new marker for kidney function on the Intensive care unit? *Intensive Care Med Exp*. 2017;5 (Suppl 2): 51-52.
- Beunders R, Meekes M, Struck J, Pickkers P. Proenkephalin, a novel biomarker for kidney function, is earlier in detecting acute kidney injury compared to creatinine. *Intensive Care Med Exp*. 2019;7 (Suppl 3): 000960.
- Marino R, Struck J, Hartmann O, Maisel AS, Rehfeldt M, Magrini L, Melander O, Bergmann A, Di Somma S. Diagnostic and short-term prognostic utility of plasma pro-enkephalin (pro-ENK) for acute kidney injury in patients admitted with sepsis in the emergency

department. *J Nephrol.* 2015;28: 717-24.

9. Hollinger A, Wittebole X, Francois B, Pickkers P, Antonelli M, Gayat E, Chousterman BG, Lascarrou JB, Dugernier T, Di Somma S, Struck J, Bergmann A, Beishuizen A, Constatin JM, Damoiseil C, Deye N, Gaudry S, Huberlant V, Marx G, Mercier E, Oueslati H, Hartmann O, Sonnevile R, Laterre PF, Mebazaa A, Legrand M. Proenkephalin A 119-159 (Penkid) is an early biomarker of septic acute kidney injury: The Kidney in Sepsis and Septic Shock (Kid-SSS) study. *Kidney Int Rep.* 2018;3: 1424-33.
10. Gayat E, Cariou A, Deye N, Vieillard-Baron A, Jaber S, Damoiseil C, Lu Q, Monnet X, Rennuit I, Azoulay E, Leone M, Oueslati H, Guidet B, Friedman D, Tesniere A, Sonnevile R, Montravers P, Pili-Floury S, Lefrant JY, Duranteau J, Laterre PF, Brechot N, Chevreul K, Michel M, Cholley B, Legrand M, Launay JM, Vicaut E, Singer M, Resche-Rigon M, Mebazaa A. Determinants of long-term outcome in ICU survivors: results from the FROG-ICU study. *Crit Care.* 2018;22: 8.
11. Ng LL, Sandhu JK, Narayan H, Quinn PA, Squire IB, Davies JE, Bergmann A, Maisel A, Jones DJ. Proenkephalin and prognosis after acute myocardial infarction. *J Am Coll Cardiol.* 2014;63: 280-9.
12. Shah KS, Taub P, Patel M, Rehfeldt M, Struck J, Clopton P, Mehta RL, Maisel AS. Proenkephalin predicts acute kidney injury in cardiac surgery patients. *Clin Nephrol.* 2015;83: 29-35.
13. Ng LL, Squire IB, Jones DJL, Cao TH, Chan DCS, Sandhu JK, Quinn PA, Davies JE, Struck J, Hartmann O, Bergmann A, Mebazaa A, Gayat E, Arrigo M, Akiyama E, Sabti Z, Lohrmann J, Twerenbold R, Herrmann T, Schumacher C, Kozhuharov N, Mueller C, and Great Network. Proenkephalin, renal dysfunction, and prognosis in patients with acute heart failure: A GREAT Network study. *J Am Coll Cardiol.* 2017;69: 56-69.
14. Breidthardt T, Jaeger C, Christ A, Klima T, Mosimann T, Twerenbold R, Boeddinghaus J, Nestelberger T, Badertscher P, Struck J, Bergmann A, Hartmann O, Kalbermatter S, Marenzi G, Mueller C. Proenkephalin for the early detection of acute kidney injury in hospitalized patients with chronic kidney disease. *Eur J Clin Invest.* 2018: e12999.
15. Emmens JE, Ter Maaten JM, Damman K, van Veldhuisen DJ, de Boer RA, Struck J, Bergmann A, Sama IE, Streng KW, Anker SD, Dickstein K, Lang CC, Metra M, Samani NJ, Ng LL, Voors AA. Proenkephalin, an opioid system surrogate, as a novel comprehensive renal marker in heart failure. *Circ Heart Fail.* 2019;12: e005544.
16. Holaday, JW. Cardiovascular effects of endogenous opiate systems. *Annu Rev Pharmacol Toxicol.* 1983;23: 541-94.
17. Sezen, SF, Kenigs VA, Kapusta DR. Renal excretory responses produced by the delta opioid agonist, BW373U86, in conscious rats. *J Pharmacol Exp Ther.* 1998;287: 238-45.
18. Denning, GM, Ackermann LW, Barna TJ, Armstrong JG, Stoll LL, Weintraub NL, Dickson EW. Proenkephalin expression and enkephalin release are widely observed in non-neuronal tissues. *Peptides.* 2008;29: 83-92.
19. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45: 942-6.
20. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15: 108-15.
21. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
22. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
23. McEnroe RJ, Dimeski G, Miller WG, et al. Supplemental Tables for Interference Testing in Clinical Chemistry. CLSI Document EP37, First Edition, 2018.
24. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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**Coefficiente de correlación,  $r = 0,92$**

## REFERENCIAS

1. Ernst A, Köhrle J, Bergmann A. Proenkephalin A 119-159, a stable proenkephalin A precursor fragment identified in human circulation. *Peptides*. 2006;27: 1835-40.
2. Donato LJ, Meeusen JW, Lieske JC, Bergmann D, Sparwasser A, Jaffe AS. Analytical performance of an immunoassay to measure proenkephalin. *Clin Biochem*. 2018;58: 72-77.
3. Beunders R, Struck J, Wu AH, Zarbock A, Di Somma S, Mehta RL, Koyner JL, Nadim MK, Maisel AS, Murray PT, Neath SX, Jaffe AS, Pickkers P. Proenkephalin (PENK) as a novel biomarker for kidney function. *J Appl Lab Med*. 2017; DOI: 10.1373/jalm.2017.023598 400-12.
4. Matsue Y, Ter Maaten JM, Struck J, Metra M, O'Connor CM, Ponikowski P, Teerlink JR, Cotter G, Davison B, Cleland JG, Givertz MM, Bloomfield DM, Dittrich HC, van Veldhuisen DJ, van der Meer P, Damman K, Voors AA. Clinical correlates and prognostic value of proenkephalin in acute and chronic heart failure. *J Card Fail*. 2017;23: 231-39.
5. Kieneker LM, Hartmann O, Struck J, Bergmann A, Gansevoort RT, Joosten MM, van den Berg E, de Boer RA, Bakker SJL. Plasma proenkephalin and poor long-term outcome in renal transplant recipients. *Transplant Direct*. 2017;3: e190.
6. Leitje G, Beunders R, van Groenendal R, Peters E, Kox M, Pickkers P. Proenkephalin, the new marker for kidney function on the Intensive care unit? *Intensive Care Med Exp*. 2017;5 (Suppl 2): 51-52.
7. Beunders R, Meekes M, Struck J, Pickkers P. Proenkephalin, a novel biomarker for kidney function, is earlier in detecting acute kidney injury compared to creatinine. *Intensive Care Med Exp*. 2019;7 (Suppl 3): 000960.
8. Marino R, Struck J, Hartmann O, Maisel AS, Rehfeldt M, Magrini L, Melander O, Bergmann A, Di Somma S. Diagnostic and short-term prognostic utility of plasma pro-enkephalin (pro-ENK) for acute kidney injury in patients admitted with sepsis in the emergency department. *J Nephrol*. 2015;28: 717-24.
9. Hollinger A, Wittebole X, Francois B, Pickkers P, Antonelli M, Gayat E, Chousterman BG, Lascarrou JB, Dugernier T, Di Somma S, Struck J,

- Bergmann A, Beishuizen A, Constatin JM, Damoiseil C, Deye N, Gaudry S, Huberlant V, Marx G, Mercier E, Oueslati H, Hartmann O, Sonnevile R, Laterre PF, Mebazaa A, Legrand M. Proenkephalin A 119-159 (Penkid) is an early biomarker of septic acute kidney injury: The Kidney in Sepsis and Septic Shock (Kid-SSS) study. *Kidney Int Rep.* 2018;3: 1424-33.
10. Gayat E, Cariou A, Deye N, Vieillard-Baron A, Jaber S, Damoiseil C, Lu Q, Monnet X, Rennuit I, Azoulay E, Leone M, Oueslati H, Guidet B, Friedman D, Tesniere A, Sonnevile R, Montravers P, Pili-Floury S, Lefrant JY, Duranteau J, Laterre PF, Brechot N, Chevreul K, Michel M, Cholley B, Legrand M, Launay JM, Vicaut E, Singer M, Resche-Rigon M, Mebazaa A. Determinants of long-term outcome in ICU survivors: results from the FROG-ICU study. *Crit Care.* 2018;22: 8.
  11. Ng LL, Sandhu JK, Narayan H, Quinn PA, Squire IB, Davies JE, Bergmann A, Maisel A, Jones DJ. Proenkephalin and prognosis after acute myocardial infarction. *J Am Coll Cardiol.* 2014;63: 280-9.
  12. Shah KS, Taub P, Patel M, Rehfeldt M, Struck J, Clopton P, Mehta RL, Maisel AS. Proenkephalin predicts acute kidney injury in cardiac surgery patients. *Clin Nephrol.* 2015;83: 29-35.
  13. Ng LL, Squire IB, Jones DJL, Cao TH, Chan DCS, Sandhu JK, Quinn PA, Davies JE, Struck J, Hartmann O, Bergmann A, Mebazaa A, Gayat E, Arrigo M, Akiyama E, Sabti Z, Lohrmann J, Twerenbold R, Herrmann T, Schumacher C, Kozhuharov N, Mueller C, and Great Network. Proenkephalin, renal dysfunction, and prognosis in patients with acute heart failure: A GREAT Network study. *J Am Coll Cardiol.* 2017;69: 56-69.
  14. Breidthardt T, Jaeger C, Christ A, Klima T, Mosimann T, Twerenbold R, Boeddinghaus J, Nestelberger T, Badertscher P, Struck J, Bergmann A, Hartmann O, Kalbermatter S, Marenzi G, Mueller C. Proenkephalin for the early detection of acute kidney injury in hospitalized patients with chronic kidney disease. *Eur J Clin Invest.* 2018: e12999.
  15. Emmens JE, Ter Maaten JM, Damman K, van Veldhuisen DJ, de Boer RA, Struck J, Bergmann A, Sama IE, Streng KW, Anker SD, Dickstein K, Lang CC, Metra M, Samani NJ, Ng LL, Voors AA. Proenkephalin, an opioid system surrogate, as a novel comprehensive renal marker in heart failure. *Circ Heart Fail.* 2019;12: e005544.
  16. Holaday, JW. Cardiovascular effects of endogenous opiate systems. *Annu Rev Pharmacol Toxicol.* 1983;23: 541-94.
  17. Sezen, SF, Kenigs VA, Kapusta DR. Renal excretory responses produced by the delta opioid agonist, BW373U86, in conscious rats. *J Pharmacol Exp Ther.* 1998;287: 238-45.
  18. Denning, GM, Ackermann LW, Barna TJ, Armstrong JG, Stoll LL, Weintraub NL, Dickson EW. Proenkephalin expression and enkephalin release are widely observed in non-neuronal tissues. *Peptides.* 2008;29: 83-92.
  19. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45: 942-6.
  20. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15: 108-15.
  21. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
  22. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
  23. McEnroe RJ, Dimeski G, Miller WG, et al. Supplemental Tables for Interference Testing in Clinical Chemistry. CLSI Document EP37, First Edition, 2018.
  24. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.