English

IB10 sphingotest® 3-in-1 Cardiac Troponin I/CK-MB/Myoglobin

For the Quantitative Determination of
Troponin I, CK-MB, Myoglobin in
Lithium Heparin Human Whole Blood and Plasma

Explanation of Symbols

	CE Mark of Conformity
	Manufacturer
	Catalog Number
	Expiry date/Use by
	Lot number
	<i>In Vitro Diagnostic</i> medical device
	Consult instructions for use
	Store between 2 °C and 8 °C
	European Authorized Representative
	Contains sufficient for <n> tests
	Do not reuse
	Serial Number
	Insert disc label side up

Nexus Dx, Inc.,
6759 Mesa Ridge Road
San Diego, California 92121 USA
Telephone: 1 (858) 410 4600
Fax: 1 (858) 410 4700

TheraGenesis GmbH
Bahnhofstrasse 5
55276 Oppenheim, Germany
Tel: +49 (0) 151 506 403 14

022-00065 REV B

IB10 sphingotest® 3-in-1 Cardiac

Troponin I/CK-MB/Myoglobin

For *In Vitro* Diagnostic Use

INTENDED USE

The IB10 sphingotest® 3-in-1 Cardiac is a rapid point-of-care (POC) immunoassay for the *in vitro* quantitative determination of Cardiac Troponin I (cTnI), CK-MB and Myoglobin, in lithium heparin whole blood and plasma. This Test is intended for use in conjunction with the Nexus IB10 Analyzer and provides quantitative results in 20 minutes.

The IB10 sphingotest® 3-in-1 Cardiac is designed for professional use only and may be used in hospital central laboratories and in alternate care settings such as emergency departments, critical care units, and other sites where near patient testing is practiced.

SUMMARY AND EXPLANATION OF THE TEST

Acute coronary syndrome (ACS), including Coronary Heart Disease (CHD) and Angina Pectoris (AP), is the leading cause of morbidity and mortality among both men and women, affecting more than 16.0 million people in the United States. The presentation of ACS is varied, with acute myocardial infarction (AMI) being the most dramatic presentation. Symptoms of ACS may develop suddenly. In 2008, there were over 800,000 new AMIs in the US.¹ On an international basis, the World Health Organization (WHO) reports over 7.2 million deaths a year from CHD. Globally it is estimated that 17.3 million deaths were due to Cardiovascular Disease in 2008.² This includes death due to both CHD and stroke.

For over 30 years, the WHO has recommended that the diagnosis of AMI be guided by positive findings on at least two of the following three criteria: (1) patient history/physical examination; (2) electrocardiogram; (3) changes in blood levels of cardiac protein markers.³ While the original cardiac marker was defined as total enzymatic creatine kinase, today's markers, as measured by immunoassay, include Myoglobin (Myo), Creatine Kinase isoform MB (CK-MB) and the cardiac striated muscle protein, Troponin I (cTnI). Patient history and physical examination are critical, but often provide insufficient information to differentiate AMI from other cardiac abnormalities. In the absence of well defined electrocardiographic ST-segment elevations, measurements of cardiac protein markers are essential for the accurate diagnosis of ACS.⁴

Several studies have suggested that a panel of all three cardiac markers enhances the reliability of earlier identification and risk stratification of patients presenting with chest pain compared to a single marker alone.^{5,6} This is due to differences in the temporal pattern of release into the bloodstream of various cardiac biomarkers following myocyte injury or cell death.

MYOGLOBIN

Myoglobin (Myo) is an oxygen-binding heme protein with a molecular weight of 17,800 daltons, typically found in cardiac and skeletal tissues. The protein is localized in the cytosol of the cell and constitutes about 2% of the total muscle protein. Because of its relatively low molecular weight and high concentration, it is among the first proteins released into the systemic blood upon cardiac myocyte cellular necrosis (cell death).⁷ Myoglobin is found in the blood at an elevated concentration approximately 2-3 hours following cardiac myocyte injury or necrosis with levels peaking after 9-12 hours.⁸

CK-MB

In the cytoplasm, creatine kinase (CK) catalyzes the reversible phosphorylation of creatine with ATP. Systolic creatine kinase exists as a dimeric molecule formed from two types of single polypeptide subunits. These two subunits combine to form three different isoenzymes of creatine kinase: CK-MM, CK-BB and CK-MB. CK-MM is predominantly expressed in skeletal muscle tissue, while CK-MB isoenzyme is most abundant in cardiac muscle tissue. CK-MB is released into the blood stream upon damage or death of cardiac tissue.^{9,10} Its elevation is noted within 4 to 8 hours following cardiac injury, reaching maximum concentrations 12-24 hours after onset of symptoms. Like Myoglobin, CK-MB is not completely cardiac specific and is found in low concentrations in striated skeletal muscle and other tissue.

cTnI

Cardiac Troponin I (cTnI) is part of the troponin complex which, together with tropomyosin, forms the main component that regulates the Ca²⁺ sensitive ATP-ase activity of actomyosin in striated muscle (skeletal and cardiac). The troponin complex consists of three subunits, troponin T (TnT), troponin I (TnI), and troponin C (TnC). Each subunit has a distinct function with TnC as the site of Ca²⁺ binding, TnT the tropomyosin binding, and TnI as the inhibitory subunit.¹¹ Different isoforms of TnI exist in the skeletal and cardiac muscles (sTnI and cTnI, respectively) with distinct immunologic epitopes that allow the production of extremely specific cardiac TnI antibodies.¹²

cTnI determinations aid in the diagnosis of myocardial infarction (MI) and patients with non-ST-segment elevation (NSTEMI) acute coronary syndrome. Elevated cTnI levels in the NSTEMI ACS subpopulation correlate with the relative risk of mortality, MI or increased probability of ischemic events requiring urgent revascularization procedures.

Following myocyte necrosis, cTnI is released into the circulation with levels exceeding the upper reference limit of normal within 4-6 hours. Peak levels are reached after 12-24 hours.¹³ This release profile into the peripheral blood is similar to CK-MB. However, CK-MB levels return to normal values in about 72 hours, while levels of cTnI may remain elevated for 5-7 days.

TnI in Acute Coronary Syndrome

Since the introduction of the first commercial assays for Troponin I in 1996, numerous clinical studies have confirmed that cTnI is the biomarker of choice in identifying cardiac myocyte ischemia and necrosis.¹⁴ Even low levels of cTnI just above the 99th percentile of levels obtained from a normal reference population indicate increased risk. To assure that the test results are interpreted with as much accuracy as possible, consensus groups now unanimously recommend that chest pain patients with low level troponin values on presentation to medical personnel, especially in the absence of confirmatory ECG evidence, be tested 2 additional times during the subsequent 12-24 hour period to determine whether cTnI is indeed rising during the first 24 hours.^{4,15-17}

PRINCIPLE

The Nexus Dx immunochemistry system combines chemistry with microfluidics and centrifugal flow to rapidly prepare a cell free plasma from whole blood that can then be moved through a channel to rehydrate, solubilize and mix with freeze dried immunoconjugates. Using a combination of active flow and capillary action, the test sample is quantitatively measured in 20 minutes with an optical signal level proportional to the analyte(s) concentration.

After addition of the patient sample, the entire test is performed within the Nexus IB10 Analyzer which provides control of the temperature of the disc, as well as the analysis process sequence, centrifugal flow, mixing, incubation time, final signal measurement, quantitation and reporting of results.

The Test disc includes a positive internal control to ensure that the test has operated properly. Each lot is calibrated to ensure that lot-to-lot variability is minimized. Lot specific calibration along with additional information such as the lot expiration date is contained on a QR code label affixed to each disc. It is recommended that external controls be tested at appropriate time intervals to confirm that the system and test lot are performing within acceptable limits.

REAGENTS

The IB10 sphingotest® 3-in-1 Cardiac contains all required reagents to evaluate the levels of cTnI, CK-MB and Myoglobin.

For the cTnI Test, reagents include: dye conjugated monoclonal anti-cTnI antibody, biotin conjugated polyclonal anti-cTnI antibodies and streptavidin immobilized at the detection area on the disc.

For the CK-MB Test, reagents include: dye conjugated monoclonal anti-CK-MB antibody, biotin conjugated goat polyclonal anti-CK-MM antibodies and streptavidin immobilized at the detection area on the disc.

For the Myoglobin Test, reagents include: dye conjugated monoclonal anti-myoglobin antibody, biotin conjugated polyclonal rabbit anti-myoglobin antibodies and streptavidin immobilized at the detection area on the disc.

MATERIALS PROVIDED

Each box contains the following:

- 10 IB10 sphingotest® 3-in-1 Cardiac discs, each individually sealed in a foil pouch with a desiccant.
- 1 Instructions for Use (IFU).

MATERIALS/EQUIPMENT REQUIRED BUT NOT PROVIDED

1. Nexus IB10 Analyzer - Model #BCA-IB10.
2. Commercially available Troponin I/CK-MB/Myoglobin controls for external Quality Control (QC). **Contact the Distributor in your area** for recommended external QC materials or related technical assistance.
3. Calibrated reusable fixed or variable volume pipette gun with high precision and accuracy capable of delivering 500 µL of whole blood or plasma.
4. Disposable pipette tips capable of accepting and delivering 500 µL of whole blood or plasma.

PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Carefully follow the Instructions for Use.
- Prior to testing controls or patient samples, ensure that the Analyzer software is updated to the latest version. Refer to the Nexus IB10 manual for specific instructions.

- Wear disposable gloves while handling samples.
- Handle samples with care. Samples and used Test discs should be treated as potentially infectious and should be discarded as biohazardous material according to local regulations.
- Thoroughly wash hands following handling.
- The result obtained from the IB10 sphingotest® 3-in-1 Cardiac does not provide a definitive diagnosis and should be interpreted by a physician in conjunction with other appropriate test results according to current medical Guidelines and patient clinical findings.
- Keep the Test in the sealed pouch until ready for use.
- Do not use the Test if the pouch is damaged or the seal is broken.
- Do not use the Test after the expiration date printed on the pouch.
- Prior to use, place the unopened pouch at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
- Always pay attention to cleanliness when handling the Test. Avoid any contamination from fingerprints or foreign substances. Do not contaminate the sample channel inlet.
- Do not drop or damage the Test disc.
- The Test disc should be inserted, with the label side up, into the Nexus IB10 Analyzer tray immediately after injecting the sample into the disc.
- Do not turn the disc upside down.
- This is a quantitative test; therefore no visual interpretation of the results should be made.

STORAGE AND STABILITY

- Store the IB10 sphingotest® 3-in-1 Cardiac between 2 and 8 °C (35 to 46 °F) until the expiration date printed on the pouch is reached.
- The IB10 sphingotest® 3-in-1 Cardiac in its sealed pouch is stable at 18 to 30 °C/64 to 86 °F for 30 days, provided the expiration date printed on the pouch is not exceeded.

SAMPLE COLLECTION AND PREPARATION

- The IB10 sphingotest® 3-in-1 Cardiac is to be run using lithium heparin whole blood or plasma samples.
- As cardiac proteins are relatively unstable, it is recommended that samples be tested as soon as possible.
- Whole blood must be tested within 48 hours of collection.
- Whole blood and plasma samples can be transported and stored up to 48 hours at ambient temperature.
- Plasma samples should be kept frozen at -20 °C (- 4 °F) or lower if longer storage is required.
- Allow samples to equilibrate to room temperature (19 to 25 °C/66 to 77 °F) prior to testing.
- For serial testing of the same patient during the 0–12 hour triage period, the same sample type (whole blood or plasma) should be used.

PROCEDURE

Nexus IB10 Analyzer

Consult the Nexus IB10 Analyzer User Manual

For Analyzer installation, start up and complete instructions for use refer to the **Nexus IB10 Analyzer User Manual**. Operator must consult the User Manual prior to use to become familiar with the proper operation and quality control procedures.

PERFORMING SYSTEM CHECK AND DISC CALIBRATION

Each time the Nexus IB10 Analyzer is turned on, a Self Check is automatically performed. The QR code on each Test disc contains information for disc calibration which the Analyzer automatically reads when running a test.

RUNNING QC WITH EXTERNAL CONTROLS

The manufacturer recommends the use of commercially available Troponin I/CK-MB/Myoglobin controls (please refer to section **Materials/Equipment Required But Not Provided**). Ensure that the Troponin I/CK-MB/Myoglobin Controls are handled and prepared according to the corresponding Instructions for Use (IFU).

1. Remove an unopened Test pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes prior to testing.
2. Open the pouch and remove the Test disc.

3. Place the Test disc on a level surface.
 4. On the Nexus IB10 Analyzer press **New Analysis**.
 5. The Analyzer will perform a general system check.
 6. Enter the external control ID manually (up to 14 characters can be used for the ID) or enter the external control ID by using the barcode scanner.
 7. Mix the external quality control vial by gently inverting the vial several times before sampling.
 8. **Testing External Control Sample on the IB10 sphingotest® 3-in-1 Cardiac.**
 - Using a precision pipette (fixed or adjusted to 500 µL) slowly draw the well mixed external quality control sample into the pipette tip.
 - Positioning the tapered pipette tip at a 45° angle, pierce the **X** on the red dot to expose the sample channel inlet.
 - Slowly express the external control sample into the inlet applying minimum, but continuous force, on the pipette plunger.
 - Express sample to the first stop on the pipette, at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
 - Press **OK** on the Nexus IB10 Analyzer display.
 - When the tray opens, insert the filled Test disc into the tray and press **Run**.
 - In 20 minutes, the Nexus IB10 Analyzer will display the results on the screen.
 - Results will print automatically (if selected during Set Up) or press **Print**.
 - When the test is complete, analyze and compare the result with the Expected Value reported in the external control IFU for the external control level as measured using the IB10 sphingotest® 3-in-1 Cardiac.
 - Remove the Test disc and discard in appropriate receptacle.
 - If the external control result is outside the expected range, refer to the Quality Control section below.
- Note: If the test run is cancelled before a Test result is displayed, the Test disc cannot be reused and should be disposed of appropriately.*

TESTING PATIENT SAMPLES ON THE NEXUS IB10 ANALYZER

1. Remove an unopened Test pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
2. Open the pouch and remove the Test disc.
3. Place the Test disc on a level surface.
4. On the Nexus IB10 Analyzer press **New Analysis**.
5. The Analyzer will perform a general system check.
6. Enter the Patient ID manually (up to 14 characters can be used for the ID) or enter the Patient ID by using the barcode scanner.
7. Mix the whole blood patient sample tube by gently inverting the tube several times before testing.
8. **Testing Patient Samples on the IB10 sphingotest® 3-in-1 Cardiac.**
 - Using a precision pipette (fixed or adjusted to 500 µL) slowly draw the well mixed patient sample into the pipette tip.
 - Positioning the tapered pipette tip at a 45° angle, pierce the **X** on the red dot to expose the sample channel inlet.
 - Slowly express the patient sample into the inlet applying minimum, but continuous force, on the pipette plunger.
 - Express sample to the first stop on the pipette, at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
 - Press **OK** on the Nexus IB10 Analyzer display.
 - When the tray opens, insert the filled Test disc into the tray and press **Run**.
 - In 20 minutes, the Nexus IB10 Analyzer will display the results on the screen.
 - Results will print automatically (if selected during Set Up) or press **Print**.
 - Remove the Test disc and discard in appropriate receptacle.

Note: If the test run is cancelled before a test result is displayed, the Test disc cannot be reused and should be disposed of appropriately.

Quality Control

EXTERNAL CONTROLS

Good laboratory practice includes the use of external controls to ensure proper test performance. It is recommended that prior to using a new lot of IB10 sphingotest® 3-in-1 Cardiac the performance of the lot should be confirmed by testing with external controls (see section **Materials/Equipment Required But Not Provided**) to ensure the test will deliver the correct result. The frequency of quality control testing should be determined according to individual laboratory standard quality control procedures. Upon confirmation of the expected results, the Test discs are

ready for use with patient samples. Controls should also be used any time the validity of test results is questionable. If external controls do not perform as expected, do not use the IB10 sphingotest® 3-in-1 Cardiac and contact the Distributor in your area for Technical Assistance.

INTERNAL CONTROL

Each IB10 sphingotest® 3-in-1 Cardiac has a built in positive procedural control. The Nexus IB10 Analyzer automatically determines the presence of this control thereby confirming that the test run has delivered a valid result. If the control does not form or if the Analyzer does not recognize it, the test result is considered "invalid" and the test must be repeated.

LIMITATIONS

The results of the IB10 sphingotest® 3-in-1 Cardiac should be used in conjunction with other available laboratory and clinical information. Universally accepted Good Clinical Practice Guidelines recommend collecting serial samples during the early triage period from 0 to > 12 hours.^{4,15-17}

Other substances and/or factors not listed, e.g. technical or procedural errors, may interfere with the test and cause inaccurate results.

Nexus Dx, Inc. offers products for their intended use. Refer to the specific product literature for the intended use statements for each product. Product claims are subject to change. Nexus Dx, Inc.'s expressed and implied warranties (inclusive of implied warranties of merchantability and fitness) are conditional upon adherence to, or observance of Nexus Dx, Inc. published directions with regard to the use of Nexus Dx, Inc. products. Under no circumstances will Nexus Dx, Inc. be liable for any indirect or consequential damages.

For Technical Assistance please contact your local Distributor.

Performance Characteristics

MEASURING RANGE

The IB10 sphingotest® 3-in-1 Cardiac has been demonstrated to provide measurable results at Troponin I/CK-MB/Myoglobin levels of:

0.05 ng/mL - 30 ng/mL - cTnI
2.0 ng/mL - 60 ng/mL - CK-MB
30.0 ng/mL - 500 ng/mL - Myoglobin

ANALYTICAL SENSITIVITY

The LoD (Limit of Detection) for each analyte was determined according to Clinical and Laboratory Standards Institute (CLSI) Guideline EP17-A.¹⁸ The proportion of false positives (α) and false negatives (β) are less than 5%.

The LoD of the IB10 sphingotest® 3-in-1 Cardiac for each analyte is:

0.05 ng/mL - cTnI
2.0 ng/mL - CK-MB
30.0 ng/mL - Myoglobin

LoQ (Limit of Quantitation) is the lowest concentration that can be reproducibly measured with a total coefficient of variation of at most 15%.

For the 3 analytes LoQ was determined to be:

0.10 ng/mL - cTnI
3.0 ng/mL - CK-MB
50.0 ng/mL - Myoglobin

CROSS-REACTIVITY AND INTERFERING SUBSTANCES

Drugs

The following drugs were tested for potential interference in the IB10 sphingotest® 3-in-1 Cardiac (Table 1). The list of drugs encompasses common prescription and over-the-counter medications, as well as those often prescribed in a cardiac patient population. The drugs were tested at concentrations recommended in the CLSI Approved Guideline EP7-A2 'Interference Testing in Clinical Chemistry,'¹⁹ or at least three times the highest concentration reported following a therapeutic dose. No significant interference with the IB10 sphingotest® 3-in-1 Cardiac measurement was observed for the drugs listed in the table below.

Table 1.

Drug	Drug	Drug
Acetaminophen	Captopril	Nifedipine
Acetylsalicylic acid (aspirin)	Digoxin	Phenytoin
Ampicillin	Dopamine	Theophylline
Ascorbic acid (vitamin C)	Erythromycin	Verapamil
Atenolol	Furosemide	
Caffeine	Methyl-DOPA	

PROTEINS, PEPTIDES AND ENDOGENOUS SUBSTANCES

The following proteins, peptides and endogenous substances were tested for potential cross-reactivity and interference in the IB10 sphingotest® 3-in-1 Cardiac at the maximum concentration of the substance indicated (Table 2). No substance demonstrated significant cross-reactivity or interference using the method provided in CLSI EP7-A2.

Table 2.

Substance	Maximum Concentration	
Proteins and Peptides		
Cardiac Troponin T	1000 ng/mL	
Cardiac Troponin C	1000 ng/mL	
Skeletal Troponin I	250 ng/mL	
CKMM	781 ng/mL	
CKBB	5000 ng/mL	
Endogenous substances		
Albumin	2.5 g/dL	
Hemoglobin	0.1 g/dL	
Rheumatoid Factor (RF)	203 IU/mL	
Creatinine	2 mg/dL	
Bilirubin	2.5 mg/dL	
Triglycerides	0.5 g/dL	
Urea	18 mg/dL	
Cholesterol	280 mg/dL	
Substance	Analyte	Maximum Concentration
Biotin	Myoglobin	400 ng/mL
	CK-MB	200 ng/mL

HOOK EFFECT

No high dose Hook Effect was observed for concentrations up to 500 ng/mL for cTnI, 200 ng/mL for CK-MB and 4000 ng/mL for Myoglobin.

PRECISION

The precision of the IB10 sphingotest® 3-in-1 Cardiac was determined using samples where cTnI, CK-MB and Myoglobin were added to normal human plasma at multiple concentrations (Table 3). The within-day and total precision were performed in two runs per day, in replicates of 4 per run at each concentration level over a 15 day period for a total number of repetitions of 120 at each concentration level. The within-run, total variances and coefficients of variation (CVs) were computed according to CLSI Guideline EP5-A2.²⁰

Table 3A.

IB10 sphingotest® 3-in-1 Cardiac - Assay Precision for cTnI

Sample	Mean (ng/mL)	Within-Run Precision		Total Precision	
		Std. dev. (ng/mL)	CV%	Std. dev. (ng/mL)	CV%
1	0.46	0.06	12.6%	0.06	12.6%
2	0.95	0.11	12.0%	0.12	13.0%
3	4.43	0.41	9.4%	0.45	10.1%

Table 3B.

IB10 sphingotest® 3-in-1 Cardiac - Assay Precision for CK-MB

Sample	Mean (ng/mL)	Within-Run Precision		Total Precision	
		Std. dev. (ng/mL)	CV%	Std. dev. (ng/mL)	CV%
1	7.13	0.64	9.0	0.73	10.3
2	13.33	1.49	11.2	1.53	11.5

Table 3C.

IB10 sphingotest® 3-in-1 Cardiac - Assay Precision for Myoglobin

Sample	Mean (ng/mL)	Within-Run Precision		Total Precision	
		Std. dev. (ng/mL)	CV%	Std. dev.(ng/mL)	CV%
1	106.1	11.4	10.8	12.8	12.1
2	202.7	25.6	12.7	25.6	12.7

WHOLE BLOOD Vs. PLASMA COMPARISON

A comparison study was performed using lithium heparin whole blood and plasma samples. When performing a Passing-Bablok regression analysis comparing the whole blood (WB) concentrations versus the corresponding plasma concentrations from the same subject samples, the following relationship was determined:

$$\text{WB} = 1.00 \text{ Plasma (95\% C.I. = [0.93-1.07])} + 0.01 \text{ ng/mL - cTnI}$$

$$\text{WB} = 1.17 \text{ Plasma (95\% C.I. = [1.09-1.29])} - 0.19 \text{ ng/mL - CK-MB}$$

$$\text{WB} = 0.92 \text{ Plasma (95\% C.I. = [0.83-1.04])} + 0.00 \text{ ng/mL - Myoglobin}$$

Expected Values

UPPER REFERENCE LIMIT

From a population of 224 individuals, the IB10 sphingotest® 3-in-1 Cardiac was used to determine the concentration upper reference limits of cTnI, CK-MB and Myoglobin. This population included apparently healthy individuals. The upper reference limits for each analyte is:

0.10 ng/mL (99th reference percentile) - cTnI

8.58 ng/mL (95th reference percentile) - CK-MB

99.84 ng/mL (95th reference percentile) - Myoglobin

Each laboratory should establish a reference range that represents the patient population that is to be evaluated at their facility.

METHOD COMPARISON

An equivalence study was performed between the IB10 sphingotest® 3-in-1 Cardiac using the Nexus IB10 Analyzer and the Ortho VITROS® Troponin I ES, Myoglobin and CK-MB Tests.

cTnI

A total of 253 samples were tested within a concentration range of cTnI of 0.05 ng/mL to 30 ng/mL. The Passing-Bablok regression was:

$$\text{IB10 sphingotest® 3-in-1 Cardiac (TnI)} = 0.80 \text{ (VITROS TnI ES)} + 0.00 \text{ ng/mL}$$

Correlation coefficient, $r = 0.92$

CK-MB

A total of 161 samples were tested within measuring range of CK-MB. The Passing-Bablok regression was:

$$\text{IB10 sphingotest® 3-in-1 Cardiac (CK-MB)} = 1.29 \text{ (VITROS CK-MB)} - 0.54 \text{ ng/mL}$$

Correlation coefficient, $r = 0.88$

MYOGLOBIN

A total of 117 samples were tested within measuring range of Myoglobin. The Passing-Bablok regression was:

$$\text{IB10 sphingotest® 3-in-1 Cardiac (Myo)} = 1.00 \text{ (VITROS MYOG)} - 5.1 \text{ ng/mL}$$

Correlation coefficient, $r = 0.90$

REFERENCES

1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. *Circulation* 2012;125:e2-e220.
2. Cardiovascular Diseases, World Health Organization factsheet. <http://www.who.int/mediacentre/factsheets/fs317/en/index.html>
3. World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria of diagnosis of ischemic heart disease. *Circulation* 1979;59:607-9.
4. Thygesen K, Alpert JS, White HD, Joint ESC/ACCF/AHA/WHF Task Force for the Redefinition of Myocardial Infarction. Universal definition of myocardial infarction. *Eur Heart J* 2007;28:2525-38; *Circulation* 2007;116:2634-53; *J Am Coll Cardiol* 2007;50:2173-95.
5. Newby LK, Storrow AB, Gibler B, et al. Bedside multimarker testing for risk stratification in chest pain patients – the chest pain evaluation by creatine kinase-MB, myoglobin and troponin I (CHECKMATE) study. *Circulation* 2001;103:1832-7.
6. Apple FS, Voss E, Lund L, et al. Cardiac troponin, CK-MB and myoglobin for the early detection of acute myocardial infarction and monitoring of reperfusion following thrombolytic therapy. *Clin Chim Acta* 1995;237:59-66.
7. Ohman EM, Casy C, Bengston JR, et al. Early detection of acute myocardial infarction: Additional diagnostic information from serum concentrations of myoglobin in patients without ST elevations. *Br Heart J* 1990;63:335-8.
8. Tucker JF, Collins RA, Anderson AJ, et al. Value of serial myoglobin levels in the early diagnosis of patients admitted for acute myocardial infarction. *Ann Emerg Med* 1994;24:704-8.
9. Cairns JA, Missirrli E, Walker WH. Usefulness of serial determinations of myoglobin and creatine kinase in serum compared for assessment of acute myocardial infarction. *Clin Chem* 1983;29:469-73.
10. Lott JA, Abbott LB. Creatine kinase isoenzymes. *Clin Lab Med* 1986;6:547-76.
11. Mehegan JP, Tobacman LS. Cooperative interaction between troponin molecules bound to the cardiac thin filament. *J Biol Chem* 1991;266:966-72.
12. Bodor GS, Porter S, Landt Y, Ladenson JH. Development of monoclonal antibody for an assay of cardiac troponin I and preliminary results in suspected cases of myocardial infarction. *Clin Chem* 1992;38:2203-14.
13. Adams JE, Bodor GS, Davila-Roman VG, et al. Cardiac troponin I: a marker with high specificity for cardiac injury. *Circulation* 1993;88:101-6.
14. Heidenreich PA, Alloggiamento T, Melsop K, et al. The prognostic value of troponin in patients with non-ST elevation acute coronary syndrome: a meta analysis. *J Am Coll Cardiol* 2001;38:478-85.
15. Morrow DA, Rifai N, Tanasjevic MJ, et al. Clinical efficacy of three assays for cardiac troponin I for risk stratification in acute coronary syndromes - a TIMI 11b substudy. *Clin Chem* 2000;46:453-60.
16. Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation* 2011;124:136-45.
17. Jaffe AS, Apple FS, Morrow A, et al. Being rational about (im)precision: a statement from the biochemistry subcommittee of the Joint European Society of Cardiology/ American College of Cardiology Foundation/ American Heart Association/World Heart Federation Task Force for the Definition of Myocardial Infarction. *Clin Chem* 2010;56:941-3.
18. Tholen DW, Linnet K, Kondratovich M, et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved

Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.

19. McEnroe, RJ, Burnitt, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
20. Tholen DW Kallner A, Kennedy JW, et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline-Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004.

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Korrelationskoeffizient, $r = 0,90$

REFERENZEN

1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. *Circulation* 2012;125:e2-e220.
2. Cardiovascular Diseases, World Health Organization factsheet. <http://www.who.int/mediacentre/factsheets/fs317/en/index.html>
3. World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria of diagnosis of ischemic heart disease. *Circulation* 1979;59:607-9.
4. Thygesen K, Alpert JS, White HD, Joint ESC/ACCF/AHA/WHF Task Force for the Redefinition of Myocardial Infarction. Universal definition of myocardial infarction. *Eur Heart J* 2007;28:2525-38; *Circulation* 2007;116:2634-53; *J Am Coll Cardiol* 2007;50:2173-95.
5. Newby LK, Storrow AB, Gibler B, et al. Bedside multimarker testing for risk stratification in chest pain patients – the chest pain evaluation by creatine kinase-MB, myoglobin and troponin I (CHECKMATE) study. *Circulation* 2001;103:1832-7.
6. Apple FS, Voss E, Lund L, et al. Cardiac troponin, CK-MB and myoglobin for the early detection of acute myocardial infarction and monitoring of reperfusion following thrombolytic therapy. *Clin Chim Acta* 1995;237:59-66.
7. Ohman EM, Casy C, Bengtson JR, et al. Early detection of acute myocardial infarction: Additional diagnostic information from serum concentrations of myoglobin in patients without ST elevations. *Br Heart J* 1990;63:335-8.
8. Tucker JF, Collins RA, Anderson AJ, et al. Value of serial myoglobin levels in the early diagnosis of patients admitted for acute myocardial infarction. *Ann Emerg Med* 1994;24:704-8.
9. Cairns JA, Missirris E, Walker WH. Usefulness of serial determinations of myoglobin and creatine kinase in serum compared for assessment of acute myocardial infarction. *Clin Chem* 1983;29:469-73.
10. Lott JA, Abbott LB. Creatine kinase isoenzymes. *Clin Lab Med* 1986;6:547-76.
11. Mehegan JP, Tobacman LS. Cooperative interaction between troponin molecules bound to the cardiac thin filament. *J Biol Chem* 1991;266:966-72.
12. Bodor GS, Porter S, Landt Y, Ladenson JH. Development of monoclonal antibody for an assay of cardiac troponin I and preliminary results in suspected cases of myocardial infarction. *Clin Chem* 1992;38:2203-14.
13. Adams JE, Bodor GS, Davila-Roman VG, et al. Cardiac troponin I: a marker with high specificity for cardiac injury. *Circulation* 1993;88:101-6.
14. Heidenreich PA, Alloggiamento T, Melsop K, et al. The prognostic value of troponin in patients with non-ST elevation acute coronary syndrome: a meta analysis. *J Am Coll Cardiol* 2001;38:478-85.
15. Morrow DA, Rifai N, Tanasjevic MJ, et al. Clinical efficacy of three assays for cardiac troponin I for risk stratification in acute coronary syndromes - a TIMI 11b substudy. *Clin Chem* 2000;46:453-60.
16. Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation* 2011;124:136-45.
17. Jaffe AS, Apple FS, Morrow A, et al. Being rational about (im)precision: a statement from the biochemistry subcommittee of the Joint European Society of Cardiology/ American College of Cardiology Foundation/ American Heart Association/World Heart Federation Task Force for the Definition of Myocardial Infarction. *Clin Chem* 2010;56:941-3.
18. Tholen DW, Linnet K, Kondratovich M, et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.
19. McEnroe, RJ, Burnett, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
20. Tholen DW, Kallner A, Kennedy JW, et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline-Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004.

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Coefficient de corrélation, $r = 0,90$

RÉFÉRENCES

1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. *Circulation* 2012;125:e2-e220.
2. Cardiovascular Diseases, World Health Organization factsheet. <http://www.who.int/mediacentre/factsheets/fs317/en/index.html>
3. World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria of diagnosis of ischemic heart disease. *Circulation* 1979;59:607-9.
4. Thygesen K, Alpert JS, White HD, Joint ESC/ACCF/AHA/WHF Task Force for the Redefinition of Myocardial Infarction. Universal definition of myocardial infarction. *Eur Heart J* 2007;28:2525-38; *Circulation* 2007;116:2634-53; *J Am Coll Cardiol* 2007;50:2173-95.
5. Newby LK, Storrow AB, Gibler B, et al. Bedside multimarker testing for risk stratification in chest pain patients – the chest pain evaluation by creatine kinase-MB, myoglobin and troponin I (CHECKMATE) study. *Circulation* 2001;103:1832-7.
6. Apple FS, Voss E, Lund L, et al. Cardiac troponin, CK-MB and myoglobin for the early detection of acute myocardial infarction and monitoring of reperfusion following thrombolytic therapy. *Clin Chim Acta* 1995;237:59-66.
7. Ohman EM, Casy C, Bengtson JR, et al. Early detection of acute myocardial infarction: Additional diagnostic information from serum concentrations of myoglobin in patients without ST elevations. *Br Heart J* 1990;63:335-8.
8. Tucker JF, Collins RA, Anderson AJ, et al. Value of serial myoglobin levels in the early diagnosis of patients admitted for acute myocardial infarction. *Ann Emerg Med* 1994;24:704-8.
9. Cairns JA, Missirrli E, Walker WH. Usefulness of serial determinations of myoglobin and creatine kinase in serum compared for assessment of acute myocardial infarction. *Clin Chem* 1983;29:469-73.
10. Lott JA, Abbott LB. Creatine kinase isoenzymes. *Clin Lab Med* 1986;6:547-76.
11. Mehegan JP, Tobacman LS. Cooperative interaction between troponin molecules bound to the cardiac thin filament. *J Biol Chem* 1991;266:966-72.
12. Bodor GS, Porter S, Landt Y, Ladenson JH. Development of monoclonal antibody for an assay of cardiac troponin I and preliminary results in suspected cases of myocardial infarction. *Clin Chem* 1992;38:2203-14.
13. Adams JE, Bodor GS, Davila-Roman VG, et al. Cardiac troponin I: a marker with high specificity for cardiac injury. *Circulation* 1993;88:101-6.
14. Heidenreich PA, Alloggiamento T, Melsop K, et al. The prognostic value of troponin in patients with non-ST elevation acute coronary syndrome: a meta analysis. *J Am Coll Cardiol* 2001;38:478-85.
15. Morrow DA, Rifai N, Tanasjevic MJ, et al. Clinical efficacy of three assays for cardiac troponin I for risk stratification in acute coronary syndromes - a TIMI 11b substudy. *Clin Chem* 2000;46:453-60.
16. Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation* 2011;124:136-45.
17. Jaffe AS, Apple FS, Morrow A, et al. Being rational about (im)precision: a statement from the biochemistry subcommittee of the Joint European Society of Cardiology/ American College of Cardiology Foundation/ American Heart Association/World Heart Federation Task Force for the Definition of Myocardial Infarction. *Clin Chem* 2010;56:941-3.
18. Tholen DW, Linnet K, Kondratovich M, et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.
19. McEnroe, RJ, Burnitt, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
20. Tholen DW, Kallner A, Kennedy JW, et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline-Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004.

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BIBLIOGRAFIA

1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. Circulation 2012;125:e2-e220.
2. Cardiovascular Diseases, World Health Organization factsheet. <http://www.who.int/mediacentre/factsheets/fs317/en/index.html>
3. World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria of diagnosis of ischemic heart disease. Circulation 1979;59:607-9.
4. Thygesen K, Alpert JS, White HD, Joint ESC/ACCF/AHA/WHF Task Force for the Redefinition of Myocardial Infarction. Universal definition of myocardial infarction. Eur Heart J 2007;28:2525-38; Circulation 2007;116:2634-53; J Am Coll Cardiol 2007;50:2173-95.
5. Newby LK, Storrow AB, Gibler B, et al. Bedside multimarker testing for risk stratification in chest pain patients – the chest pain evaluation by

- creatine kinase-MB, myoglobin and troponin I (CHECKMATE) study. *Circulation* 2001;103:1832-7.
- 6. Apple FS, Voss E, Lund L, et al. Cardiac troponin, CK-MB and myoglobin for the early detection of acute myocardial infarction and monitoring of reperfusion following thrombolytic therapy. *Clin Chim Acta* 1995;237:59-66.
 - 7. Ohman EM, Casy C, Bengston JR, et al. Early detection of acute myocardial infarction: Additional diagnostic information from serum concentrations of myoglobin in patients without ST elevations. *Br Heart J* 1990;63:335-8.
 - 8. Tucker JF, Collins RA, Anderson AJ, et al. Value of serial myoglobin levels in the early diagnosis of patients admitted for acute myocardial infarction. *Ann Emerg Med* 1994;24:704-8.
 - 9. Cairns JA, Missirrli E, Walker WH. Usefulness of serial determinations of myoglobin and creatine kinase in serum compared for assessment of acute myocardial infarction. *Clin Chem* 1983;29:469-73.
 - 10. Lott JA, Abbott LB. Creatine kinase isoenzymes. *Clin Lab Med* 1986;6:547-76.
 - 11. Mehegan JP, Tobacman LS. Cooperative interaction between troponin molecules bound to the cardiac thin filament. *J Biol Chem* 1991;266:966-72.
 - 12. Bodor GS, Porter S, Landt Y, Ladenson JH. Development of monoclonal antibody for an assay of cardiac troponin I and preliminary results in suspected cases of myocardial infarction. *Clin Chem* 1992;38:2203-14.
 - 13. Adams JE, Bodor GS, Davila-Roman VG, et al. Cardiac troponin I: a marker with high specificity for cardiac injury. *Circulation* 1993;88:101-6.
 - 14. Heidenreich PA, Alloggiamento T, Melsop K, et al. The prognostic value of troponin in patients with non-ST elevation acute coronary syndrome: a meta analysis. *J Am Coll Cardiol* 2001;38:478-85.
 - 15. Morrow DA, Rifai N, Tanasjevic MJ, et al. Clinical efficacy of three assays for cardiac troponin I for risk stratification in acute coronary syndromes - a TIMI 11b substudy. *Clin Chem* 2000;46:453-60.
 - 16. Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation* 2011;124:136-45.
 - 17. Jaffe AS, Apple FS, Morrow A, et al. Being rational about (im)precision: a statement from the biochemistry subcommittee of the Joint European Society of Cardiology/ American College of Cardiology Foundation/ American Heart Association/World Heart Federation Task Force for the Definition of Myocardial Infarction. *Clin Chem* 2010;56:941-3.
 - 18. Tholen DW, Linnet K, Kondratovich M, et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.
 - 19. McEnroe, RJ, Burnett, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
 - 20. Tholen DW, Kallner A, Kennedy JW, et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline-Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004.

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IB10 sphingotest® 3-in-1 Cardiac (Mio) = 1,00 (VITROS MYOG) – 5,1 ng/ml

Coeficiente de correlación, $r = 0,90$

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1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. *Circulation* 2012;125:e2-e220.
2. Cardiovascular Diseases, World Health Organization factsheet. <http://www.who.int/mediacentre/factsheets/fs317/en/index.html>
3. World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria of diagnosis of ischemic heart disease. *Circulation* 1979;59:607-9.
4. Thygesen K, Alpert JS, White HD, Joint ESC/ACCF/AHA/WHF Task Force for the Redefinition of Myocardial Infarction. Universal definition of myocardial infarction. *Eur Heart J* 2007;28:2525-38; *Circulation* 2007;116:2634-53; *J Am Coll Cardiol* 2007;50:2173-95.
5. Newby LK, Storrow AB, Gibler B, et al. Bedside multimarker testing for risk stratification in chest pain patients – the chest pain evaluation by creatine kinase-MB, myoglobin and troponin I (CHECKMATE) study. *Circulation* 2001;103:1832-7.
6. Apple FS, Voss E, Lund L, et al. Cardiac troponin, CK-MB and myoglobin for the early detection of acute myocardial infarction and monitoring of reperfusion following thrombolytic therapy. *Clin Chim Acta* 1995;237:59-66.
7. Ohman EM, Casy C, Bengston JR, et al. Early detection of acute myocardial infarction: Additional diagnostic information from serum concentrations of myoglobin in patients without ST elevations. *Br Heart J* 1990;63:335-8.
8. Tucker JF, Collins RA, Anderson AJ, et al. Value of serial myoglobin levels in the early diagnosis of patients admitted for acute myocardial infarction. *Ann Emerg Med* 1994;24:704-8.
9. Cairns JA, Missirris E, Walker WH. Usefulness of serial determinations of myoglobin and creatine kinase in serum compared for assessment of acute myocardial infarction. *Clin Chem* 1983;29:469-73.
10. Lott JA, Abbott LB. Creatine kinase isoenzymes. *Clin Lab Med* 1986;6:547-76.
11. Mehegan JP, Tobacman LS. Cooperative interaction between troponin molecules bound to the cardiac thin filament. *J Biol Chem* 1991;266:966-72.
12. Bodor GS, Porter S, Landt Y, Ladenson JH. Development of monoclonal antibody for an assay of cardiac troponin I and preliminary results in suspected cases of myocardial infarction. *Clin Chem* 1992;38:2203-14.
13. Adams JE, Bodor GS, Davila-Roman VG, et al. Cardiac troponin I: a marker with high specificity for cardiac injury. *Circulation* 1993;88:101-6.
14. Heidenreich PA, Alloggiamento T, Melsop K, et al. The prognostic value of troponin in patients with non-ST elevation acute coronary syndrome: a meta analysis. *J Am Coll Cardiol* 2001;38:478-85.
15. Morrow DA, Rifai N, Tanasjevic MJ, et al. Clinical efficacy of three assays for cardiac troponin I for risk stratification in acute coronary syndromes - a TIMI 11b substudy. *Clin Chem* 2000;46:453-60.
16. Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation* 2011;124:136-45.
17. Jaffe AS, Apple FS, Morrow A, et al. Being rational about (im)precision: a statement from the biochemistry subcommittee of the Joint European Society of Cardiology/ American College of Cardiology Foundation/ American Heart Association/World Heart Federation Task Force for the Definition of Myocardial Infarction. *Clin Chem* 2010;56:941-3.
18. Tholen DW, Linnet K, Kondratovich M, et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.
19. McEnroe, RJ, Burnett, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
20. Tholen DW, Kallner A, Kennedy JW, et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline-Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004.

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1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. *Circulation* 2012;125:e2-e220.
2. Cardiovascular Diseases, World Health Organization factsheet. <http://www.who.int/mediacentre/factsheets/fs317/en/index.html>
3. World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria of diagnosis of ischemic heart disease. *Circulation* 1979;59:607-9.
4. Thygesen K, Alpert JS, White HD, Joint ESC/ACCF/AHA/WHF Task Force for the Redefinition of Myocardial Infarction. Universal definition of myocardial infarction. *Eur Heart J* 2007;28:2525-38; *Circulation* 2007;116:2634-53; *J Am Coll Cardiol* 2007;50:2173-95.
5. Newby LK, Storrow AB, Gibler B, et al. Bedside multimarker testing for risk stratification in chest pain patients – the chest pain evaluation by creatine kinase-MB, myoglobin and troponin I (CHECKMATE) study. *Circulation* 2001;103:1832-7.
6. Apple FS, Voss E, Lund L, et al. Cardiac troponin, CK-MB and myoglobin for the early detection of acute myocardial infarction and monitoring of reperfusion following thrombolytic therapy. *Clin Chim Acta* 1995;237:59-66.
7. Ohman EM, Casy C, Bengtson JR, et al. Early detection of acute myocardial infarction: Additional diagnostic information from serum concentrations of myoglobin in patients without ST elevations. *Br Heart J* 1990;63:335-8.
8. Tucker JF, Collins RA, Anderson AJ, et al. Value of serial myoglobin levels in the early diagnosis of patients admitted for acute myocardial infarction. *Ann Emerg Med* 1994;24:704-8.
9. Cairns JA, Missirlis E, Walker WH. Usefulness of serial determinations of myoglobin and creatine kinase in serum compared for assessment of acute myocardial infarction. *Clin Chem* 1983;29:469-73.
10. Lott JA, Abbott LB. Creatine kinase isoenzymes. *Clin Lab Med* 1986;6:547-76.
11. Mehegan JP, Tobacman LS. Cooperative interaction between troponin molecules bound to the cardiac thin filament. *J Biol Chem* 1991;266:966-72.
12. Bodor GS, Porter S, Landt Y, Ladenson JH. Development of monoclonal antibody for an assay of cardiac troponin I and preliminary results in suspected cases of myocardial infarction. *Clin Chem* 1992;38:2203-14.
13. Adams JE, Bodor GS, Davila-Roman VG, et al. Cardiac troponin I: a marker with high specificity for cardiac injury. *Circulation* 1993;88:101-6.
14. Heidenreich PA, Alloggiamento T, Melsop K, et al. The prognostic value of troponin in patients with non-ST elevation acute coronary syndrome: a meta analysis. *J Am Coll Cardiol* 2001;38:478-85.
15. Morrow DA, Rifai N, Tanasjevic MJ, et al. Clinical efficacy of three assays for cardiac troponin I for risk stratification in acute coronary syndromes - a TIMI 11b substudy. *Clin Chem* 2000;46:453-60.
16. Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation* 2011;124:136-45.
17. Jaffe AS, Apple FS, Morrow A, et al. Being rational about (im)precision: a statement from the biochemistry subcommittee of the Joint European Society of Cardiology/ American College of Cardiology Foundation/ American Heart Association/World Heart Federation Task Force for the Definition of Myocardial Infarction. *Clin Chem* 2010;56:941-3.

18. Tholen DW, Linnet K, Kondratovich M, et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.
19. McEnroe, RJ, Burnitt, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
20. Tholen DW Kallner A, Kennedy JW, et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline-Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004.