English

IB10 sphingotest® NT-proBNP

For the Quantitative Determination of NT-proBNP in
EDTA and Lithium Heparin Human Whole Blood and Plasma

Explanation of Symbols

	CE Mark of Conformity
	Manufacturer
	Catalog Number
	Expiry date/Use by
	Lot number
	<i>In Vitro</i> Diagnostic medical device
	Consult instructions for use
	Store between 2 °C and 8 °C
	European Authorized Representative
	Contains sufficient for <n> tests
	Do not reuse
	Serial Number
	Insert disc label side up

Nexus Dx, Inc.,
6759 Mesa Ridge Road
San Diego, California 92121 USA
Telephone: 1 (858) 410 4600
Fax: 1 (858) 410 4700

TheraGenesis GmbH
Bahnhofstrasse 5
55276 Oppenheim, Germany
Tel: +49 (0) 151 506 403 14

022-00064 REV B

IB10 sphingotest® NT-proBNP

For *In Vitro* Diagnostic Use

INTENDED USE

The IB10 sphingotest® NT-proBNP is a rapid test for the *in vitro* quantitative determination of N-terminal pro-brain natriuretic peptide (NT-proBNP) in human whole blood or plasma using EDTA or lithium heparin as the anticoagulant. The NT-proBNP Test is intended for use in conjunction with the Nexus IB10 Analyzer and provides quantitative results in 20 minutes.

This Test is designed for professional use only and may be used in hospital central laboratories and in alternate care settings such as emergency departments, critical care units, and other sites where near patient testing is practiced.

NT-proBNP determinations aid in the diagnosis and management of Congestive Heart Failure (CHF).

SUMMARY AND EXPLANATION OF THE TEST

Congestive Heart Failure (CHF) is a debilitating disease currently afflicting more than 6 million people in the United States with approximately 280,000 associated deaths each year.¹ In the United States, it is the most rapidly growing cardiovascular disease and it has been estimated that approximately 20 million more Americans may have asymptomatic cardiac impairment.²

The class of cardiac neurohormones was first described by de Bold *et al.*^{3,4} This family of structurally similar but genetically distinct molecules includes BNP, atrial natriuretic peptide (ANP) and C-type natriuretic peptide (CNP). These three natriuretic peptides (NPs) are synthesized as high molecular weight precursors. NT-proBNP is a cleavage product of pro-brain natriuretic peptide (proBNP). The precursor molecule proBNP is a peptide consisting of 108 amino acids.^{5,6} During intracellular peptide maturation, cleavage of this precursor molecule by an endopeptidase results in formation of the biologically active BNP peptide and the biologically inactive 76 amino acid N-terminal fragment, NT-proBNP.⁷

The NPs possess potent diuretic, natriuretic and vasodilatory properties and have been reported as valuable diagnostic and prognostic markers in cardiovascular disease, particularly for patients in New York Heart Association (NYHA) classes I-IV CHF.^{8,9} In particular, measurement of plasma concentrations of NT-proBNP has utility as a valuable tool for aiding in the diagnosis and the assessment of severity of patients with CHF.^{10,11}

PRINCIPLE

The Nexus Dx immunochemistry system combines chemistry with microfluidics and centrifugal flow to rapidly prepare a cell free plasma from whole blood that can then be moved through a channel to rehydrate, solubilize and mix with freeze dried immunoconjugates. Using a combination of active flow and capillary action, the test sample is quantitatively measured in 20 minutes with an optical signal level proportional to the analyte(s) concentration.

After addition of the patient sample, the entire test is performed within the Nexus IB10 Analyzer which provides control of the temperature of the disc, as well as the sequence, centrifugal flow, mixing, incubation time, final signal measurement, quantitation and reporting of results. The Test disc includes a positive internal control to ensure that the test has operated properly. Each lot is calibrated to ensure that lot-to-lot variability is minimized. Lot specific calibration along with additional information such as the lot expiration date is contained on a QR code label affixed to each disc. It is recommended that external controls be tested at appropriate time intervals to confirm that the system and test lot are performing within acceptable limits.

REAGENTS

The IB10 sphingotest® NT-proBNP contains all required reagents to evaluate the level of NT-proBNP including biotin conjugated polyclonal anti-NT- proBNP antibody, dye conjugated monoclonal anti-NT-proBNP antibody and streptavidin immobilized at the detection area on the disc.

MATERIALS PROVIDED

Each box contains the following:

- IB10 sphingotest® NT-proBNP discs, each individually sealed in a foil pouch with a desiccant.
- Instructions for Use (IFU).

MATERIALS/EQUIPMENT REQUIRED BUT NOT PROVIDED

1. Nexus IB10 Analyzer – Model #BCA-IB10.
2. Commercially available NT-proBNP Controls for external Quality Control (QC). **Contact the Distributor in your area** for recommended external QC materials or related technical assistance.
3. Calibrated reusable fixed or variable volume pipette gun with high precision and accuracy capable of delivering 500 µL of whole blood or plasma.
4. Disposable pipette tips capable of accepting and delivering 500 µL of whole blood or plasma.

PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Carefully follow the Instructions for Use.
- Prior to testing controls or patient samples, ensure that the Analyzer software is updated to the latest version. Refer to the Nexus IB10 manual for specific instructions.
- Wear disposable gloves while handling samples.
- Handle samples with care. Samples and used Test discs should be treated as if potentially infectious and should be discarded as biohazardous materials according to local regulations.
- Thoroughly wash hands following handling.
- The result obtained from the IB10 sphingotest® NT-proBNP does not provide a definitive diagnosis and should be interpreted by the physician in conjunction with other appropriate test results and patient clinical findings.
- Keep the Test in the sealed pouch until ready for use.
- Do not use the Test if the pouch is damaged or the seal is broken.
- Do not use the Test after the expiration date printed on the pouch.
- Prior to use, place the unopened pouch at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
- Always pay attention to cleanliness when handling the Test. Avoid any contamination from fingerprints or foreign substances. Do not contaminate the sample channel inlet.
- Do not drop or damage the disc.
- The Test disc should be inserted with the label side up into the Nexus IB10 Analyzer tray immediately after injecting the sample into the disc.
- Do not turn the disc upside down.
- This is a quantitative test; therefore no visual interpretation of the result should be made.

STORAGE AND STABILITY

- Store the IB10 sphingotest® NT-proBNP disc between 2 and 8 °C (35 to 46 °F) until the expiration date printed on the pouch is reached.
- The IB10 sphingotest® NT-proBNP in its sealed pouch is stable at 18 to 30 °C (64 to 86 °F) for 30 days, provided the expiration date printed on the pouch is not exceeded.

SAMPLE COLLECTION AND PREPARATION

- The IB10 sphingotest® NT-proBNP is to be run using EDTA or lithium heparin whole blood or plasma samples.
- As cardiac proteins are relatively unstable, it is recommended that sample be tested as soon as possible; however, whole blood must be tested within 24 hours of collection.
- Whole blood samples can be transported and stored up to 24 hours at ambient temperature.
- Plasma samples should be kept frozen at -20 °C (-4 °F) or lower if longer storage is required.
- Allow samples to equilibrate to room temperature (19 to 25 °C/66 to 77 °F) prior to testing.
- For serial testing of the same patient, the same sample type (EDTA or lithium heparin whole blood or plasma) should be used.

PROCEDURE

Nexus IB10 Analyzer

 Consult the Nexus IB10 Analyzer User Manual

For Analyzer installation, start up and complete instructions for use refer to the **Nexus IB10 Analyzer User Manual**. Operator must consult the User Manual prior to use to become familiar with the proper operation and quality control procedures.

PERFORMING SYSTEM CHECK AND DISC CALIBRATION

Each time the Nexus IB10 Analyzer is turned on, a Self-Check is automatically performed. The QR code on each test disc contains information for disc calibration which the Analyzer automatically reads when running a test.

RUNNING QC WITH EXTERNAL CONTROLS

The manufacturer recommends the use of commercially available NT-proBNP Controls (please refer to section **Materials/Equipment Required But Not Provided**). Ensure that the NT-proBNP controls are handled and prepared according to the corresponding Instructions for Use (IFU).

1. Remove an unopened pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes prior to testing.
 2. Open the pouch and remove the Test disc.
 3. Place the Test disc on a level surface.
 4. On the Nexus IB10 Analyzer press **New Analysis**.
 5. The Analyzer will perform a general system check.
 6. Enter the external control ID manually (up to 20 characters can be used for the ID) or enter the external control ID by using the barcode scanner.
 7. Mix the external quality control vial by gently inverting the vial several times before sampling.
- 8. Testing External Control Sample on the IB10 sphingotest® NT-proBNP.**
- Using a precision pipette (fixed or adjusted to 500 µL) slowly draw well-mixed external quality control sample into the pipette tip.
 - Positioning the tapered pipette tip at a 45° angle, pierce the X on the red dot to expose the sample channel inlet.
 - Slowly express the external control sample into the inlet applying minimum, but continuous force, on the pipette plunger.
 - Express sample to the first stop on the pipette, at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
 - Press **QC** on the Nexus IB10 Analyzer display.
 - When the tray opens, insert the filled Test disc into the tray and press **Run**.
 - The tray will close and perform a disc validity check.
 - A screen appears for selecting quality control materials (please refer to the **Quality Control Settings** section of the **Nexus IB10 Analyzer User Manual** for how to update quality control material [external controls]).
 - Select the quality control material to be tested.
 - Press **OK** on the Nexus IB10 Analyzer display.
 - In 20 minutes, the Nexus IB10 Analyzer will display the results on the screen.
 - Results will print automatically (if selected during Set Up) or press **Print**.
 - When the test is complete, analyze and compare the result with the Expected Value reported in the external control IFU for the external control level as measured using the IB10 sphingotest® NT-proBNP.
 - Remove the Test disc and discard in appropriate receptacle.
 - If the external control result is outside the expected range, refer to the Quality Control section below.

Note: If the test run is cancelled before a Test result is displayed, the Test disc cannot be reused and should be disposed of appropriately.

TESTING PATIENT SAMPLES ON THE NEXUS IB10 ANALYZER

1. Remove an unopened Test pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
 2. Open the pouch and remove the Test disc.
 3. Place the Test disc on a level surface.
 4. On the Nexus IB10 Analyzer press **New Analysis**.
 5. The Analyzer will perform a general system check.
 6. Enter Patient ID manually (up to 20 characters can be used for the ID) or enter the Patient ID by using the barcode scanner.
 7. Mix the whole blood patient sample tube by gently inverting the tube several times before testing.
- 8. Testing Patient Sample on the IB10 sphingotest® NT-proBNP.**
- Using a precision pipette (fixed or adjusted to 500 µL), slowly draw the well-mixed patient sample into the pipette tip.
 - Positioning the tapered pipette tip at a 45° angle, pierce the X on the red dot to expose the sample channel inlet.
 - Slowly express the patient sample into the inlet applying minimum, but continuous force, on the pipette plunger.
 - Express sample to the first stop on the pipette, at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
 - Press **OK** on the Nexus IB10 Analyzer display.
 - When the tray opens, insert the filled Test disc into the tray and press **Run**.
 - In 20 minutes, the Nexus IB10 Analyzer will display the results on the screen.
 - Results will print automatically (if selected during Set Up) or press **Print**.
 - Remove the Test disc and discard in appropriate receptacle.

Note: If the test run is cancelled before a test result is displayed, the Test disc cannot be reused and should be disposed of appropriately.

INTERPRETATION OF RESULTS

The range of NT-proBNP concentration reported by the Nexus IB10 Analyzer is 30 pg/mL to 5,000 pg/mL. Results below or above this range will be shown as “< 30 pg/mL” or “> 5,000 pg/mL,” respectively.

- **Recommended Decision Threshold Values:**
 - Patients under 75 years of age: 125 pg/mL
 - Patients 75 years of age and older: 450 pg/mL
- NT-proBNP results less than or equal to the decision threshold values are considered normal values, representative of patients without CHF.
- Results greater than the recommended decision threshold values are considered abnormal and are suggestive of patients with CHF.
- NT-proBNP results greater than 5,000 pg/mL are considered very high values for NT-proBNP and exceed the upper limits of the IB10 sphingotest® NT-proBNP.

Quality Control

EXTERNAL CONTROLS

Good laboratory practice includes the use of external controls to ensure proper test performance. It is recommended that prior to using a new lot or shipment of IB10 sphingotest® NT-proBNP, the performance of the lot should be confirmed by testing with external controls (see section **Materials/Equipment Required But Not Provided**) to ensure the test will deliver the correct result. The frequency of quality control testing should be determined according to individual laboratory standard quality control procedures. Upon confirmation of the expected results, the test discs are ready for use with patient samples. Controls should also be used any time the validity of test results is questionable. If external controls do not perform as expected, do not use the IB10 sphingotest® NT-proBNP and contact the Distributor in your area for Technical Assistance.

INTERNAL CONTROL

The IB10 sphingotest® NT-proBNP has a built in positive procedural control. The Nexus IB10 Analyzer automatically determines the presence of this control thereby confirming that the test run has delivered a valid result. If the control does not form or if it is not recognized by the Analyzer, the test result is considered “invalid” and the test must be repeated.

LIMITATIONS

The results of the IB10 sphingotest® NT-proBNP should be used in conjunction with other available laboratory and clinical information. Evaluations of lithium heparin plasma samples showed a mean 13% increase in observed values compared to EDTA plasma samples. A comparison of 48 paired EDTA and lithium heparin plasma samples yielded a Passing-Bablok regression line with the following equation:

$$\text{Lithium heparin plasma} = 1.16 \text{ (EDTA plasma)} - 2.38 \text{ pg/mL}$$

Correlation coefficient, $r = 0.96$

It is not recommended that lithium heparin and EDTA anticoagulated samples from the same patient be used interchangeably with this Test. Reference ranges should be established for the evaluation of lithium heparin samples and for EDTA samples.

Other substances and/or factors not listed, e.g. technical or procedural errors, may interfere with the test and cause inaccurate results.

Nexus Dx, Inc. offer products for their intended use. Refer to the specific product literature for the intended use statements for each product. Product claims are subject to change. Nexus Dx, Inc.'s expressed and implied warranties (inclusive of implied warranties of merchantability and fitness) are conditional upon adherence to, or observance of Nexus Dx, Inc. published directions with regard to the use of Nexus Dx, Inc. products. Under no circumstances will Nexus Dx, Inc. be liable for any indirect or consequential damages.

For Technical Assistance please contact the Distributor in your area.

The IB10 sphingotest® NT-proBNP is manufactured under license from Roche Diagnostics GmbH.

Performance Characteristics

MEASURING RANGE

The IB10 sphingotest® NT-proBNP has been demonstrated to provide measurable results at NT-proBNP levels from 30 pg/mL to 5,000 pg/mL.

ANALYTICAL SENSITIVITY

The LoD (Limit of Detection) of the IB10 sphingotest® NT-proBNP is 30 pg/mL, determined according to the Clinical and Laboratory Standards Institute (CLSI, formerly NCCLS) guideline EP17-A2.¹² The proportion of false positives (α) and false negatives (β) are less than 5% (Limit of Blank = 18 pg/mL). LoQ (Limit of Quantitation) is the lowest NT-proBNP concentration that can be reproducibly measured with a total coefficient of variation of at most 15%. It was determined to be 50 pg/mL.

CROSS-REACTIVITY AND INTERFERING SUBSTANCES

Drugs

The following drugs were tested for potential interference in the IB10 sphingotest® NT-proBNP (Table 1). The list of drugs encompasses common prescription and over-the-counter medications, as well as those often prescribed in a cardiac patient population. The drugs were tested at concentrations as recommended in the CLSI Approved Guideline EP7-A2 'Interference Testing in Clinical Chemistry,'¹³ or at least three times the highest concentration reported following a therapeutic dose. No significant interference with the IB10 sphingotest® NT-proBNP was observed for the drugs listed in the table below.

Table 1.

Drug	Drug	Drug
Acetaminophen	Caffeine	Methyl-DOPA
Acetylsalicylic acid (aspirin)	Captopril	Nifedipine
Allopurinol	Digoxin	Phenytoin
Ampicillin	Dopamine	Theophylline
Ascorbic acid (vitamin C)	Erythromycin	Verapamil
Atenolol	Furosemide	

PROTEINS, PEPTIDES AND ENDOGENOUS SUBSTANCES

The following proteins, peptides and endogenous substances were tested for potential cross-reactivity and interference in the IB10 sphingotest® NT-proBNP at the maximum concentration of substance indicated (Table 2). No substance demonstrated significant cross-reactivity or interference using the method provided in CLSI EP7-A2.

Table 2.

Substance	Maximum Concentration
Proteins and Peptides	
BNP-32	3.5 µg/mL
Atrial Natriuretic Polypeptide (ANP) (1-28)	3.1 µg/mL
NT-ProANP (1-30)	3.5 µg/mL
NT-Pro-ANP (31-67)	1.0 µg/mL
NT-Pro-ANP (79-98)	1.0 µg/mL
C-type Natriuretic Peptide-53 (CNP-53)	2.2 µg/mL
Endothelin I	20 pg/mL
Endogenous substances	
Albumin	5 g/dL
Hemoglobin	1.4 g/dL
Rheumatoid Factor (RF)	1500 IU/mL
Creatinine	2 mg/dL
Bilirubin	60 mg/dL
Triglycerides	5 g/dL
Urea	18 mg/dL
Cholesterol	250 mg/dL
Biotin	400 ng/mL

HOOK EFFECT

No high dose Hook Effect was observed for NT-proBNP concentrations up to 300,000 pg/mL.

PRECISION

The precision of IB10 sphingotest® NT-proBNP was determined using samples where recombinant NT-proBNP was added to normal human plasma at three concentrations (Table 3). The within-run and total precision was performed in two runs per day, in replicates of 5 per run at each concentration level for a 20 day period for a total number of repetitions of 200 at each concentration level. The within-run and total coefficients of variation (CVs) were computed according to CLSI guideline EP05-A3.¹⁴

Table 3.

Sample	Mean level (pg/mL)	Within-Run Precision		Total Precision	
		Std. dev. (pg/mL)	CV (%)	Std. dev. (pg/mL)	CV (%)
1	130.07	16.65	12.8	17.95	13.8
2	434.19	51.23	11.8	53.84	12.4
3	1615.8	184.20	11.4	187.43	11.6

WHOLE BLOOD Vs. PLASMA COMPARISON

A comparison study was performed using whole blood and plasma samples. When performing a Passing-Bablok regression analysis comparing the whole blood concentrations versus the corresponding plasma concentrations from the same subject samples, a slope of 1.01 (95% C.I. = [0.97-1.06]), an intercept of 2.72 pg/mL, and a correlation coefficient of $r = 0.99$ were obtained.

Expected Values

RECOMMENDED DECISION THRESHOLD VALUES

From calibration based on the reference Roche Elecsys® proBNP II assay, the recommended Decision Threshold Values for IB10 sphingotest® NT-proBNP are as follows:

- Patients under 75 years of age: 125 pg/mL
- Patients 75 years of age and older: 450 pg/mL

Each laboratory should establish a reference range that represents the patient population that is to be evaluated.

METHOD COMPARISON

An equivalence study was performed between the IB10 sphingotest® NT-proBNP using the Nexus IB10 Analyzer and the Roche Elecsys® proBNP II assay. A total of 321 samples were tested within the range of 30–5000 pg/mL. The Passing-Bablok regression was:

$$\text{IB10 sphingotest® NT-proBNP} = 0.97 \text{ (Roche Elecsys® proBNP II assay)} - 11.81 \text{ pg/mL}$$

Correlation coefficient, $r = 0.96$

REFERENCES

1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. Circulation 2012;125: e2-e220.
2. Consensus Recommendations for the Management of Chronic Heart Failure. Am J Cardiol 1999; 83(2A):1A-38A.
3. de Bold AJ. Heart atria granularity effects of changes in water-electrolyte balance. Proc Soc Exp Biol Med 1979;161:508-11.
4. de Bold AJ, Borenstein HB, Veress AT, et al. A rapid and potent natriuretic response to intravenous injection of atrial myocardial extracts in rats. Life Sci 1981;28:89-94.
5. Hall, C. Essential biochemistry and physiology of (NT-pro) BNP. Eur J Heart Fail 2004;6:257-60.
6. Gotze JP, Kastrup J. Plasma pro-brain natriuretic peptides are strong biochemical markers in clinical cardiology. Scan J Clin Lab Invest 2001;61, Suppl 234:47-51.
7. The Criteria Committee of the New York Heart Association. Nomenclature and Criteria for Diagnosis of Diseases of the Heart and Great Vessels. 9th ed. Boston, Mass: Little, Brown and Co; 1994:253-256.
8. Costello-Boerrigter LC, Burnett JC. The prognostic value of N-terminal proB-type natriuretic peptide. Nat Clin Pract Cardiovasc Med 2005;2:194-201.

9. Cowie MR, Jourdain P, Maisel A, et al. Clinical applications of B-type natriuretic peptide (BNP) testing. *Eur Heart J* 2003;24:1710-18.
10. Mair J, Hammerer-Lercher A, Puschendorf B. The Impact of cardiac natriuretic peptide determination on the diagnosis and management of heart failure. *Clin Chem Lab Med* 2001;39:571-88.
11. McDonagh TA, Holmer S, Raymond I, et al. NT-proBNP and the diagnosis of heart failure: a pooled analysis of three European epidemiological studies. *Eur J Heart Fail* 2004;6:269-73.
12. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
13. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
14. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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Korrelationskoeffizient, $r = 0,96$

REFERENZEN

1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. Circulation 2012;125: e2-e220.
2. Consensus Recommendations for the Management of Chronic Heart Failure. Am J Cardiol 1999; 83(2A):1A-38A.
3. de Bold AJ. Heart atria granularity effects of changes in water-electrolyte balance. Proc Soc Exp Biol Med 1979;161:508-11.
4. de Bold AJ, Borenstein HB, Veress AT, et al. A rapid and potent natriuretic response to intravenous injection of atrial myocardial extracts in rats. Life Sci 1981;28:89-94.
5. Hall, C. Essential biochemistry and physiology of (NT-pro) BNP. Eur J Heart Fail 2004;6:257-60.
6. Gotze JP, Kastrup J. Plasma pro-brain natriuretic peptides are strong biochemical markers in clinical cardiology. Scan J Clin Lab Invest 2001;61, Suppl 234:47-51.
7. The Criteria Committee of the New York Heart Association. Nomenclature and Criteria for Diagnosis of Diseases of the Heart and Great Vessels. 9th ed. Boston, Mass: Little, Brown and Co; 1994:253-256.
8. Costello-Boerrigter LC, Burnett JC. The prognostic value of N-terminal proB-type natriuretic peptide. Nat Clin Pract Cardiovasc Med 2005;2:194-201.
9. Cowie MR, Jourdain P, Maisel A, et al. Clinical applications of B-type natriuretic peptide (BNP) testing. Eur Heart J 2003;24:1710-18.
10. Mair J, Hammerer-Lercher A, Puschendorf B. The Impact of cardiac natriuretic peptide determination on the diagnosis and management of heart failure. Clin Chem Lab Med 2001;39:571-88.
11. McDonagh TA, Holmer S, Raymond I, et al. NT-proBNP and the diagnosis of heart failure: a pooled analysis of three European epidemiological studies. Eur J Heart Fail 2004;6:269-73.
12. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
13. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
14. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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Coefficient de corrélation, $r = 0,96$

RÉFÉRENCES

1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. Circulation 2012;125: e2-e220.
2. Consensus Recommendations for the Management of Chronic Heart Failure. Am J Cardiol 1999; 83(2A):1A-38A.
3. de Bold AJ. Heart atria granularity effects of changes in water-electrolyte balance. Proc Soc Exp Biol Med 1979;161:508-11.
4. de Bold AJ, Borenstein HB, Veress AT, et al. A rapid and potent natriuretic response to intravenous injection of atrial myocardial extracts in rats. Life Sci 1981;28:89-94.
5. Hall, C. Essential biochemistry and physiology of (NT-pro) BNP. Eur J Heart Fail 2004;6:257-60.
6. Gotze JP, Kastrup J. Plasma pro-brain natriuretic peptides are strong biochemical markers in clinical cardiology. Scan J Clin Lab Invest 2001;61, Suppl 234:47-51.
7. The Criteria Committee of the New York Heart Association. Nomenclature and Criteria for Diagnosis of Diseases of the Heart and Great Vessels. 9th ed. Boston, Mass: Little, Brown and Co; 1994:253-256.
8. Costello-Boerigter LC, Burnett JC. The prognostic value of N-terminal proB-type natriuretic peptide. Nat Clin Pract Cardiovasc Med 2005;2:194-201.
9. Cowie MR, Jourdain P, Maisel A, et al. Clinical applications of B-type natriuretic peptide (BNP) testing. Eur Heart J 2003;24:1710-18.
10. Mair J, Hammerer-Lercher A, Puschendorf B. The Impact of cardiac natriuretic peptide determination on the diagnosis and management of heart failure. Clin Chem Lab Med 2001;39:571-88.
11. McDonagh TA, Holmer S, Raymond I, et al. NT-proBNP and the diagnosis of heart failure: a pooled analysis of three European epidemiological studies. Eur J Heart Fail 2004;6:269-73.
12. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
13. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
14. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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RIFERIMENTI

1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. *Circulation* 2012;125: e2-e220.
2. Consensus Recommendations for the Management of Chronic Heart Failure. *Am J Cardiol* 1999; 83(2A):1A-38A.
3. de Bold AJ. Heart atria granularity effects of changes in water-electrolyte balance. *Proc Soc Exp Biol Med* 1979;161:508-11.
4. de Bold AJ, Borenstein HB, Veress AT, et al. A rapid and potent natriuretic response to intravenous injection of atrial myocardial extracts in rats. *Life Sci* 1981;28:89-94.
5. Hall, C. Essential biochemistry and physiology of (NT-pro) BNP. *Eur J Heart Fail* 2004;6:257-60.
6. Gotze JP, Kastrup J. Plasma pro-brain natriuretic peptides are strong biochemical markers in clinical cardiology. *Scan J Clin Lab Invest* 2001;61, Suppl 234:47-51.
7. The Criteria Committee of the New York Heart Association. Nomenclature and Criteria for Diagnosis of Diseases of the Heart and Great Vessels. 9th ed. Boston, Mass: Little, Brown and Co; 1994:253-256.
8. Costello-Boerrigter LC, Burnett JC. The prognostic value of N-terminal proB-type natriuretic peptide. *Nat Clin Pract Cardiovasc Med* 2005;2:194-201.
9. Cowie MR, Jourdain P, Maisel A, et al. Clinical applications of B-type natriuretic peptide (BNP) testing. *Eur Heart J* 2003;24:1710-18.
10. Mair J, Hammerer-Lercher A, Puschendorf B. The Impact of cardiac natriuretic peptide determination on the diagnosis and management of heart failure. *Clin Chem Lab Med* 2001;39:571-88.
11. McDonagh TA, Holmer S, Raymond I, et al. NT-proBNP and the diagnosis of heart failure: a pooled analysis of three European epidemiological studies. *Eur J Heart Fail* 2004;6:269-73.
12. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
13. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
14. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. Circulation 2012;125: e2-e220.
2. Consensus Recommendations for the Management of Chronic Heart Failure. Am J Cardiol 1999; 83(2A):1A-38A.
3. de Bold AJ. Heart atria granularity effects of changes in water-electrolyte balance. Proc Soc Exp Biol Med 1979;161:508-11.
4. de Bold AJ, Borenstein HB, Veress AT, et al. A rapid and potent natriuretic response to intravenous injection of atrial myocardial extracts in rats. Life Sci 1981;28:89-94.

5. Hall, C. Essential biochemistry and physiology of (NT-pro) BNP. *Eur J Heart Fail* 2004;6:257-60.
6. Gotze JP, Kastrup J. Plasma pro-brain natriuretic peptides are strong biochemical markers in clinical cardiology. *Scan J Clin Lab Invest* 2001;61, Suppl 234:47-51.
7. The Criteria Committee of the New York Heart Association. Nomenclature and Criteria for Diagnosis of Diseases of the Heart and Great Vessels. 9th ed. Boston, Mass: Little, Brown and Co; 1994:253-256.
8. Costello-Boerrigter LC, Burnett JC. The prognostic value of N-terminal proB-type natriuretic peptide. *Nat Clin Pract Cardiovasc Med* 2005;2:194-201.
9. Cowie MR, Jourdain P, Maisel A, et al. Clinical applications of B-type natriuretic peptide (BNP) testing. *Eur Heart J* 2003;24:1710-18.
10. Mair J, Hammerer-Lercher A, Puschendorf B. The Impact of cardiac natriuretic peptide determination on the diagnosis and management of heart failure. *Clin Chem Lab Med* 2001;39:571-88.
11. McDonagh TA, Holmer S, Raymond I, et al. NT-proBNP and the diagnosis of heart failure: a pooled analysis of three European epidemiological studies. *Eur J Heart Fail* 2004;6:269-73.
12. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
13. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
14. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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