English

IB10 sphingotest® Troponin-99

For the Quantitative Determination of Troponin-99 in
Lithium Heparin Human Whole Blood and Plasma

Explanation of Symbols

	CE Mark of Conformity
	Manufacturer
	Catalog Number
	Expiry date/Use by
	Lot number
	<i>In Vitro</i> Diagnostic medical device
	Consult instructions for use
	Store between 2 °C and 8 °C
	European Authorized Representative
	Contains sufficient for <n> tests
	Do not reuse
	Serial Number
	Insert disc label side up

Nexus Dx, Inc.,
6759 Mesa Ridge Road
San Diego, California 92121 USA
Telephone: 1 (858) 410 4600
Fax: 1 (858) 410 4700

TheraGenesis GmbH
Bahnhofstrasse 5
55276 Oppenheim, Germany
Tel: +49 (0) 151 506 403 14

022-00067 REV B

IB10 sphingotest® Troponin-99

For *In Vitro* Diagnostic Use

INTENDED USE

The IB10 sphingotest® Troponin-99 is a rapid point-of-care (POC) immunoassay for the *in vitro* quantitative determination of Cardiac Troponin-99 (cTnI) in lithium heparin whole blood and plasma. The Troponin-99 Test is intended for use in conjunction with the Nexus IB10 Analyzer and provides quantitative results in 20 minutes.

This Test is designed for professional use only and may be used in hospital central laboratories and in alternate care settings such as emergency departments, critical care units, and other sites where near patient testing is practiced.

cTnI determinations aid in the diagnosis of myocardial infarction (MI) and patients with non-ST-segment elevation (NSTEMI) acute coronary syndromes (ACS). Elevated cTnI levels in the NSTEMI ACS subpopulation correlate with the relative risk of mortality, MI or increased probability of ischemic events, requiring urgent revascularization procedures.

SUMMARY AND EXPLANATION OF THE TEST

Acute coronary syndrome (ACS) including Coronary Heart Disease (CHD) and Angina Pectoris (AP) is the leading cause of morbidity and mortality among both men and women, affecting more than 16.0 million people in the United States. The presentation of ACS is varied, with acute myocardial infarction (AMI) being the most dramatic presentation. Symptoms of ACS may develop suddenly. In 2008, there were over 800,000 new AMIs in the US.¹ On an international basis, the World Health Organization (WHO) reports over 7.2 million deaths a year from CHD. Globally it is estimated that 17.3 million deaths were due to Cardiovascular Disease in 2008.² This includes death due to both CHD and stroke.

For over 30 years, the WHO has recommended that the diagnosis of AMI be guided by positive findings on at least two of the following three criteria: (1) patient history/physical examination; (2) electrocardiogram; (3) changes in blood levels of cardiac protein markers.³ Previously the primary marker of choice was creatine kinase (MB isoenzyme), but research over the past 20 years has confirmed that cardiac striated muscle proteins of the troponin complex (cTnI or cTnT) are far more accurate and sensitive in differential diagnosis. Patient history and physical examination are critical, but often provide insufficient information to differentiate AMI from other cardiac abnormalities. In the absence of well defined electrocardiographic ST-segment elevations, troponin measurements are essential for the accurate diagnosis of ACS.⁴ The utility of determining the serum levels of the cardiac isoforms of TnI has been validated in numerous multi-center trials and cTnI is now considered the biomarker of choice in the initial diagnosis of AMI.⁵

The temporal release of cTnI into the serum has been investigated and compared to those of the other established cardiac markers, such as creatine kinase MB (CK-MB), and myoglobin.⁶ Following myocyte necrosis, cTnI is released into the circulation with levels exceeding the upper reference limit of normal within 4-6 hours, and peak levels are reached after 12-24 hours.⁷ This early release profile is similar to CK-MB. However, CK-MB levels return to normal values in about 72 hours, while levels of cTnI remain elevated for 5-7 days.

cTnI IN ACUTE CORONARY SYNDROME

Since the introduction of the first commercial assays for Troponin-99 in 1996, numerous clinical studies have confirmed that cTnI is the biomarker of choice in identifying cardiac myocyte ischemia and necrosis.⁸ Even low levels of cTnI just above the 99th percentile of levels obtained from a normal reference population indicate increased risk. To assure that the test results are interpreted with as much accuracy as possible, consensus groups now unanimously recommend that chest pain patients with low level troponin values on presentation to medical personnel, especially in the absence of confirmatory ECG evidence, be tested 2 additional times during the subsequent 12-24 hour period to determine whether cTnI is indeed rising during the first 24 hours.^{4,9-11}

PRINCIPLE

The Nexus Dx immunochemistry system combines chemistry with microfluidics and centrifugal flow to rapidly prepare a cell free plasma from whole blood that can then be moved through a channel to rehydrate, solubilize and mix with freeze dried immunoconjugates. Using a combination of active flow and capillary action, the test sample is quantitatively measured in 20 minutes with an optical signal level proportional to the analyte(s) concentration.

After addition of the patient sample, the entire test is performed within the Nexus IB10 Analyzer which provides control of the temperature of the disc, as well as the sequence, centrifugal flow, mixing, incubation time, final signal measurement, quantitation and reporting of results. The Test disc includes a positive internal control to ensure that the Test disc has operated properly. Each lot is calibrated to ensure that lot-to-lot variability is minimized. Lot specific calibration along with additional information such as the lot expiration date is contained on a QR code label affixed to each disc. It is recommended that external controls be tested at appropriate time intervals to confirm that the reagent lot is performing within acceptable limits.

REAGENTS

The IB10 sphingotest® Troponin-99 contains all required reagents to evaluate the level of cTnI including dye conjugated monoclonal anti-cTnI antibody, biotin conjugated polyclonal anti-cTnI antibody and streptavidin immobilized at the detection area on the disc.

MATERIALS PROVIDED

Each box contains the following:

- 10 IB10 sphingotest® Troponin-99 discs, each individually sealed in a foil pouch with a desiccant.
- 1 Instructions for Use (IFU).

MATERIALS/EQUIPMENT REQUIRED BUT NOT PROVIDED

1. Nexus IB10 Analyzer - Model #BCA-IB10
2. Commercially available Troponin-99 Controls for external Quality Control (QC). **Contact the Distributor in your area** for recommended external QC materials or related technical assistance.
3. Calibrated reusable fixed or variable volume pipette gun with high precision and accuracy capable of delivering 500 µL of whole blood or plasma.
4. Disposable pipette tips capable of accepting and delivering 500 µL of whole blood or plasma.

PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Carefully follow the Instructions for Use.
- Prior to testing controls or patient samples, ensure that the Analyzer software is updated to the latest version. Refer to the Nexus IB10 manual for specific instructions.
- Wear disposable gloves while handling samples.
- Handle samples with care. Samples and used Test discs should be treated as potentially infectious and should be discarded as biohazardous material according to local regulations.
- Thoroughly wash hands following handling.
- The result obtained from the IB10 sphingotest® Troponin-99 does not provide a definitive diagnosis and should be interpreted by a physician in conjunction with other laboratory test results according to current medical Guidelines and patient clinical findings.
- Keep the Test in the sealed pouch until ready for use.
- Do not use the Test if the pouch is damaged or the seal is broken.
- Do not use the Test after the expiration date printed on the pouch.
- Prior to use, place the unopened pouch at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
- Always pay attention to cleanliness when handling the Test. Avoid any contamination from fingerprints or foreign substances. Do not contaminate the sample channel inlet.
- Do not drop or damage the Test disc.
- The Test disc should be inserted with the label side up, into the Nexus IB10 Analyzer tray immediately after injecting the sample into the disc.
- Do not turn the disc upside down.
- This is a quantitative test; therefore no visual interpretation of the results should be made.

STORAGE AND STABILITY

- Store the IB10 sphingotest® Troponin-99 between 2 and 8 °C (35 to 46 °F) until the expiration date printed on the pouch is reached.
- The IB10 sphingotest® Troponin-99 in its sealed pouch is stable at 18 to 30 °C/64 to 86 °F for 30 days, provided the expiration date printed on the pouch is not exceeded.

SAMPLE COLLECTION AND PREPARATION

- The IB10 sphingotest® Troponin-99 is to be run using lithium heparin whole blood or plasma samples.
- As cardiac proteins are relatively unstable, it is recommended that fresh samples be tested as soon as possible.
- Whole blood must be tested within 48 hours of collection.

- Whole blood and plasma samples can be transported and stored up to 48 hours at ambient temperature.
- Plasma samples should be kept frozen at -20 °C (-4 °F) or lower if longer storage is required.
- Allow samples to equilibrate to room temperature (19 to 25 °C/66 to 77 °F) prior to testing.
- For serial testing of the same patient during the 0–12 hour triage period, the same sample type (whole blood or plasma) should be used.

PROCEDURE

Nexus IB10 Analyzer

Consult the Nexus IB10 Analyzer User Manual

For Analyzer installation, start up and complete instructions for use refer to the **Nexus IB10 Analyzer User Manual**. Operator must consult the User Manual prior to use to become familiar with the proper operation and quality control procedures.

PERFORMING SYSTEM CHECK AND DISC CALIBRATION

Each time the Nexus IB10 Analyzer is turned on, a Self Check is automatically performed. The QR code on each Test disc contains information for disc calibration which the Analyzer automatically reads when running a test.

RUNNING QC WITH EXTERNAL CONTROLS

The manufacturer recommends the use of commercially available Troponin-99 Controls (please refer to section **Materials/Equipment Required But Not Provided**). Ensure that the Troponin-99 Controls are handled and prepared according to the corresponding Instructions for Use (IFU).

1. Remove an unopened Test pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes prior to testing.
 2. Open the pouch and remove the Test disc.
 3. Place the Test disc on a level surface.
 4. On the Nexus IB10 Analyzer press **New Analysis**.
 5. The Analyzer will perform a general system check.
 6. Enter the external control ID manually (up to 14 characters can be used for the ID) or enter the external control ID by using the barcode scanner.
 7. Mix external quality control vial by gently inverting the vial several times before sampling.
 8. **Testing External Control Sample on the IB10 sphingotest® Troponin-99.**
 - Using a precision pipette (fixed or adjusted to 500 µL) slowly draw well mixed external quality control sample into the pipette tip.
 - Positioning the tapered pipette tip at a 45° angle, pierce the X on the red dot to expose the sample channel inlet.
 - Slowly express the external control sample into the inlet applying minimum, but continuous force, on the pipette plunger.
 - Express sample to the first stop on the pipette, at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
 - Press **OK** on the Nexus IB10 Analyzer display.
 - When the tray opens, insert the filled Test disc into the tray and press **Run**.
 - In 20 minutes, the Nexus IB10 Analyzer will display the results on the screen.
 - Results will print automatically (if selected during Set Up) or press **Print**.
 - When the test is complete, analyze and compare the result with the expected value reported in the external control IFU for the external control level as measured using the IB10 sphingotest® Troponin-99.
 - Remove the Test disc and discard in appropriate receptacle.
 - If the external control result is outside the expected range, refer to the Quality Control Section below.
- Note: If the test run is cancelled before a test result is displayed, the Test disc cannot be reused and should be disposed of appropriately.*

TESTING PATIENT SAMPLES ON THE NEXUS IB10 ANALYZER

1. Remove an unopened Test pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
2. Open the pouch and remove the Test disc.
3. Place the Test disc on a level surface.
4. On the Nexus IB10 Analyzer press **New Analysis**.

5. The Analyzer will perform a general system check.
6. Enter the Patient ID manually (up to 14 characters can be used for the ID) or enter the Patient ID by using the barcode scanner.
7. Mix the whole blood patient sample tube by gently inverting the tube several times before testing.

8. Testing Patient Sample on the IB10 sphingotest® Troponin-99.

- Using a precision pipette (fixed or adjusted to 500 µL) slowly draw well mixed patient sample into the pipette tip.
- Positioning the tapered pipette at a 45° angle, pierce the X on the red dot to expose the sample channel inlet.
- Slowly express the patient sample into the inlet applying minimum, but continuous force, on the pipette plunger.
- Express sample to the first stop on the pipette at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
- Press **OK** on the Nexus IB10 Analyzer display.
- When the tray opens, insert the filled Test disc into the tray and press **Run**.
- In 20 minutes, the Nexus IB10 Analyzer will display the results on the screen.
- Results will print automatically (if selected during Set Up) or press **Print**.
- Remove the Test disc and discard in appropriate receptacle.

Note: If the test run is cancelled before a test result is displayed, the Test disc cannot be reused and should be disposed of appropriately.

INTERPRETATION OF RESULTS

The range of Troponin-99 concentration reported by the Nexus IB10 Analyzer is 0.05 ng/mL to 30.0 ng/mL. Results below or above this range will be shown as “< 0.05 ng/mL” or “> 30.0 ng/mL” respectively.

Quality Control

EXTERNAL CONTROLS

Good laboratory practice includes the use of external controls to ensure proper test performance. It is recommended that prior to using a new lot of IB10 sphingotest® Troponin-99, the performance of the lot should be confirmed by testing with external controls (see section **Materials/Equipment Required But Not Provided**) to ensure the Test will deliver the correct result. The frequency of quality control testing should be determined according to individual laboratory standard quality control procedures. Upon confirmation of the expected results, the Test discs are ready for use with patient samples. Controls should also be used any time the validity of test results is questionable. If external controls do not perform as expected, do not use the IB10 sphingotest® Troponin-99 and contact the Distributor in your area for Technical Assistance.

INTERNAL CONTROL

Each IB10 sphingotest® Troponin-99 has a built in positive procedural control. The Nexus IB10 Analyzer automatically determines the presence of this control thereby confirming that the test run has delivered a valid result. If the control does not form or is not recognized by the Analyzer, the test result is considered “invalid” and the test must be repeated.

LIMITATIONS

The results of the IB10 sphingotest® Troponin-99 should be used in conjunction with other available laboratory and clinical information.

Universally accepted Good Clinical Practice Guidelines recommend collecting serial samples during the early triage period from 0 to > 12 hours.^{4,9-11}

Other substances and/or factors not listed, e.g. technical or procedural errors, may interfere with the test and cause inaccurate results.

Nexus Dx, Inc. offers products for their intended use. Refer to the specific product literature for the intended use statements for each product. Product claims are subject to change. Nexus Dx, Inc. expressed and implied warranties (inclusive of implied warranties of merchantability and fitness) are conditional upon adherence to, or observance of Nexus Dx, Inc. published directions with regard to the use of Nexus Dx, Inc. products. Under no circumstances will Nexus Dx, Inc. be liable for any indirect or consequential damages.

For Technical Assistance please contact your local Distributor.

Performance Characteristics

MEASURING RANGE

The IB10 sphingotest® Troponin-99 has been demonstrated to provide measurable results at Troponin-99 levels of 0.05 ng/mL to 30 ng/mL.

ANALYTICAL SENSITIVITY

The LoD (Limit of Detection) of the IB10 sphingotest® Troponin-99 is 0.05 ng/mL, determined according to Clinical and Laboratory Standards Institute (CLSI, formerly NCCLS) Guideline EP17-A.¹² The proportions of false positives (α) less than 5% and false negatives (β) are less than 5% (limit of blank = 0.02 ng/mL).

LoQ (Limit of Quantitation) is the lowest cTnI concentration that can be reproducibly measured with a total coefficient of variation of at most 15%. It was determined to be 0.1 ng/mL.

CROSS-REACTIVITY AND INTERFERING SUBSTANCES

Drugs

The following drugs were tested for potential interference in the IB10 sphingotest® Troponin-99 (Table 1). The list of drugs encompasses common prescription and over-the-counter medications, as well as those often prescribed in a cardiac patient population. The drugs were tested at concentrations as recommended in the CLSI Approved Guideline EP7-A2 "Interference Testing in Clinical Chemistry,"¹³ or at least three times the highest concentration reported following a therapeutic dose. No significant interference with the IB10 sphingotest® Troponin-99 measurement was observed for the drugs listed in the table below.

Table 1.

Drug	Drug	Drug
Acetaminophen	Caffeine	Methyl-DOPA
Acetylsalicylic acid (aspirin)	Captopril	Nifedipine
Allopurinol	Digoxin	Phenytoin
Ampicillin	Dopamine	Theophylline
Ascorbic acid (vitamin C)	Erythromycin	Verapamil
Atenolol	Furosemide	

PROTEINS, PEPTIDES AND ENDOGENOUS SUBSTANCES

The following proteins, peptides and endogenous substances were tested for potential cross-reactivity in the IB10 sphingotest® Troponin-99 at the maximum concentration of substance indicated (Table 2). No substance demonstrated significant cross-reactivity using the method provided in CLSI EP7-A2.

Table 2.

Substance	Maximum Concentration
Cardiac Troponin T	1000 ng/mL
Cardiac Troponin C	1000 ng/mL
Skeletal Troponin-99	250 ng/mL
Endogenous substances	
Albumin	2.5 g/dL
Hemoglobin	0.1 g/dL
Rheumatoid Factor (RF)	203 IU/mL
Creatinine	2 mg/dL
Bilirubin	2.5 mg/dL
Triglycerides	0.5 g/dL
Urea	18 mg/dL
Cholesterol	280 mg/dL
Biotin	500 ng/mL

HOOK EFFECT

No high dose Hook Effect was observed for cTnI concentrations up to 500 ng/mL.

PRECISION

The precision of IB10 sphingotest® Troponin-99 was determined using samples where cTnI was added to normal human plasma at three concentrations (Table 3). The within-day and total precision were performed in two runs per day, in replicates of 4 per run at each concentration level over a 15 day period for a total number of repetitions of 120 at each concentration level. The within-run, total variances and coefficients of variation (CVs) were computed according to CLSI guideline EP5-A2.¹⁴

Table 3.

Sample	Mean level (ng/mL)	Within-Run Precision		Total Precision	
		Std. dev. (ng/mL)	CV (%)	Std. dev. (ng/mL)	CV (%)
1	0.46	0.06	12.6	0.06	12.6
2	0.95	0.11	12.0	0.12	13.0
3	4.43	0.41	9.4	0.45	10.1

WHOLE BLOOD Vs. PLASMA CORRELATION

A comparison study was performed using whole blood and plasma samples. When performing a Passing-Bablok regression analysis comparing the whole blood concentrations versus the corresponding plasma concentrations from the same subject samples, a slope of 1.00 (95% C.I. = [.93-1.07]) and an intercept of 0.01 ng/mL were found.

Expected Values

UPPER REFERENCE LIMIT

From a population of 224 individuals, the IB10 sphingotest® Troponin-99 was used to determine the concentration upper reference limits of cTnI. This population included apparently healthy individuals. The 99th percentile upper reference limit as determined with the samples is 0.10 ng/mL cTnI. Each laboratory should establish a reference range that represents the patient population that is to be evaluated at their facility.

METHOD COMPARISON

An equivalence study was performed between the IB10 sphingotest® Troponin-99 using the Nexus IB10 Analyzer, and the Ortho VITROS® Troponin I ES Test system. A total of 253 samples were tested within a concentration range of cTnI of 0.05 ng/mL to 30 ng/mL. The Passing-Bablok regression was:

$$\text{IB10 sphingotest® Troponin-99} = 0.80 \text{ (VITROS TnI ES Test)} + 0.00 \text{ ng/mL}$$

Correlation coefficient, $r = 0.92$

REFERENCES

1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. Circulation 2012;125:e2-e220.
2. Cardiovascular Diseases, World Health Organization factsheet. <http://www.who.int/mediacentre/factsheets/fs317/en/index.html>
3. World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria of diagnosis of ischemic heart disease. Circulation 1979;59:607-9.
4. Thygesen K, Alpert JS, White HD, Joint ESC/ACCF/AHA/WHF Task force for the Redefinition of Myocardial Infarction. Universal definition of myocardial infarction. Eur Heart J 2007;28:2525-38; Circulation 2007;116:2634-53; J Am Coll Cardiol 2007;50:2173-95.
5. Larue C, Calzolari C, Bertinchant J-P, et al. Cardiac specific immunoenzymatic assay of Troponin-99 in the early phase of acute myocardial infarction. Clin Chem 1993;39:972-9.
6. Apple FS, Voss E, Lund L, et al. Cardiac troponin CK-MB and myoglobin for the early detection of acute myocardial infarction and monitoring of reperfusion following thrombolytic therapy. Clin Chim Acta 1995;237:59-66.
7. Adams JE, Bodor GS, Davila-Roman VG, et al. Cardiac troponin-99: a marker with high specificity for cardiac injury. Circulation 1993;88:101-6.
8. Heidenreich PA, Alloggiamento T, Melsop K, et al. The prognostic value of troponin-99 in patients with non-ST elevation acute coronary syndrome: a meta analysis. J Am Coll Cardiol 2001;38:478-85.

9. Morrow DA, Rifai N, Tanasjevic MJ, et al. Clinical efficacy of three assays for cardiac troponin-I for risk stratification in acute coronary syndromes- a TIMI 11b substudy. *Clin Chem* 2000;46:453-60.
10. Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation* 2011;124:136-45.
11. Jaffe AS, Apple FS, Morrow A, et al. Being rational about (im)precision: a statement from the biochemistry subcommittee of the Joint European Society of Cardiology/ American College of Cardiology Foundation/ American Heart Association/World Heart Federation Task Force for the Definition of Myocardial Infarction. *Clin Chem* 2010;56:941-3.
12. Tholen DW, Linnet K, Kondratovich M et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.
13. McEnroe, RJ, Burnitt, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
14. Tholen DW, Kallner A, Kennedy JW et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline – Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004.

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Korrelationskoeffizient, r = 0,92

REFERENZEN

1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. *Circulation* 2012;125:e2-e220.
2. Cardiovascular Diseases, World Health Organization factsheet. <http://www.who.int/mediacentre/factsheets/fs317/en/index.html>
3. World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria of diagnosis of ischemic heart disease. *Circulation* 1979;59:607-9.
4. Thygesen K, Alpert JS, White HD, Joint ESC/ACCF/AHA/WHF Task force for the Redefinition of Myocardial Infarction. Universal definition of myocardial infarction. *Eur Heart J* 2007;28:2525-38; *Circulation* 2007;116:2634-53; *J Am Coll Cardiol* 2007;50:2173-95.
5. Larue C, Calzolari C, Bertinchant J-P, et al. Cardiac specific immunoenzymatic assay of Troponin-99 in the early phase of acute myocardial infarction. *Clin Chem* 1993;39:972-9.
6. Apple FS, Voss E, Lund L, et al. Cardiac troponin CK-MB and myoglobin for the early detection of acute myocardial infarction and monitoring of reperfusion following thrombolytic therapy. *Clin Chim Acta* 1995;237:59-66.
7. Adams JE, Bodor GS, Davila-Roman VG, et al. Cardiac troponin-99: a marker with high specificity for cardiac injury. *Circulation* 1993;88:101-6.
8. Heidenreich PA, Alloggiamento T, Melsop K, et al. The prognostic value of troponin-99 in patients with non-ST elevation acute coronary syndrome: a meta analysis. *J Am Coll Cardiol* 2001;38:478-85.
9. Morrow DA, Rifai N, Tanasijevic MJ, et al. Clinical efficacy of three assays for cardiac troponin-99 for risk stratification in acute coronary syndromes- a TIMI 11b substudy. *Clin Chem* 2000;46:453-60.
10. Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation* 2011;124:136-45.
11. Jaffe AS, Apple FS, Morrow A, et al. Being rational about (im)precision: a statement from the biochemistry subcommittee of the Joint European Society of Cardiology/ American College of Cardiology Foundation/ American Heart Association/World Heart Federation Task Force for the Definition of Myocardial Infarction. *Clin Chem* 2010;56:941-3.
12. Tholen DW, Linnet K, Kondratovich M et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.
13. McEnroe, RJ, Burnitt, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
14. Tholen DW Kallner A, Kennedy JW et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline – Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004.

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Coefficient de corrélation, $r = 0,92$

RÉFÉRENCES

1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. *Circulation* 2012;125:e2-e220.
2. Cardiovascular Diseases, World Health Organization factsheet. <http://www.who.int/mediacentre/factsheets/fs317/en/index.html>
3. World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria of diagnosis of ischemic heart disease. *Circulation* 1979;59:607-9.
4. Thygesen K, Alpert JS, White HD, Joint ESC/ACCF/AHA/WHF Task force for the Redefinition of Myocardial Infarction. Universal definition of myocardial infarction. *Eur Heart J* 2007;28:2525-38; *Circulation* 2007;116:2634-53; *J Am Coll Cardiol* 2007;50:2173-95.
5. Larue C, Calzolari C, Bertinchant J-P, et al. Cardiac specific immunoenzymatic assay of Troponin-99 in the early phase of acute myocardial infarction. *Clin Chem* 1993;39:972-9.
6. Apple FS, Voss E, Lund L, et al. Cardiac troponin CK-MB and myoglobin for the early detection of acute myocardial infarction and monitoring of reperfusion following thrombolytic therapy. *Clin Chim Acta* 1995;237:59-66.
7. Adams JE, Bodor GS, Davila-Roman VG, et al. Cardiac troponin-99: a marker with high specificity for cardiac injury. *Circulation* 1993;88:101-6.
8. Heidenreich PA, Alloggiamento T, Melsop K, et al. The prognostic value of troponin-99 in patients with non-ST elevation acute coronary syndrome: a meta analysis. *J Am Coll Cardiol* 2001;38:478-85.
9. Morrow DA, Rifai N, Tanasjevic MJ, et al. Clinical efficacy of three assays for cardiac troponin-99 for risk stratification in acute coronary syndromes- a TIMI 11b substudy. *Clin Chem* 2000;46:453-60.
10. Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation* 2011;124:136-45.
11. Jaffe AS, Apple FS, Morrow A, et al. Being rational about (im)precision: a statement from the biochemistry subcommittee of the Joint European Society of Cardiology/ American College of Cardiology Foundation/ American Heart Association/World Heart Federation Task Force for the Definition of Myocardial Infarction. *Clin Chem* 2010;56:941-3.
12. Tholen DW, Linnet K, Kondratovich M et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.
13. McEnroe, RJ, Burnett, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
14. Tholen DW, Kallner A, Kennedy JW et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline – Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004.

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BIBLIOGRAFIA

1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. *Circulation* 2012;125:e2-e220.
2. Cardiovascular Diseases, World Health Organization factsheet. <http://www.who.int/mediacentre/factsheets/fs317/en/index.html>
3. World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria of diagnosis of ischemic heart disease. *Circulation* 1979;59:607-9.
4. Thygesen K, Alpert JS, White HD, Joint ESC/ACCF/AHA/WHF Task force for the Redefinition of Myocardial Infarction. Universal definition of myocardial infarction. *Eur Heart J* 2007;28:2525-38; *Circulation* 2007;116:2634-53; *J Am Coll Cardiol* 2007;50:2173-95.
5. Larue C, Calzolari C, Bertinchant J-P, et al. Cardiac specific immunoenzymatic assay of Troponin-99 in the early phase of acute myocardial infarction. *Clin Chem* 1993;39:972-9.
6. Apple FS, Voss E, Lund L, et al. Cardiac troponin CK-MB and myoglobin for the early detection of acute myocardial infarction and monitoring of reperfusion following thrombolytic therapy. *Clin Chim Acta* 1995;237:59-66.
7. Adams JE, Bodor GS, Davila-Roman VG, et al. Cardiac troponin-99: a marker with high specificity for cardiac injury. *Circulation* 1993;88:101-6.
8. Heidenreich PA, Alloggiamento T, Melsop K, et al. The prognostic value of troponin-99 in patients with non-ST elevation acute coronary syndrome: a meta analysis. *J Am Coll Cardiol* 2001;38:478-85.
9. Morrow DA, Rifai N, Tanasjevic MJ, et al. Clinical efficacy of three assays for cardiac troponin-99 for risk stratification in acute coronary syndromes- a TIMI 11b substudy. *Clin Chem* 2000;46:453-60.
10. Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation* 2011;124:136-45.
11. Jaffe AS, Apple FS, Morrow A, et al. Being rational about (im)precision: a statement from the biochemistry subcommittee of the Joint European Society of Cardiology/ American College of Cardiology Foundation/ American Heart Association/World Heart Federation Task Force for the Definition of Myocardial Infarction. *Clin Chem* 2010;56:941-3.
12. Tholen DW, Linnet K, Kondratovich M et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.
13. McEnroe, RJ, Burnitt, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
14. Tholen DW Kallner A, Kennedy JW et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline – Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004.

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1. Roger VL, Go AS, Lloyd-Jones DM, on behalf of the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Heart disease and stroke statistics—2012 update: a report from the American Heart Association. Circulation 2012;125:e2-e220.
2. Cardiovascular Diseases, World Health Organization factsheet. <http://www.who.int/mediacentre/factsheets/fs317/en/index.html>
3. World Health Organization. Report of the Joint International Society and Federation of Cardiology/World Health Organization Task Force on Standardization of Clinical Nomenclature. Nomenclature and criteria of diagnosis of ischemic heart disease. Circulation 1979;59:607-9.
4. Thygesen K, Alpert JS, White HD, Joint ESC/ACCF/AHA/WHF Task force for the Redefinition of Myocardial Infarction. Universal definition of myocardial infarction. Eur Heart J 2007;28:2525-38; Circulation 2007;116:2634-53; J Am Coll Cardiol 2007;50:2173-95.
5. Larue C, Calzolari C, Bertinchant J-P, et al. Cardiac specific immunoenzymatic assay of Troponin-99 in the early phase of acute myocardial infarction. Clin Chem 1993;39:972-9.
6. Apple FS, Voss E, Lund L, et al. Cardiac troponin CK-MB and myoglobin for the early detection of acute myocardial infarction and monitoring of reperfusion following thrombolytic therapy. Clin Chim Acta 1995;237:59-66.
7. Adams JE, Bodor GS, Davila-Roman VG, et al. Cardiac troponin-99: a marker with high specificity for cardiac injury. Circulation 1993;88:101-6.
8. Heidenreich PA, Alloggiamento T, Melsop K, et al. The prognostic value of troponin-99 in patients with non-ST elevation acute coronary

- syndrome: a meta analysis. *J Am Coll Cardiol* 2001;38:478-85.
- 9. Morrow DA, Rifai N, Tanasjevic MJ, et al. Clinical efficacy of three assays for cardiac troponin-99 for risk stratification in acute coronary syndromes- a TIMI 11b substudy. *Clin Chem* 2000;46:453-60.
 - 10. Reichlin T, Irfan A, Twerenbold R, et al. Utility of absolute and relative changes in cardiac troponin concentrations in the early diagnosis of acute myocardial infarction. *Circulation* 2011;124:136-45.
 - 11. Jaffe AS, Apple FS, Morrow A, et al. Being rational about (im)precision: a statement from the biochemistry subcommittee of the Joint European Society of Cardiology/ American College of Cardiology Foundation/ American Heart Association/World Heart Federation Task Force for the Definition of Myocardial Infarction. *Clin Chem* 2010;56:941-3.
 - 12. Tholen DW, Linnet K, Kondratovich M et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.
 - 13. McEnroe, RJ, Burnett, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
 - 14. Tholen DW, Kallner A, Kennedy JW et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline – Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004.