English

IB10 sphingotest® beta-hCG

For the Quantitative Determination of Human Chorionic Gonadotropin in
Lithium Heparin Human Blood and Plasma

Explanation of Symbols

	CE Mark of Conformity
	Manufacturer
	Catalog Number
	Expiry date/Use by
	Lot number
	<i>In Vitro</i> Diagnostic medical device
	Consult instructions for use
	Store between 2 °C and 8 °C
	European Authorized Representative
	Contains sufficient for <n> tests
	Do not reuse
	Serial Number
	Insert disc label side up

Nexus Dx, Inc.,
6759 Mesa Ridge Road
San Diego, California 92121 USA
Telephone: 1 (858) 410 4600
Fax: 1 (858) 410 4700

TheraGenesis GmbH
Bahnhofstrasse 5
55276 Oppenheim, Germany
Tel: +49 (0) 151 506 403 14

022-00066 REV B

IB10 sphingotest® beta-hCG

For *In Vitro* Diagnostic Use

INTENDED USE

The IB10 sphingotest® beta-hCG is a rapid point-of-care (POC) immunoassay for the *in vitro* quantitative determination of human chorionic gonadotropin (hCG) in lithium heparin whole blood and plasma. The IB10 sphingotest® beta-hCG is intended for use in conjunction with the Nexus IB10 Analyzer and provides quantitative results in 20 minutes.

This Test detects total hCG, measuring both the intact hCG molecule and its free beta subunit. The quantitative measurement of hCG aids in the early detection of pregnancy.

This Test is designed for professional use only and may be used in hospital central laboratories and in alternate care settings such as emergency departments, critical care units, and other sites where near patient testing is practiced.

SUMMARY AND EXPLANATION OF THE TEST

hCG is a glycoprotein hormone composed of 2 dissimilar subunits designated as α and β . It is produced by trophoblast cells of the placenta shortly after implantation of a fertilized ovum into the uterine wall.^{1,2} hCG is the primary analyte used for pregnancy confirmation and monitoring due to its rapid rise in both blood and urine soon after conception.

The detection of hCG within 3–4 weeks of the last menstrual period (LMP) is the most reliable indicator for the confirmation of pregnancy. During a normal pregnancy, levels of hCG in the blood vary but are approximately 25–50 mIU/mL in the week after conception and rise exponentially doubling every 1.5–3 days during the first six weeks.³ hCG levels continue to rise through the end of the first trimester, followed by a slow decline as the pregnancy reaches full-term (~ 40 weeks).

In addition to its synthesis during pregnancy, elevated hCG levels are also associated with trophoblastic disease and certain non-trophoblastic neoplasms. The possibility of having these diseases must be eliminated before a diagnosis of pregnancy can be made.^{4–5} A small amount of hCG is normally produced by the pituitary gland in conjunction with the structurally similar glycoprotein hormones, which include luteinizing hormone (LH), follicle-stimulating hormone (FSH) and thyroid-stimulating hormone (TSH).^{6,7}

PRINCIPLE

The Nexus Dx immunochemistry system combines chemistry with microfluidics and centrifugal flow to rapidly prepare a cell free plasma from whole blood that can then be moved through a channel to rehydrate, solubilize and mix with freeze dried immunoconjugates. Using a combination of active flow and capillary action, the test sample is quantitatively measured in 20 minutes with an optical signal level proportional to the analyte(s) concentration.

After addition of the patient sample, the entire test is performed within the Nexus IB10 Analyzer which provides control of the temperature of the disc, as well as the reaction sequence, centrifugal flow, mixing, incubation time, final signal measurement, quantitation and reporting of results. The Test disc includes a positive internal control to ensure that the Test disc has operated properly. Each lot is calibrated to provide accurate analyte concentration.

Lot specific calibration along with additional information such as the lot expiration date is contained on a QR code label affixed to each disc. It is recommended that external controls be tested at appropriate time intervals to confirm that the system and test lot are performing within acceptable limits.

REAGENTS

The IB10 sphingotest® beta-hCG contains all required reagents to evaluate the level of hCG including dye conjugated monoclonal anti-hCG antibody, biotin conjugated monoclonal anti-hCG antibody and streptavidin immobilized at the detection area on the disc.

MATERIALS PROVIDED

Each box contains the following:

- 10 IB10 sphingotest® beta-hCG discs, each individually sealed in a foil pouch with a desiccant.
- 1 Instructions for Use (IFU).

MATERIALS/EQUIPMENT REQUIRED BUT NOT PROVIDED

1. Nexus IB10 Analyzer - Model #BCA-IB10.
2. Commercially available hCG Controls for external Quality Control (QC). **Contact the Distributor in your area** for recommended external QC materials or related technical assistance.
3. Calibrated reusable fixed or variable volume pipette gun with high precision and accuracy capable of delivering 500 μ L of whole blood or plasma.
4. Disposable pipette tips capable of accepting and delivering 500 μ L of whole blood or plasma.

PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Carefully follow the Instructions for Use.
- Prior to testing controls or patient samples, ensure that the Analyzer software is updated to the latest version. (Refer to the Nexus IB10 manual for specific instructions).
- Wear disposable gloves while handling samples.
- Handle samples with care. Samples and used Test discs should be treated as potentially infectious and should be discarded as biohazardous material according to local regulations.
- Thoroughly wash hands following handling.
- The result obtained from the IB10 sphingotest® beta-hCG does not provide a definitive diagnosis and should be interpreted by a physician in conjunction with other laboratory test results and according to current medical guidelines and patient clinical findings.
- Keep the Test disc in the sealed pouch until ready for use.
- Do not use the Test disc if the pouch is damaged or the seal is broken.
- Do not use the Test disc after the expiration date printed on the pouch.
- Prior to use, place the unopened pouch at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
- Always pay attention to cleanliness when handling the Test. Avoid any contamination from fingerprints or foreign substances. Do not contaminate the sample channel inlet.
- Do not drop or damage the Test disc.
- The Test disc should be inserted with the label side up, into the Nexus IB10 Analyzer tray immediately after injecting the sample into the disc.
- Do not turn the disc upside down.
- This is a quantitative test; therefore no visual interpretation of the results should be made.

STORAGE AND STABILITY

- Store the IB10 sphingotest® beta-hCG between 2 and 8 °C (35 to 46 °F) until the expiration date printed on the pouch is reached.
- The IB10 sphingotest® beta-hCG in its sealed pouch is stable at 18 to 30 °C/64 to 86 °F for 30 days, provided the expiration date printed on the pouch is not exceeded.

SAMPLE COLLECTION AND PREPARATION

- The IB10 sphingotest® beta-hCG is to be run using lithium heparin whole blood or plasma samples.
- It is recommended that samples be tested as soon as possible.
- Whole blood must be tested within 24 hours of collection.
- Plasma samples should be kept frozen at -20 °C (- 4 °F) or lower if longer storage is required.
- Allow samples to equilibrate to room temperature (19 to 25 °C/66 to 77 °F) prior to testing.

PROCEDURE

Nexus IB10 Analyzer

 **Consult the Nexus IB10 Analyzer User Manual**

For Analyzer installation, start up and complete instructions for use refer to the **Nexus IB10 Analyzer User Manual**. Operator must consult the User Manual prior to use to become familiar with the proper operation and quality control procedures.

PERFORMING SYSTEM CHECK AND DISC CALIBRATION

Each time the Nexus IB10 Analyzer is turned on, a Self Check is automatically performed. The QR code on each Test disc contains information for disc calibration which the Analyzer automatically reads when running a test.

RUNNING QC WITH EXTERNAL CONTROLS

The manufacturer recommends the use of commercially available hCG Controls (please refer to section **Materials/Equipment Required But Not Provided**). Ensure that the external controls are handled and prepared according to the corresponding Instructions for Use (IFU).

1. Remove an unopened Test pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes prior to testing.
2. Open the pouch and remove the Test disc.
3. Place the Test disc on a level surface.
4. On the Nexus IB10 Analyzer press **New Analysis**.
5. The Analyzer will perform a general system check.
6. Enter the external control ID manually (up to 14 characters can be used for the ID) or enter the external control ID by using the barcode scanner.
7. Mix the external quality control vial by gently inverting the vial several times before sampling.

8. Testing External Control Sample on the IB10 sphingotest® beta-hCG.

- Using a precision pipette (fixed or adjusted to 500 µL) slowly draw well-mixed external quality control sample into the pipette tip.
- Positioning the tapered pipette tip at a 45° angle, pierce the **X** on the red dot to expose the sample channel inlet.
- Slowly express the external control sample into the inlet applying minimum, but continuous force, on the pipette plunger.
- Express sample to the first stop on the pipette, at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
- Press **QC** on the Nexus IB10 Analyzer display.
- When the tray opens, insert the filled Test disc into the tray and press **Run**.
- The tray will close and perform a disc validity check.
- A screen appears for selecting quality control materials (please refer to the **Quality Control Settings** section of the **Nexus IB10 Analyzer User Manual** for how to update quality control material [external controls]).
- Select the quality control material to be tested.
- Press **OK** on the Nexus IB10 Analyzer display.
- In 20 minutes, the Nexus IB10 Analyzer will display the results on the screen.
- Results will print automatically (if selected during Set Up) or press **Print**.
- When the test is complete, analyze and compare the result with the expected value reported in the external control IFU for the external control level as measured using the IB10 sphingotest® beta-hCG.
- Remove the Test disc and discard in appropriate receptacle.
- If the external control result is outside the expected range, refer to the Quality Control section below.

Note: If the test run is cancelled before a Test result is displayed, the Test disc cannot be reused and should be disposed of appropriately.

TESTING PATIENT SAMPLES ON THE NEXUS IB10 ANALYZER

1. Remove an unopened Test pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
2. Open the pouch and remove the Test disc.
3. Place the Test disc on a level surface.
4. On the Nexus IB10 Analyzer press **New Analysis**.
5. The Analyzer will perform a general system check.
6. Enter the Patient ID manually (up to 14 characters can be used for the ID) or enter the Patient ID by using the barcode scanner.
7. Mix the whole blood patient sample by gently inverting the tube several times before testing.

8. Testing Patient Sample on the IB10 sphingotest® beta-hCG.

- Using a precision pipette (fixed or adjusted to 500 µL) slowly draw well-mixed patient sample into the pipette tip.
- Positioning the tapered pipette tip at a 45° angle, pierce the **X** on the red dot to expose the sample channel inlet.
- Slowly express the patient sample into the inlet applying minimum, but continuous force, on the pipette plunger.
- Express sample to the first stop on the pipette, at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
- Press **OK** on the Nexus IB10 Analyzer display.
- When the tray opens, insert the filled Test disc into the tray and press **Run**.
- In 20 minutes, the Nexus IB10 Analyzer will display the results on the screen.
- Results will print automatically (if selected during Set Up) or press **Print**.
- Remove the Test disc and discard in appropriate receptacle.

Note: If the test run is cancelled before a test result is displayed, the Test disc cannot be reused and should be disposed of appropriately.

TRACEABILITY OF CALIBRATION

Calibration of the IB10 sphingotest® beta-hCG is traceable by in-house reference calibrators assigned by the 5th WHO International Standard for hCG (07/364), (National Institute of Biological Standards and Controls -NIBSC).

INTERPRETATION OF RESULTS

The range of hCG concentration reported by the Nexus IB10 Analyzer is 4.0 mIU/mL to 400 mIU/mL. Results below or above this range will be shown as “< 4.0 mIU/mL” or “> 400 mIU/mL” respectively.

- **Recommended Decision Threshold Values:**

hCG results greater than or equal to 25 mIU/mL are considered positive. See section - **Limitations** for additional guidance in interpreting test results.

Quality Control

EXTERNAL CONTROLS

Good laboratory practice includes the use of external controls to ensure proper test performance. It is recommended that prior to using a new lot of IB10 sphingotest® beta-hCG, the performance of the lot should be confirmed by testing with external controls (see section **Materials/Equipment Required But Not Provided**) to ensure the Test will deliver the correct result. The frequency of quality control testing should be determined according to individual laboratory standard quality control procedures. Upon confirmation of the expected results, the Test discs are ready for use with patient samples. Controls should also be used any time the validity of test results is questionable. If external controls do not perform as expected, do not use the IB10 sphingotest® and contact the Distributor in your area for Technical Assistance.

INTERNAL CONTROL

Each IB10 sphingotest® beta-hCG has a built in positive procedural control. The Nexus IB10 Analyzer automatically determines the presence of this control thereby confirming that the test run has delivered a valid result. If the control does not form or if it is not recognized by the Analyzer, the test result is considered “invalid” and the test must be repeated.

LIMITATIONS

- For diagnostic purposes, the test result should always be evaluated and interpreted in conjunction with the patient’s medical history, symptoms and other clinical information. Ectopic pregnancy cannot be distinguished from normal pregnancy by hCG measurements alone.⁸
- If the hCG level is inconsistent with clinical evidence, the results should be confirmed by an alternative method. This may include a urine hCG test or repeating blood or plasma measurements on a different test system.⁹ The absence of urinary hCG may suggest a falsely elevated blood or plasma result.
- As with any assay employing mouse antibodies, the possibility exists for interference by heterophilic antibodies in the sample.^{10,11} Specimens from individuals who have been regularly exposed to animals or who have been treated with animal products may contain these antibodies.
- Falsely elevated or falsely depressed hCG results might be caused by the presence of heterophilic antibodies, such as human anti-mouse antibodies (HAMA), nonspecific protein binding or hCG-like substances.
- A high dose hook effect was not observed up to 700,000 mIU/mL.
- Elevated hCG levels have been reported for the patients with trophoblastic or non-trophoblastic neoplasms.
- Persistent low but elevated levels of hCG may occur in apparently healthy, non-pregnant subjects.¹²
- hCG levels less than 25 mIU/mL do not exclude pregnancy. Because the concentration of hCG rises rapidly during the early stage of pregnancy, further testing with samples obtained 48 hours later is recommended if pregnancy is suspected.¹³
- Positive hCG levels may remain detectable for several weeks after delivery or abortion.¹⁴
- Specimens which test positive during the initial days after conception may be negative later due to natural termination of the pregnancy.
- Peri- and postmenopausal women may elicit low but elevated hCG levels unrelated to pregnancy.¹⁵

Other substances and/or factors not listed, e.g. technical or procedural errors, may interfere with the test and cause inaccurate results.

Nexus Dx, Inc. offers products for their intended use. Refer to the specific product literature for the intended use statements for each product.

Product claims are subject to change. Nexus Dx, Inc.’s expressed and implied warranties (inclusive of implied warranties of merchantability and fitness) are conditional upon adherence to, or observance of Nexus Dx, Inc. published directions with regard to the use of Nexus Dx, Inc. products. Under no circumstances will Nexus Dx, Inc. be liable for any indirect or consequential damages.

For Technical Assistance please contact your local Distributor.

EXPECTED VALUES

Each laboratory should establish its own expected values that represents the population that is to be evaluated at their facility. As a guide, the following ranges were determined. The IB10 sphingotest® beta-hCG was used to determine the concentration upper reference limit of hCG in human plasma samples collected from apparently healthy, non-pregnant individuals. The 95th percentile upper reference limit as determined with the samples is 5.42 mIU/mL hCG.

Reference Group	N	Median (mIU/mL)	95 th percentile (mIU/mL)
Male	150	0.00	4.75
Non-pregnant female (all ages)	248	1.56	5.75
Female age ≤ 50	150	0.43	4.15
Female age > 50	98	3.26	7.91
Total	398	0.48	5.42

Representative hCG ranges during normal pregnancy based on Last Menstrual Period (LMP) are summarized below. Because other clinical reference citations may show different values, each laboratory should establish its own reference ranges.

After LMP (weeks)	hCG Range (mIU/mL) ¹⁶
4	5-100
5	200-3000
6	10,000-80,000
7-14	90,000-500,000
15-26	5,000-8,000
27-40	3,000-15,000

Performance Characteristics

MEASURING RANGE

The IB10 sphingotest® beta-hCG provides accurate measurements of hCG within 4.0 mIU/mL – 400 mIU/mL interval.

ANALYTICAL SENSITIVITY

The LoD (Limit of Detection) of the IB10 sphingotest® beta-hCG is 4.0 mIU/mL, and was determined according to Clinical and Laboratory Standards Institute (CLSI) Guideline EP17-A.¹⁷ The proportion of false positives (α) and false negatives (β) are less than 5% (limit of blank = 1.70 mIU/mL).

LoQ (Limit of Quantitation) is the lowest hCG concentration that can be reproducibly measured with a total coefficient of variation of at most 15%. It was determined to be 5.0 mIU/mL.

CROSS-REACTIVITY AND INTERFERING SUBSTANCES

Potential Interferents

The following potentially interfering substances were tested for cross-reactivity and interference in the IB10 sphingotest® beta-hCG at the maximum concentration of substance indicated (Table 1). No substance demonstrated significant cross-reactivity or interference using the method provided in the CLSI EP7-A2.¹⁸

Table 1.

Potential Interferents	Maximum Concentration
FSH	400 mIU/mL
LH	400 mIU/mL
TSH	250 µIU/mL
Albumin	5 g/dL

Potential Interferents	Maximum Concentration
Hemoglobin	0.1 g/dL
Rheumatoid Factor (RF)	220 IU/mL
Creatinine	2 mg/dL
Bilirubin	2.5 mg/dL
Triglycerides (Mixtures)	0.5 g/dL
Urea	18 mg/dL
Cholesterol	280 mg/dL
Biotin	400 ng/mL

HOOK EFFECT

No high dose Hook Effect was observed for hCG concentrations up to 700,000 mIU/mL.

PRECISION

The precision of the IB10 sphingotest® beta-hCG was determined using samples where hCG was added to normal human plasma at three concentrations (Table 2). The within-day and total precision measurements were performed in two runs per day, in replicates of 5 per run at each concentration level over a 20 day period for a total number of repetitions of 200 at each concentration level. The within-run and total precision were computed according to the CLSI guideline EP5-A2.¹⁹

Table 2.

Sample	Mean (mIU/mL)	Within-Run Precision		Total Precision	
		Std. dev. (mIU/mL)	CV%	Std. dev. (mIU/mL)	CV%
1	97.0	1.30	13.4	1.34	13.8
2	24.31	3.03	12.4	3.08	12.7
3	190.52	25.46	13.4	25.96	13.6

WHOLE BLOOD Vs. PLASMA COMPARISON

A comparison study was performed using whole blood and plasma samples. When performing a Passing-Bablok regression analysis comparing the whole blood concentrations versus the corresponding plasma concentrations from the same subject samples (N=30), a slope of 1.05 (95% C.I. = [1.00-1.09]) and an intercept of -0.51 mIU/mL were obtained.

METHOD COMPARISON

An equivalence study was performed between the IB10 sphingotest® beta-hCG using the Nexus IB10 Analyzer, and the Ortho VITROS® Total β-hCG II Test system. A total of 215 samples were tested within a concentration range of hCG of 4.0 mIU/mL to 400 mIU/mL. The Passing-Bablok regression was:

$$\text{IB10 sphingotest® beta-hCG} = 0.96 \text{ (VITROS Total } \beta\text{-hCG Test)} + 1.75 \text{ mIU/mL}$$

Rank correlation coefficient = 0.99

REFERENCES

1. Braunstein GD, Rasor J, Adler D, et al. Serum human chorionic gonadotropin levels throughout normal pregnancy. Am J Obstet Gynecol 1976;126:678-81.
2. Brody S, Carlstrom G. Immuno-assay of human chorionic gonadotropin in normal and pathologic pregnancy. J Clin Endocrinol Metab 1962;22:564-74.
3. Vaitukaitis JL, Ross GT, Braunstein GT, Rayford PL. Gonadotropins and their subunits:basic and clinical studies Recent Prog Horm Res 1976;32:289-331.
4. Braunstein GD, Vaitukaitis JL, Carbone PP, Ross GT. Ectopic production of human chorionic gonadotropin by neoplasm. Ann Intern Med

- 1973;78:39-45.
- 5. Stenman UH, Alfthan H, Hotakainen K. Human chorionic gonadotropin in cancer. *Clin Biochem* 2004;37:549-61.
 - 6. Stenman UH, Alfthan H, Ranta T, et al. Serum levels of human chorionic gonadotropin in nonpregnant women and men are modulated by gonadotropin-releasing hormone and sex steroids. *Clin Endocrinol Metab* 1987;64:730-6.
 - 7. Birken S, Maydelman Y, Gawinowicz MA, et al. Isolation and characterization of human pituitary chorionic gonadotropin. *Endocrinology* 1996;137:1402-11.
 - 8. Kadar N, Caldwell BV, Romero R. A method of screening for ectopic pregnancy and its indications. *Obstet Gynecol* 1981;58:162-6.
 - 9. Cole LA. Phantom hCG and phantom choriocarcinoma. *Gynecol Oncol* 1998;71:325-9.
 - 10. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay interference. *J Clin Immunoassay* 1992;15:108-115.
 - 11. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45:942-956.
 - 12. Borkowski A, Muquardt C. Human chorionic gonadotropin in the plasma of normal, nonpregnant subjects. *N Engl J Med* 1979;301:298-302.
 - 13. Lenton EA, Neal LM, Sulaiman R. Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy. *Fertil Steril* 1982;37:773-8.
 - 14. Steier JA, Bergsjo P, Myking OL. Human chorionic gonadotropin in maternal plasma after induced abortion, spontaneous abortion and removed ectopic pregnancy. *Obstet Gynecol* 1984;64:391-4.
 - 15. Snyder JA, Haymond S, Parvin CA, et al. Diagnostic consideration in the measurement of human chorionic gonadotropin in aging women. *Clin Chem* 2005;51:1830-5.
 - 16. Wu, HBA. *Tietz clinical guide to laboratory tests*, 4th Ed. Philadelphia, 2006, Saunders-Elsevier, pp. 252-259.
 - 17. Tholen DW, Linnet K, Kondratovich M, et al. NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*, NCCLS document EP17-A, Volume 24 Number 35, 2004.
 - 18. McEnroe, RJ, Burnitt, MF, Powers, DM, et al. *CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition*. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
 - 19. Tholen DW, Kallner A, Kennedy JW, et al. NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline – Second Edition*. NCCLS document EP5-A2, Volume 25 Number 27, 2004

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Rang Korrelationskoeffizient = 0,99

REFERENZEN

1. Braunstein GD, Rasor J, Adler D, et al. Serum human chorionic gonadotropin levels throughout normal pregnancy. Am J Obstet Gynecol 1976;126:678-81.
2. Brody S, Carlstrom G. Immuno-assay of human chorionic gonadotropin in normal and pathologic pregnancy. J Clin Endocrinol Metab 1962;22:564-74.
3. Vaitukaitis JL, Ross GT, Braunstein GT, Rayford PL. Gonadotropins and their subunits:basic and clinical studies Recent Prog Horm Res 1976;32:289-331.
4. Braunstein GD, Vaitukaitis JL, Carbone PP, Ross GT. Ectopic production of human chorionic gonadotropin by neoplasm. Ann Intern Med 1973;78:39-45.
5. Stenman UH, Alfthan H, Hotakainen K. Human chorionic gonadotropin in cancer. Clin Biochem 2004;37:549-61.
6. Stenman UH, Alfthan H, Ranta T, et al. Serum levels of human chorionic gonadotropin in nonpregnant women and men are modulated by gonadotropin-releasing hormone and sex steroids. Clin Endocrinol Metab 1987;64:730-6.
7. Birken S, Maydelman Y, Gawinowicz MA, et al. Isolation and characterization of human pituitary chorionic gonadotropin. Endocrinology 1996;137:1402-11.
8. Kadar N, Caldwell BV, Romero R. A method of screening for ectopic pregnancy and its indications. Obstet Gynecol 1981;58:162-6.
9. Cole LA. Phantom hCG and phantom choriocarcinoma. Gynecol Oncol 1998;71:325-9.
10. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay interference. J Clin Immunoassay 1992;15:108-115.
11. Kricka LJ. Human anti-animal antibody interference in immunological assays. Clin Chem 1999;45:942-956.
12. Borkowski A, Muquardt C. Human chorionic gonadotropin in the plasma of normal, nonpregnant subjects. N Engl J Med 1979;301:298-302.
13. Lenton EA, Neal LM, Sulaiman R. Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy. Fertil Steril 1982;37:773-8.
14. Steier JA, Bergsjo P, Myking OL. Human chorionic gonadotropin in maternal plasma after induced abortion, spontaneous abortion and

- removed ectopic pregnancy. *Obstet Gynecol* 1984;64:391-4.
- 15. Snyder JA, Haymond S, Parvin CA, et al. Diagnostic consideration in the measurement of human chorionic gonadotropin in aging women. *Clin Chem* 2005;51:1830-5.
 - 16. Wu, HBA. Tietz clinical guide to laboratory tests, 4th Ed. Philadelphia, 2006, Saunders-Elsevier, pp. 252-259.
 - 17. Tholen DW, Linnet K, Kondratovich M, et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.
 - 18. McEnroe, RJ, Burnitt, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
 - 19. Tholen DW Kallner A, Kennedy JW et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline – Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004..

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RÉFÉRENCES

- Braunstein GD, Rasor J, Adler D, et al. Serum human chorionic gonadotropin levels throughout normal pregnancy. *Am J Obstet Gynecol* 1976;126:678-81.
- Brody S, Carlstrom G. Immuno-assay of human chorionic gonadotropin in normal and pathologic pregnancy. *J Clin Endocrinol Metab* 1962;22:564-74.
- Vaitukaitis JL, Ross GT, Braunstein GT, Rayford PL. Gonadotropins and their subunits:basic and clinical studies *Recent Prog Horm Res* 1976;32:289-331.
- Braunstein GD, Vaitukaitis JL, Carbone PP, Ross GT. Ectopic production of human chorionic gonadotropin by neoplasm. *Ann Intern Med* 1973;78:39-45.
- Stenman UH, Alfthan H, Hotakainen K. Human chorionic gonadotropin in cancer. *Clin Biochem* 2004;37:549-61.
- Stenman UH, Alfthan H, Ranta T, et al. Serum levels of human chorionic gonadotropin in nonpregnant women and men are modulated by gonadotropin-releasing hormone and sex steroids. *Clin Endocrinol Metab* 1987;64:730-6.
- Birken S, Maydelman Y, Gawinowicz MA, et al. Isolation and characterization of human pituitary chorionic gonadotropin. *Endocrinology* 1996;137:1402-11.
- Kadar N, Caldwell BV, Romero R. A method of screening for ectopic pregnancy and its indications. *Obstet Gynecol* 1981;58:162-6.
- Cole LA. Phantom hCG and phantom choriocarcinoma. *Gynecol Oncol* 1998;71:325-9.
- Levinson SS. The nature of heterophilic antibodies and their role in immunoassay interference. *J Clin Immunoassay* 1992;15:108-115.
- Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45:942-956.
- Borkowski A, Muquardt C. Human chorionic gonadotropin in the plasma of normal, nonpregnant subjects. *N Engl J Med* 1979;301:298-302.
- Lenton EA, Neal LM, Sulaiman R. Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy. *Fertil Steril* 1982;37:773-8.
- Steier JA, Bergsjo P, Myking OL. Human chorionic gonadotropin in maternal plasma after induced abortion, spontaneous abortion and removed ectopic pregnancy. *Obstet Gynecol* 1984;64:391-4.
- Snyder JA, Haymond S, Parvin CA, et al. Diagnostic consideration in the measurement of human chorionic gonadotropin in aging women. *Clin Chem* 2005;51:1830-5.
- Wu, HBA. Tietz clinical guide to laboratory tests, 4th Ed. Philadelphia, 2006, Saunders-Elsevier, pp. 252-259.
- Tholen DW, Linnet K, Kondratovich M, et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.

18. McEnroe, RJ, Burnitt, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
19. Tholen DW Kallner A, Kennedy JW et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline – Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004..

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RIFERIMENTI

1. Braunstein GD, Rasor J, Adler D, et al. Serum human chorionic gonadotropin levels throughout normal pregnancy. *Am J Obstet Gynecol* 1976;126:678-81.
2. Brody S, Carlstrom G. Immuno-assay of human chorionic gonadotropin in normal and pathologic pregnancy. *J Clin Endocrinol Metab* 1962;22:564-74.
3. Vaitukaitis JL, Ross GT, Braunstein GT, Rayford PL. Gonadotropins and their subunits:basic and clinical studies *Recent Prog Horm Res* 1976;32:289-331.
4. Braunstein GD, Vaitukaitis JL, Carbone PP, Ross GT. Ectopic production of human chorionic gonadotropin by neoplasm. *Ann Intern Med* 1973;78:39-45.
5. Stenman UH, Alfthan H, Hotakainen K. Human chorionic gonadotropin in cancer. *Clin Biochem* 2004;37:549-61.
6. Stenman UH, Alfthan H, Ranta T, et al. Serum levels of human chorionic gonadotropin in nonpregnant women and men are modulated by gonadotropin-releasing hormone and sex steroids. *Clin Endocrinol Metab* 1987;64:730-6.
7. Birken S, Maydelman Y, Gawinowicz MA, et al. Isolation and characterization of human pituitary chorionic gonadotropin. *Endocrinology* 1996;137:1402-11.
8. Kadar N, Caldwell BV, Romero R. A method of screening for ectopic pregnancy and its indications. *Obstet Gynecol* 1981;58:162-6.
9. Cole LA. Phantom hCG and phantom choriocarcinoma. *Gynecol Oncol* 1998;71:325-9.
10. Levinson SS. The nature of heterophilic antibodies and their role of in immunoassay interference. *J Clin Immunoassay* 1992;15:108-115.
11. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45:942-956.
12. Borkowski A, Muquardt C. Human chorionic gonadotropin in the plasma of normal, nonpregnant subjects. *N Engl J Med* 1979;301:298-302.
13. Lenton EA, Neal LM, Sulaiman R. Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy. *Fertil Steril* 1982;37:773-8.
14. Steier JA, Bergsjo P, Myking OL. Human chorionic gonadotropin in maternal plasma after induced abortion, spontaneous abortion and removed ectopic pregnancy. *Obstet Gynecol* 1984;64:391-4.
15. Snyder JA, Haymond S, Parvin CA, et al. Diagnostic consideration in the measurement of human chorionic gonadotropin in aging women. *Clin Chem* 2005;51:1830-5.
16. Wu, HBA. *Tietz clinical guide to laboratory tests*, 4th Ed. Philadelphia, 2006, Saunders-Elsevier, pp. 252-259.
17. Tholen DW, Linnet K, Kondratovich M, et al. NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*, NCCLS document EP17-A, Volume 24 Number 35, 2004.
18. McEnroe, RJ, Burnitt, MF, Powers, DM, et al. *CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition*. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
19. Tholen DW Kallner A, Kennedy JW et al. NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline – Second Edition*. NCCLS document EP5-A2, Volume 25 Number 27, 2004.

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REFERÊNCIAS

- Braunstein GD, Rasor J, Adler D, et al. Serum human chorionic gonadotropin levels throughout normal pregnancy. *Am J Obstet Gynecol* 1976;126:678-81.
- Brody S, Carlstrom G. Immuno-assay of human chorionic gonadotropin in normal and pathologic pregnancy. *J Clin Endocrinol Metab* 1962;22:564-74.
- Vaitukaitis JL, Ross GT, Braunstein GT, Rayford PL. Gonadotropins and their subunits:basic and clinical studies *Recent Prog Horm Res* 1976;32:289-331.
- Braunstein GD, Vaitukaitis JL, Carbone PP, Ross GT. Ectopic production of human chorionic gonadotropin by neoplasm. *Ann Intern Med* 1973;78:39-45.
- Stenman UH, Alfthan H, Hotakainen K. Human chorionic gonadotropin in cancer. *Clin Biochem* 2004;37:549-61.
- Stenman UH, Alfthan H, Ranta T, et al. Serum levels of human chorionic gonadotropin in nonpregnant women and men are modulated by gonadotropin-releasing hormone and sex steroids. *Clin Endocrinol Metab* 1987;64:730-6.
- Birken S, Maydelman Y, Gawinowicz MA, et al. Isolation and characterization of human pituitary chorionic gonadotropin. *Endocrinology* 1996;137:1402-11.
- Kadar N, Caldwell BV, Romero R. A method of screening for ectopic pregnancy and its indications. *Obstet Gynecol* 1981;58:162-6.
- Cole LA. Phantom hCG and phantom choriocarcinoma. *Gynecol Oncol* 1998;71:325-9.
- Levinson SS. The nature of heterophilic antibodies and their role in immunoassay interference. *J Clin Immunoassay* 1992;15:108-115.
- Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45:942-956.
- Borkowski A, Muquardt C. Human chorionic gonadotropin in the plasma of normal, nonpregnant subjects. *N Engl J Med* 1979;301:298-302.
- Lenton EA, Neal LM, Sulaiman R. Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy. *Fertil Steril* 1982;37:773-8.
- Steier JA, Bergsjo P, Myking OL. Human chorionic gonadotropin in maternal plasma after induced abortion, spontaneous abortion and removed ectopic pregnancy. *Obstet Gynecol* 1984;64:391-4.
- Snyder JA, Haymond S, Parvin CA, et al. Diagnostic consideration in the measurement of human chorionic gonadotropin in aging women. *Clin Chem* 2005;51:1830-5.
- Wu, HBA. Tietz clinical guide to laboratory tests, 4th Ed. Philadelphia, 2006, Saunders-Elsevier, pp. 252-259.
- Tholen DW, Linnet K, Kondratovich M, et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation;

- Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.
- 18. McEnroe, RJ, Burnitt, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
 - 19. Tholen DW Kallner A, Kennedy JW et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline – Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004.

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Coeficiente de correlación, $r = 0,99$

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1. Braunstein GD, Rasor J, Adler D, et al. Serum human chorionic gonadotropin levels throughout normal pregnancy. Am J Obstet Gynecol 1976;126:678-81.
2. Brody S, Carlstrom G. Immuno-assay of human chorionic gonadotropin in normal and pathologic pregnancy. J Clin Endocrinol Metab 1962;22:564-74.
3. Vaitukaitis JL, Ross GT, Braunstein GT, Rayford PL. Gonadotropins and their subunits:basic and clinical studies Recent Prog Horm Res 1976;32:289-331.
4. Braunstein GD, Vaitukaitis JL, Carbone PP, Ross GT. Ectopic production of human chorionic gonadotropin by neoplasm. Ann Intern Med 1973;78:39-45.
5. Stenman UH, Alfthan H, Hotakainen K. Human chorionic gonadotropin in cancer. Clin Biochem 2004;37:549-61.
6. Stenman UH, Alfthan H, Ranta T, et al. Serum levels of human chorionic gonadotropin in nonpregnant women and men are modulated by gonadotropin-releasing hormone and sex steroids. Clin Endocrinol Metab 1987;64:730-6.
7. Birken S, Maydelman Y, Gawinowicz MA, et al. Isolation and characterization of human pituitary chorionic gonadotropin. Endocrinology 1996;137:1402-11.
8. Kadar N, Caldwell BV, Romero R. A method of screening for ectopic pregnancy and its indications. Obstet Gynecol 1981;58:162-6.
9. Cole LA. Phantom hCG and phantom choriocarcinoma. Gynecol Oncol 1998;71:325-9.
10. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay interference. J Clin Immunoassay 1992;15:108-115.
11. Kricka LJ. Human anti-animal antibody interference in immunological assays. Clin Chem 1999;45:942-956.
12. Borkowski A, Muquardt C. Human chorionic gonadotropin in the plasma of normal, nonpregnant subjects. N Engl J Med 1979;301:298-302.
13. Lenton EA, Neal LM, Sulaiman R. Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy. Fertil Steril 1982;37:773-8.

14. Steier JA, Bergsjo P, Myking OL. Human chorionic gonadotropin in maternal plasma after induced abortion, spontaneous abortion and removed ectopic pregnancy. *Obstet Gynecol* 1984;64:391-4.
15. Snyder JA, Haymond S, Parvin CA, et al. Diagnostic consideration in the measurement of human chorionic gonadotropin in aging women. *Clin Chem* 2005;51:1830-5.
16. Wu, HBA. Tietz clinical guide to laboratory tests, 4th Ed. Philadelphia, 2006, Saunders-Elsevier, pp. 252-259.
17. Tholen DW, Linnet K, Kondratovich M, et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.
18. McEnroe, RJ, Burnitt, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
19. Tholen DW, Kallner A, Kennedy JW et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline – Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004.

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1. Braunstein GD, Rasor J, Adler D, et al. Serum human chorionic gonadotropin levels throughout normal pregnancy. Am J Obstet Gynecol 1976;126:678-81.
2. Brody S, Carlstrom G. Immuno-assay of human chorionic gonadotropin in normal and pathologic pregnancy. J Clin Endocrinol Metab 1962;22:564-74.
3. Vaitukaitis JL, Ross GT, Braunstein GT, Rayford PL. Gonadotropins and their subunits:basic and clinical studies Recent Prog Horm Res 1976;32:289-331.
4. Braunstein GD, Vaitukaitis JL, Carbone PP, Ross GT. Ectopic production of human chorionic gonadotropin by neoplasm. Ann Intern Med 1973;78:39-45.
5. Stenman UH, Alfthan H, Hotakainen K. Human chorionic gonadotropin in cancer. Clin Biochem 2004;37:549-61.
6. Stenman UH, Alfthan H, Ranta T, et al. Serum levels of human chorionic gonadotropin in nonpregnant women and men are modulated by gonadotropin-releasing hormone and sex steroids. Clin Endocrinol Metab 1987;64:730-6.
7. Birken S, Maydelman Y, Gawinowicz MA, et al. Isolation and characterization of human pituitary chorionic gonadotropin. Endocrinology 1996;137:1402-11.
8. Kadar N, Caldwell BV, Romero R. A method of screening for ectopic pregnancy and its indications. Obstet Gynecol 1981;58:162-6.
9. Cole LA. Phantom hCG and phantom choriocarcinoma. Gynecol Oncol 1998;71:325-9.
10. Levinson SS. The nature of heterophilic antibodies and their role of in immunoassay interference. J Clin Immunoassay 1992;15:108-115.
11. Kricka LJ. Human anti-animal antibody interference in immunological assays. Clin Chem 1999;45:942-956.
12. Borkowski A, Muquardt C. Human chorionic gonadotropin in the plasma of normal, nonpregnant subjects. N Engl J Med 1979;301:298-302.
13. Lenton EA, Neal LM, Sulaiman R. Plasma concentrations of human chorionic gonadotropin from the time of implantation until the second week of pregnancy. Fertil Steril 1982;37:773-8.
14. Steier JA, Bergsjo P, Myking OL. Human chorionic gonadotropin in maternal plasma after induced abortion, spontaneous abortion and removed ectopic pregnancy. Obstet Gynecol 1984;64:391-4.

15. Snyder JA, Haymond S, Parvin CA, et al. Diagnostic consideration in the measurement of human chorionic gonadotropin in aging women. *Clin Chem* 2005;51:1830-5.
16. Wu, HBA. Tietz clinical guide to laboratory tests, 4th Ed. Philadelphia, 2006, Saunders-Elsevier, pp. 252-259.
17. Tholen DW, Linnet K, Kondratovich M, et al. NCCLS. Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline, NCCLS document EP17-A, Volume 24 Number 35, 2004.
18. McEnroe, RJ, Burnitt, MF, Powers, DM, et al. CLSI Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
19. Tholen DW Kallner A, Kennedy JW et al. NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline – Second Edition. NCCLS document EP5-A2, Volume 25 Number 27, 2004.