English

IB10 sphingotest® bio-ADM®

For the Quantitative Determination of
Bioactive Adrenomedullin in EDTA Whole Blood and Plasma

Explanation of Symbols

	CE Mark of Conformity
	Manufacturer
	Catalog Number
	Expiry date/Use by
	Lot number
	<i>In Vitro</i> Diagnostic medical device
	Consult instructions for use
	Store between 2 °C and 8 °C
	European Authorized Representative
	Contains sufficient for <n> tests
	Do not reuse
	Serial Number
	Insert disc label side up

Nexus Dx, Inc.,
6759 Mesa Ridge Road
San Diego, California 92121 USA
Telephone: 1 (858) 410 4600
Fax: 1 (858) 410 4700

TheraGenesis GmbH
Bahnhofstrasse 5
55276 Oppenheim, Germany
Tel: +49 (0) 151 506 403 14

022-00076 REV D

IB10 sphingotest® bio-ADM®

For *In Vitro* Diagnostic Use

INTENDED USE

The IB10 sphingotest® bio-ADM® is a rapid point-of-care (POC) immunoassay for the *in vitro* quantitative determination of human amidated adrenomedullin peptide (1-52), in the following referred as bioactive adrenomedullin (bio-ADM®), in human EDTA whole blood and plasma. This Test is intended for use in conjunction with the Nexus IB10 Analyzer and provides quantitative results within 20 minutes.

The IB10 sphingotest® bio-ADM® is designed for professional use only and may be used in hospital central laboratories and in alternate care settings such as emergency departments, critical care units, and other sites where near patient testing is practiced.

SUMMARY AND EXPLANATION OF THE TEST

The IB10 sphingotest® bio-ADM® quantitatively measures the level of bioactive adrenomedullin (bio-ADM®) circulating in the bloodstream. Measuring bio-ADM® using the Nexus IB10 enables the assessment of the endothelial barrier up to 48 hours before the symptoms become visible. Regular assessment of the bio-ADM® levels allows for the monitoring of critically ill patients². Elevated bio-ADM® blood levels predict both blood pressure drop resulting in shock as well as leaky vessels leading to the formation of edema³. Decreasing levels of bio-ADM® reflect an improvement of the endothelial function, which is closely associated with the patient's clinical condition. Using the IB10 sphingotest® bio-ADM® assay, clinicians can monitor levels of bio-ADM® to identify treatment responders and non-responders in order to improve patient management, including discharge decisions. Clinicians should use the IB10 sphingotest® bio-ADM® results in conjunction with the patient's other laboratory findings and clinical signs, and interpret the concrete patient values in the context of the patient's clinical situation.

BIOACTIVE ADRENOMEDULLIN

In critical care settings, one of the main causes of organ failure and ultimately mortality is loss of endothelial function, which is associated with leakage of blood vessels. Although the symptoms of loss of endothelial functions are well-known, there are currently no simple, blood-based detection methods established for monitoring the worsening and improvement of endothelial function. bio-ADM® has been identified as a controlling hormone of the endothelial barrier, the interior wall protecting the blood vessels¹. In certain conditions such as septic shock⁵⁻⁶, acute heart failure⁷⁻⁹ or cardiogenic shock¹⁰, the endothelial barrier becomes leaky, and additional bio-ADM® is produced to re-seal the barrier. However, bio-ADM® has a second function. It also vasodilates the blood vessels, resulting in a dangerous blood pressure drop, which leads to shock and may ultimately escalate into multiple organ failure.

PRINCIPLE

The Nexus IB10 immunochemistry system combines chemistry with microfluidics and centrifugal flow to rapidly prepare a cell free plasma from whole blood that can then be moved through a channel to rehydrate, solubilize and mix with freeze dried immunoconjugates. Using a combination of active flow and capillary action, the test sample is quantitatively measured within 20 minutes with an optical signal level proportional to the analyte(s) concentration.

After addition of the patient sample, the entire test is performed within the Nexus IB10 Analyzer, which provides control of the temperature of the disc, as well as the reaction sequence, centrifugal flow, mixing, incubation time, final signal measurement, quantitation and reporting of results. The Test disc includes a positive internal control to ensure that the Test has operated properly. Each lot is calibrated to provide accurate analyte concentration. Lot specific calibration along with additional information such as the lot expiration date is contained on a QR code label affixed to each disc.

REAGENTS

The IB10 sphingotest® bio-ADM® contains all required reagents to evaluate the level of bioactive adrenomedullin including dye-conjugated monoclonal anti-bioactive adrenomedullin antibody, biotin-conjugated monoclonal anti-bioactive adrenomedullin antibody and streptavidin immobilized at the detection area on the disc.

MATERIALS PROVIDED

Each box contains the following:

- 10 IB10 sphingotest® bio-ADM® discs, each individually sealed in a foil pouch with a desiccant.
- Instructions for Use (IFU).

MATERIALS/EQUIPMENT REQUIRED BUT NOT PROVIDED

1. Nexus IB10 Analyzer - Model #BCA-IB10.
2. Calibrated reusable fixed or variable volume pipette gun with high precision and accuracy capable of delivering 500 µL of whole blood or plasma.
3. Disposable pipette tips capable of accepting and delivering 500 µL of whole blood or plasma.

PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Carefully follow the Instructions for Use.
- Prior to testing patient samples, ensure that the Analyzer software is updated to the latest version. (Refer to the Nexus IB10 manual for specific instructions).
- Wear disposable gloves while handling samples.
- Handle samples with care. Samples and used test discs should be treated as potentially infectious and should be discarded as biohazardous material according to local regulations.
- Thoroughly wash hands following handling.
- The result obtained from the IB10 sphingotest® bio-ADM® does not provide a definitive diagnosis and should be interpreted by a physician in conjunction with other laboratory test results and according to current medical guidelines and patient clinical findings.
- Keep the test disc in the sealed pouch until ready for use.
- Do not use the test disc if the pouch is damaged or the seal is broken.
- Do not use the test disc after the expiration date printed on the pouch.
- Prior to use, place the unopened pouch at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
- Always pay attention to cleanliness when handling the test disc. Avoid any contamination from fingerprints or foreign substances. Do not contaminate the sample channel inlet.
- Do not drop or damage the test disc.
- The test disc should be inserted with the label side up, into the Nexus IB10 Analyzer tray immediately after injecting the sample into the disc.
- Do not turn the disc upside down.
- This is a quantitative test; therefore no visual interpretation of the results should be made.

STORAGE AND STABILITY

- Store the IB10 sphingotest® bio-ADM® disc between 2 and 8 °C (35 and 46 °F) until the expiration date printed on the pouch is reached.
- The IB10 sphingotest® bio-ADM® in its sealed pouch is stable between 18 and 30 °C (64 and 86 °F) for 14 days, provided the expiration date printed on the pouch is not exceeded.

SAMPLE COLLECTION AND PREPARATION

- The IB10 sphingotest® bio-ADM® is to be run using EDTA whole blood or plasma samples.
- It is recommended that samples be tested as soon as possible after collection.
- Plasma samples should be kept frozen at -20 °C (- 4 °F) or lower if longer storage is required.
- Allow samples to equilibrate to room temperature (19 to 25 °C/66 to 77 °F) prior to testing.

PROCEDURE

Nexus IB10 Analyzer

 Consult the Nexus IB10 Analyzer User Manual

For Analyzer installation, start up and complete instructions for use, refer to the **Nexus IB10 User Manual**. Operator must consult the User Manual prior to use to become familiar with the proper operation and quality control procedures.

PERFORMING SYSTEM CHECK AND DISC CALIBRATION

Each time the Nexus IB10 Analyzer is turned on, a Self Check is automatically performed. The QR code on each test disc contains information for disc calibration which the Analyzer automatically reads when running a test.

TESTING PATIENT SAMPLES ON THE NEXUS IB10 ANALYZER

1. Remove an unopened Test pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
2. Open the pouch and remove the Test disc.
3. Place the Test disc on a level surface.
4. On the Nexus IB10 Analyzer press **New Analysis**.
5. The Analyzer will perform a general system check.
6. Enter the Patient ID manually (up to 20 characters can be used for the ID) or enter the Patient ID by using the barcode scanner.
7. Mix the whole blood patient sample tube by gently inverting the tube several times before testing.

8. Testing Patient Samples on the IB10 sphingotest® bio-ADM®

- Using a precision pipette (fixed or adjusted to 500 µL), slowly draw well-mixed patient sample into the pipette tip.
- Positioning the tapered pipette tip at a 45° angle, pierce the X on the red dot to expose the sample channel inlet.
- Slowly express the patient sample into the inlet applying minimum, but continuous force, on the pipette plunger.
- Express sample to the first stop on the pipette, at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
- Press **OK** on the Nexus IB10 Analyzer display.
- When the tray opens, insert the filled Test disc into the tray and press **Run**.
- Within 20 minutes, the Nexus IB10 Analyzer will display the results on the screen.
- Results will print automatically (if selected during Set Up) or press **Print**.
- Remove the test disc and discard in appropriate receptacle.

Note: If the test run is cancelled before a test result is displayed, the test disc cannot be reused and should be disposed of appropriately.

INTERPRETATION OF RESULTS

The range of Bioactive Adrenomedullin concentrations reported by the Nexus IB10 is 45 pg/mL to 500 pg/mL. Results below or above this range will be shown as “<45 pg/mL” or “>500 pg/mL”, respectively.

Quality Control

INTERNAL CONTROL

The IB10 sphingotest® bio-ADM® has a built in positive procedural control. The Nexus IB10 Analyzer automatically determines the presence of this control thereby confirming that the test run has delivered a valid result. If the control does not form or if it is not recognized by the Analyzer, the test result is considered “invalid” and the test must be repeated.

LIMITATIONS

The test result should be evaluated and interpreted in conjunction with the patient’s medical history, symptoms and other clinical information.

As with any assay employing mouse antibodies, the possibility exists for interference by heterophilic antibodies in the sample.^{12,13} Specimens from individuals who have been regularly exposed to animals or who have been treated with animal products may contain these antibodies.

Other substances and/or factors not listed, e.g. technical or procedural error, may interfere with the test and cause inaccurate results.

Nexus Dx, Inc. offers products for their intended use. Refer to the specific product literature for the intended use statements for each product. Product claims are subject to change. Nexus Dx, Inc.’s expressed and implied warranties (inclusive of implied warranties of merchantability and fitness) are conditional upon adherence to, or observance of Nexus Dx, Inc. published directions with regard to the use of Nexus Dx, Inc. products. Under no circumstances will Nexus Dx, Inc. be liable for any indirect or consequential damages.

For Technical Assistance please contact your local Distributor.

Performance Characteristics

MEASURING RANGE

The IB10 sphingotest® bio-ADM® has been demonstrated to provide measurable results at Bioactive Adrenomedullin levels from 45 pg/mL to 500 pg/mL.

ANALYTICAL SENSITIVITY

The LoD (Limit of Detection) of the IB10 sphingotest® bio-ADM® is 45 pg/mL, determined according to Clinical and Laboratory Standards Institute (CLSI) Approved Guideline EP17-A2.¹⁴ LoQ (Limit of Quantitation) is the lowest Bioactive Adrenomedullin concentration that can be reproducibly measured with a total coefficient of variation of at most 20%. It was determined to be 45 pg/mL.

INTERFERING SUBSTANCES AND CROSS-REACTIVITY

The following substances were tested for potential interference with Bioactive Adrenomedullin measurement in the IB10 sphingotest® bio-ADM® based on CLSI Approved Guideline EP7-A2 and EP37 (Table 1).¹⁵⁻¹⁶ No significant interference was observed at the maximum concentration for substances listed in the table below.

Table 1.

Potential Interferents	Maximum Concentration
Bilirubin (Conjugated)	8 mg/dL
Bilirubin (Free)	20 mg/dL
Biotin	1200 ng/mL
Cholesterol	400 mg/dL
Citrate	109 mmol/L
Creatinine	15 mg/dL
Glucose	1000 mg/dL
Hemoglobin	500 mg/dL
Unfractionated Heparin	330 U/dL
Fractionated Heparin	3000 U/L
Human anti-mouse antibody (HAMA)	711 ng/mL
Rheumatoid factor (RF)	1500 IU/mL
Triglycerides	750 mg/dL
Urea	120 mg/dL

HOOK EFFECT

No high dose hook effect is observed for Bioactive Adrenomedullin concentrations up to 100,000 pg/mL.

PRECISION

The precision of Bioactive Adrenomedullin measurement of the IB10 sphingotest® bio-ADM® was determined using samples where Bioactive Adrenomedullin was added to normal human plasma at two concentrations (Table 2). The within-run and total precision was performed in two runs per day, in replicates of 5 per run at each concentration level for a 10 day period for a total number of repetitions of 100 at each concentration level. The within-run and total precision were computed according to the CLSI Approved Guideline EP05-A3.¹⁷

Table 2.

Sample	Test Lot	Mean (pg/mL)	Within-Run Precision		Total Precision	
			Std. dev. (pg/mL)	CV (%)	Std. dev. (pg/mL)	CV (%)
1	A	58.5	8.6	14.7%	8.8	15.0%
	B	53.5	8.9	16.7%	9.0	16.8%
2	A	149.4	10.9	7.3%	11.9	8.0%
	B	150.6	15.7	10.4%	16.1	10.7%

WHOLE BLOOD Vs. PLASMA COMPARISON

A comparison study was performed using matching EDTA whole blood and EDTA plasma samples. When performing a Passing-Bablok regression analysis comparing the whole blood concentrations versus the corresponding plasma concentrations from the same subject samples (N=43), a slope of 1.03 (95% C.I. = [0.96-1.10]), an intercept of -7.0 pg/mL, and a correlation coefficient of $r = 0.99$ were obtained.

METHOD COMPARISON

An equivalence study was performed between the IB10 sphingotest® bio-ADM® using the Nexus IB10 Analyzer, and the sphingotest® bio-ADM® immunoassay (ILMA) assay. A total of 102 samples were tested within a concentration range of Bioactive Adrenomedullin of 42.4 pg/mL to 553.6 pg/mL. The Passing-Bablok regression was:

$$\text{IB10 sphingotest® bio-ADM®} = 0.96 \times (\text{sphingotest® bio-ADM® ILMA}) + 3.6 \text{ pg/mL}$$

Correlation coefficient, rho = 0.90

REFERENCES

1. Geven C et al. Vascular effects of adrenomedullin and the anti-adrenomedullin antibody Adrecizumab in sepsis; Shock. 2018 Aug;50(2):132-140. doi: 10.1097
2. Mebazaa A et al. Circulating adrenomedullin estimates survival and reversibility of organ failure in sepsis: the prospective observational multinational Adrenomedullin and Outcome in Sepsis and Septic Shock-1 (AdrenOSS-1) study; Crit Care. 2018 Dec 21;22(1):354. doi: 10.1186/s13054-018-2243-2
3. Caiaroni P et al. Circulating biologically active adrenomedullin (bio-ADM) predicts hemodynamic support requirement and mortality during sepsis; Chest. 2017 Aug;152(2):312-320. doi: 10.1016/j.chest.2017.03.035
4. Marino R et al. Plasma adrenomedullin is associated with short-term mortality and vasopressor requirement in patients admitted with sepsis; Crit Care. 2014 Feb 17;18(1):R34. doi: 10.1186/cc13731
5. Gayat E et al. Determinants of long-term outcome in ICU survivors results from the FROG-ICU study; Crit Care. 2018 Jan 18;22(1):8. doi: 10.1186/s13054-017-1922-8
6. Kim H et al. Circulating biologically active adrenomedullin predicts organ failure and mortality in sepsis; Ann Lab Med. 2019 Sep;39(5):454-463. doi: 10.3343/alm.2019.39.5.454
7. Pandhi P et al. Clinical value of pre-discharge bio-adrenomedullin as a marker of residual congestion and high risk of heart failure hospital readmission; Eur J Heart Fail. 2020 Apr;22(4):683-691. doi: 10.1002/ejhf.1693
8. Ter Maaten JM et al. Bio-adrenomedullin as a marker of congestion in patients with new-onset and worsening heart failure; Eur J Heart Fail. 2019 Jun;21(6):732-743. doi: 10.1002/ejhf.1437
9. Voors AA et al. Adrenomedullin in heart failure pathophysiology and therapeutic application; Eur J Heart Fail. 2019 Feb;21(2):163-171. doi: 10.1002/ejhf.1366
10. Tolppanen H et al. Adrenomedullin - a marker of impaired hemodynamics, organ dysfunction, and poor prognosis in cardiogenic shock; Ann Intensive Care. 2017 Dec;7(1):6. doi: 10.1186/s13613-016-0229-2
11. Weber et al. Sandwich immunoassay for bioactive plasma adrenomedullin; J Applied Laboratory Medicine. 2017 Sep; 2(2):222-223. doi: 10.1373/jalm.2017.023655
12. Kricka LJ. Human anti-animal antibody interference in immunological assays. Clin Chem 1999;45: 942-6.

13. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15: 108-15.
14. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
15. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
16. McEnroe RJ, Dimeski G, Miller WG, et al. Supplemental Tables for Interference Testing in Clinical Chemistry. CLSI Document EP37, First Edition, 2018.
17. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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Korrelationskoeffizient, $p = 0,90$

REFERENZEN

1. Geven C et al. Vascular effects of adrenomedullin and the anti-adrenomedullin antibody Adrecizumab in sepsis; Shock. 2018 Aug;50(2):132-140. doi: 10,1097
2. Mebazaa A et al. Circulating adrenomedullin estimates survival and reversibility of organ failure in sepsis: the prospective observational

- multinational Adrenomedullin and Outcome in Sepsis and Septic Shock-1 (AdrenOSS-1) study; Crit Care. 2018 Dec 21;22(1):354. doi: 10.1186/s13054-018-2243-2
- 3. Caironi P et al. Circulating biologically active adrenomedullin (bio-ADM) predicts hemodynamic support requirement and mortality during sepsis; Chest. 2017 Aug;152(2):312-320. doi: 10.1016/j.chest.2017.03.035
 - 4. Marino R et al. Plasma adrenomedullin is associated with short-term mortality and vasopressor requirement in patients admitted with sepsis; Crit Care. 2014 Feb 17;18(1):R34. doi: 10.1186/cc13731
 - 5. Gayat E et al. Determinants of long-term outcome in ICU survivors results from the FROG-ICU study; Crit Care. 2018 Jan 18;22(1):8. doi: 10.1186/s13054-017-1922-8
 - 6. Kim H et al. Circulating biologically active adrenomedullin predicts organ failure and mortality in sepsis; Ann Lab Med. 2019 Sep;39(5):454-463. doi: 10.3343/alm.2019.39.5.454
 - 7. Pandhi P et al. Clinical value of pre-discharge bio-adrenomedullin as a marker of residual congestion and high risk of heart failure hospital readmission; Eur J Heart Fail. 2020 Apr;22(4):683-691. doi: 10.1002/ejhf.1693
 - 8. Ter Maaten JM et al. Bio-adrenomedullin as a marker of congestion in patients with new-onset and worsening heart failure; Eur J Heart Fail. 2019 Jun;21(6):732-743. doi: 10.1002/ejhf.1437
 - 9. Voors AA et al. Adrenomedullin in heart failure pathophysiology and therapeutic application; Eur J Heart Fail. 2019 Feb;21(2):163-171. doi: 10.1002/ejhf.1366
 - 10. Tolppanen H et al. Adrenomedullin - a marker of impaired hemodynamics, organ dysfunction, and poor prognosis in cardiogenic shock; Ann Intensive Care. 2017 Dec;7(1):6. doi: 10.1186/s13613-016-0229-2
 - 11. Weber et al. Sandwich immunoassay for bioactive plasma adrenomedullin; J Applied Laboratory Medicine. 2017 Sep; 2(2):222-223. doi: 10.1373/jalm.2017.023655
 - 12. Kricka LJ. Human anti-animal antibody interference in immunological assays. Clin Chem 1999;45: 942-6.
 - 13. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. J Clin Immunoassay 1992;15: 108-15.
 - 14. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
 - 15. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
 - 16. McEnroe RJ, Dimeski G, Miller WG, et al. Supplemental Tables for Interference Testing in Clinical Chemistry. CLSI Document EP37, First Edition, 2018.
 - 17. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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Coefficient de corrélation, rho = 0,90

BIBLIOGRAPHIE

1. Geven C et al. Vascular effects of adrenomedullin and the anti-adrenomedullin antibody Adrecizumab in sepsis; Shock. 2018 Aug;50(2):132-140. doi: 10.1097
2. Mebazaa A et al. Circulating adrenomedullin estimates survival and reversibility of organ failure in sepsis: the prospective observational multinational Adrenomedullin and Outcome in Sepsis and Septic Shock-1 (AdrenOSS-1) study; Crit Care. 2018 Dec 21;22(1):354. doi: 10.1186/s13054-018-2243-2
3. Caironi P et al. Circulating biologically active adrenomedullin (bio-ADM) predicts hemodynamic support requirement and mortality during sepsis; Chest. 2017 Aug;152(2):312-320. doi: 10.1016/j.chest.2017.03.035
4. Marino R et al. Plasma adrenomedullin is associated with short-term mortality and vasopressor requirement in patients admitted with

- sepsis; Crit Care. 2014 Feb 17;18(1):R34. doi: 10.1186/cc13731
- 5. Gayat E et al. Determinants of long-term outcome in ICU survivors results from the FROG-ICU study; Crit Care. 2018 Jan 18;22(1):8. doi: 10.1186/s13054-017-1922-8
 - 6. Kim H et al. Circulating biologically active adrenomedullin predicts organ failure and mortality in sepsis; Ann Lab Med. 2019 Sep;39(5):454-463. doi: 10.3343/alm.2019.39.5.454
 - 7. Pandhi P et al. Clinical value of pre-discharge bio-adrenomedullin as a marker of residual congestion and high risk of heart failure hospital readmission; Eur J Heart Fail. 2020 Apr;22(4):683-691. doi: 10.1002/ejhf.1693
 - 8. Ter Maaten JM et al. Bio-adrenomedullin as a marker of congestion in patients with new-onset and worsening heart failure; Eur J Heart Fail. 2019 Jun;21(6):732-743. doi: 10.1002/ejhf.1437
 - 9. Voors AA et al. Adrenomedullin in heart failure pathophysiology and therapeutic application; Eur J Heart Fail. 2019 Feb;21(2):163-171. doi: 10.1002/ejhf.1366
 - 10. Tolppanen H et al. Adrenomedullin - a marker of impaired hemodynamics, organ dysfunction, and poor prognosis in cardiogenic shock; Ann Intensive Care. 2017 Dec;7(1):6. doi: 10.1186/s13613-016-0229-2
 - 11. Weber et al. Sandwich immunoassay for bioactive plasma adrenomedullin; J Applied Laboratory Medicine. 2017 Sep; 2(2):222-223. doi: 10.1373/jalm.2017.023655
 - 12. Kricka LJ. Human anti-animal antibody interference in immunological assays. Clin Chem 1999;45: 942-6.
 - 13. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. J Clin Immunoassay 1992;15: 108-15.
 - 14. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
 - 15. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
 - 16. McEnroe RJ, Dimeski G, Miller WG, et al. Supplemental Tables for Interference Testing in Clinical Chemistry. CLSI Document EP37, First Edition, 2018.
 - 17. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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BIBLIOGRAFIA

1. Geven C et al. Vascular effects of adrenomedullin and the anti-adrenomedullin antibody Adrecizumab in sepsis; Shock. 2018 Aug;50(2):132-140. doi: 10.1097
2. Mebazaa A et al. Circulating adrenomedullin estimates survival and reversibility of organ failure in sepsis: the prospective observational multinational Adrenomedullin and Outcome in Sepsis and Septic Shock-1 (AdrenOSS-1) study; Crit Care. 2018 Dec 21;22(1):354. doi: 10.1186/s13054-018-2243-2
3. Caironi P et al. Circulating biologically active adrenomedullin (bio-ADM) predicts hemodynamic support requirement and mortality during sepsis; Chest. 2017 Aug;152(2):312-320. doi: 10.1016/j.chest.2017.03.035
4. Marino R et al. Plasma adrenomedullin is associated with short-term mortality and vasopressor requirement in patients admitted with sepsis; Crit Care. 2014 Feb 17;18(1):R34. doi: 10.1186/cc13731
5. Gayat E et al. Determinants of long-term outcome in ICU survivors results from the FROG-ICU study; Crit Care. 2018 Jan 18;22(1):8. doi: 10.1186/s13054-017-1922-8
6. Kim H et al. Circulating biologically active adrenomedullin predicts organ failure and mortality in sepsis; Ann Lab Med. 2019 Sep;39(5):454-463. doi: 10.3343/alm.2019.39.5.454
7. Pandhi P et al. Clinical value of pre-discharge bio-adrenomedullin as a marker of residual congestion and high risk of heart failure hospital

- readmission; Eur J Heart Fail. 2020 Apr;22(4):683-691. doi: 10.1002/ejhf.1693
- 8. Ter Maaten JM et al. Bio-adrenomedullin as a marker of congestion in patients with new-onset and worsening heart failure; Eur J Heart Fail. 2019 Jun;21(6):732-743. doi: 10.1002/ejhf.1437
 - 9. Voors AA et al. Adrenomedullin in heart failure pathophysiology and therapeutic application; Eur J Heart Fail. 2019 Feb;21(2):163-171. doi: 10.1002/ejhf.1366
 - 10. Tolppanen H et al. Adrenomedullin - a marker of impaired hemodynamics, organ dysfunction, and poor prognosis in cardiogenic shock; Ann Intensive Care. 2017 Dec;7(1):6. doi: 10.1186/s13613-016-0229-2
 - 11. Weber et al. Sandwich immunoassay for bioactive plasma adrenomedullin; J Applied Laboratory Medicine. 2017 Sep; 2(2):222-223. doi: 10.1373/jalm.2017.023655
 - 12. Kricka LJ. Human anti-animal antibody interference in immunological assays. Clin Chem 1999;45: 942-6.
 - 13. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. J Clin Immunoassay 1992;15: 108-15.
 - 14. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
 - 15. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
 - 16. McEnroe RJ, Dimeski G, Miller WG, et al. Supplemental Tables for Interference Testing in Clinical Chemistry. CLSI Document EP37, First Edition, 2018.
 - 17. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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REFERÊNCIAS

1. Geven C et al. Vascular effects of adrenomedullin and the anti-adrenomedullin antibody Adrecizumab in sepsis; Shock. 2018 Aug;50(2):132-140. doi: 10.1097
2. Mebazaa A et al. Circulating adrenomedullin estimates survival and reversibility of organ failure in sepsis: the prospective observational multinational Adrenomedullin and Outcome in Sepsis and Septic Shock-1 (AdrenOSS-1) study; Crit Care. 2018 Dec 21;22(1):354. doi: 10.1186/s13054-018-2243-2
3. Caironi P et al. Circulating biologically active adrenomedullin (bio-ADM) predicts hemodynamic support requirement and mortality during sepsis; Chest. 2017 Aug;152(2):312-320. doi: 10.1016/j.chest.2017.03.035
4. Marino R et al. Plasma adrenomedullin is associated with short-term mortality and vasopressor requirement in patients admitted with sepsis; Crit Care. 2014 Feb 17;18(1):R34. doi: 10.1186/cc13731
5. Gayat E et al. Determinants of long-term outcome in ICU survivors results from the FROG-ICU study; Crit Care. 2018 Jan 18;22(1):8. doi: 10.1186/s13054-017-1922-8
6. Kim H et al. Circulating biologically active adrenomedullin predicts organ failure and mortality in sepsis; Ann Lab Med. 2019 Sep;39(5):454-463. doi: 10.3343/alm.2019.39.5.454

7. Pandhi P et al. Clinical value of pre-discharge bio-adrenomedullin as a marker of residual congestion and high risk of heart failure hospital readmission; Eur J Heart Fail. 2020 Apr;22(4):683-691. doi: 10.1002/ejhf.1693
8. Ter Maaten JM et al. Bio-adrenomedullin as a marker of congestion in patients with new-onset and worsening heart failure; Eur J Heart Fail. 2019 Jun;21(6):732-743. doi: 10.1002/ejhf.1437
9. Voors AA et al. Adrenomedullin in heart failure pathophysiology and therapeutic application; Eur J Heart Fail. 2019 Feb;21(2):163-171. doi: 10.1002/ejhf.1366
10. Tolppanen H et al. Adrenomedullin - a marker of impaired hemodynamics, organ dysfunction, and poor prognosis in cardiogenic shock; Ann Intensive Care. 2017 Dec;7(1):6. doi: 10.1186/s13613-016-0229-2
11. Weber et al. Sandwich immunoassay for bioactive plasma adrenomedullin; J Applied Laboratory Medicine. 2017 Sep; 2(2):222-223. doi: 10.1373/jalm.2017.023655
12. Kricka LJ. Human anti-animal antibody interference in immunological assays. Clin Chem 1999;45: 942-6.
13. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. J Clin Immunoassay 1992;15: 108-15.
14. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
15. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
16. McEnroe RJ, Dimeski G, Miller WG, et al. Supplemental Tables for Interference Testing in Clinical Chemistry. CLSI Document EP37, First Edition, 2018.
17. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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Coeficiente de correlación, rho = 0,90

REFERENCIAS

1. Geven C et al. Vascular effects of adrenomedullin and the anti-adrenomedullin antibody Adrecizumab in sepsis; Shock. 2018 Aug;50(2):132-140. doi: 10,1097
2. Mebazaa A et al. Circulating adrenomedullin estimates survival and reversibility of organ failure in sepsis: the prospective observational multinational Adrenomedullin and Outcome in Sepsis and Septic Shock-1 (AdrenOSS-1) study; Crit Care. 2018 Dec 21;22(1):354. doi: 10.1186/s13054-018-2243-2
3. Caironi P et al. Circulating biologically active adrenomedullin (bio-ADM) predicts hemodynamic support requirement and mortality during sepsis; Chest. 2017 Aug;152(2):312-320. doi: 10.1016/j.chest.2017.03.035
4. Marino R et al. Plasma adrenomedullin is associated with short-term mortality and vasopressor requirement in patients admitted with sepsis; Crit Care. 2014 Feb 17;18(1):R34. doi: 10.1186/cc13731

5. Gayat E et al. Determinants of long-term outcome in ICU survivors results from the FROG-ICU study; Crit Care. 2018 Jan 18;22(1):8. doi: 10.1186/s13054-017-1922-8
6. Kim H et al. Circulating biologically active adrenomedullin predicts organ failure and mortality in sepsis; Ann Lab Med. 2019 Sep;39(5):454-463. doi: 10.3343/alm.2019.39.5.454
7. Pandhi P et al. Clinical value of pre-discharge bio-adrenomedullin as a marker of residual congestion and high risk of heart failure hospital readmission; Eur J Heart Fail. 2020 Apr;22(4):683-691. doi: 10.1002/ejhf.1693
8. Ter Maaten JM et al. Bio-adrenomedullin as a marker of congestion in patients with new-onset and worsening heart failure; Eur J Heart Fail. 2019 Jun;21(6):732-743. doi: 10.1002/ejhf.1437
9. Voors AA et al. Adrenomedullin in heart failure pathophysiology and therapeutic application; Eur J Heart Fail. 2019 Feb;21(2):163-171. doi: 10.1002/ejhf.1366
10. Tolppanen H et al. Adrenomedullin - a marker of impaired hemodynamics, organ dysfunction, and poor prognosis in cardiogenic shock; Ann Intensive Care. 2017 Dec;7(1):6. doi: 10.1186/s13613-016-0229-2
11. Weber et al. Sandwich immunoassay for bioactive plasma adrenomedullin; J Applied Laboratory Medicine. 2017 Sep; 2(2):222-223. doi: 10.1373/jalm.2017.023655
12. Kricka LJ. Human anti-animal antibody interference in immunological assays. Clin Chem 1999;45: 942-6.
13. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. J Clin Immunoassay 1992;15: 108-15.
14. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
15. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
16. McEnroe RJ, Dimeski G, Miller WG, et al. Supplemental Tables for Interference Testing in Clinical Chemistry. CLSI Document EP37, First Edition, 2018.
17. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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1. Geven C et al. Vascular effects of adrenomedullin and the anti-adrenomedullin antibody Adrecizumab in sepsis; Shock. 2018 Aug;50(2):132-140. doi: 10.1097
2. Mebazaa A et al. Circulating adrenomedullin estimates survival and reversibility of organ failure in sepsis: the prospective observational multinational Adrenomedullin and Outcome in Sepsis and Septic Shock-1 (AdrenOSS-1) study; Crit Care. 2018 Dec 21;22(1):354. doi: 10.1186/s13054-018-2243-2
3. Caironi P et al. Circulating biologically active adrenomedullin (bio-ADM) predicts hemodynamic support requirement and mortality during sepsis; Chest. 2017 Aug;152(2):312-320. doi: 10.1016/j.chest.2017.03.035
4. Marino R et al. Plasma adrenomedullin is associated with short-term mortality and vasopressor requirement in patients admitted with sepsis; Crit Care. 2014 Feb 17;18(1):R34. doi: 10.1186/cc13731
5. Gayat E et al. Determinants of long-term outcome in ICU survivors results from the FROG-ICU study; Crit Care. 2018 Jan 18;22(1):8. doi: 10.1186/s13054-017-1922-8
6. Kim H et al. Circulating biologically active adrenomedullin predicts organ failure and mortality in sepsis; Ann Lab Med. 2019 Sep;39(5):454-463. doi: 10.3343/alm.2019.39.5.454
7. Pandhi P et al. Clinical value of pre-discharge bio-adrenomedullin as a marker of residual congestion and high risk of heart failure hospital readmission; Eur J Heart Fail. 2020 Apr;22(4):683-691. doi: 10.1002/ejhf.1693
8. Ter Maaten JM et al. Bio-adrenomedullin as a marker of congestion in patients with new-onset and worsening heart failure; Eur J Heart Fail. 2019 Jun;21(6):732-743. doi: 10.1002/ejhf.1437
9. Voors AA et al. Adrenomedullin in heart failure pathophysiology and therapeutic application; Eur J Heart Fail. 2019 Feb;21(2):163-171. doi: 10.1002/ejhf.1366
10. Tolppanen H et al. Adrenomedullin - a marker of impaired hemodynamics, organ dysfunction, and poor prognosis in cardiogenic shock; Ann Intensive Care. 2017 Dec;7(1):6. doi: 10.1186/s13613-016-0229-2
11. Weber et al. Sandwich immunoassay for bioactive plasma adrenomedullin; J Applied Laboratory Medicine. 2017 Sep; 2(2):222-223. doi: 10.1373/jalm.2017.023655
12. Kricka LJ. Human anti-animal antibody interference in immunological assays. Clin Chem 1999;45: 942-6.
13. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. J Clin Immunoassay 1992;15: 108-15.
14. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
15. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
16. McEnroe RJ, Dimeski G, Miller WG, et al. Supplemental Tables for Interference Testing in Clinical Chemistry. CLSI Document EP37, First Edition, 2018.
17. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.