
English

IB10 sphingotest® PCT
 For the Quantitative Determination of
 Procalcitonin in EDTA Whole Blood and Plasma

For U.S. Export Only

Explanation of Symbols

	CE Mark of Conformity
	Manufacturer
	Catalog Number
	Expiry date/Use by
	Lot number
	<i>In Vitro</i> Diagnostic medical device
	Consult instructions for use
	Store between 2 °C and 8 °C
	European Authorized Representative
	Contains sufficient for <n> test
	Do not reuse
	Serial Number
	Insert disc label side up

Nexus Dx, Inc.,
 6759 Mesa Ridge Road
 San Diego, California 92121 USA
 Telephone: 1 (858) 410 4600
 Fax: 1 (858) 410 4700

TheraGenesis GmbH
 Bahnhofstrasse 5
 55276 Oppenheim, Germany
 Tel: +49 (0) 151 506 403 14

022-00084 REV B

IB10 sphingotest® PCT

For *In Vitro* Diagnostic Use

INTENDED USE

The IB10 sphingotest® PCT is a rapid point-of-care (POC) immunoassay for the *in vitro* quantitative determination of Procalcitonin (PCT) in the concentration range of 0.3 µg/L to 10.00 µg/L in EDTA whole blood and plasma. This PCT Test is intended for use in conjunction with the Nexus IB10 Analyzer and provides quantitative results in 20 minutes.

The IB10 sphingotest® PCT is intended to be used in conjunction with clinical evaluations and other laboratory findings to diagnose sepsis, to assess the degree of severity, and to monitor the change of sepsis severity.

This assay is neither suitable nor intended for other applications and/or uses.

The IB10 sphingotest® PCT is designed for professional use only and may be used in hospital central laboratories and in alternate care settings such as emergency departments, critical care units, and other POC sites where near patient testing is practiced.

SUMMARY AND EXPLANATION OF THE TEST

Sepsis is the most common cause of death in intensive care units (ICU) with a mortality rate up to 50% depending on severity. The earlier the sepsis is identified and treated, the better the prognosis.¹ Early diagnosis and treatment of sepsis is cost effective and reduces the number of hospital and critical care bed days for patients.

Procalcitonin (PCT), a precursor of the peptide hormone, calcitonin (CT), is produced by neuroendocrine cells successively, which are enzymatically cleaved into immature calcitonin, katalcalcin, and the N-terminal region.² During sepsis due to bacterial infection, the final step in the synthesis of calcitonin is inhibited, resulting in an elevation of procalcitonin levels in multiple tissues.³

PCT levels in sepsis are often greater than 1 µg/L and reach values of 10 µg/L or even higher in individual patients with severe sepsis and septic shock. As the septic infection resolves, the PCT levels also return to ranges < 0.5 µg/L, with a half-life of 24 hours. Consequently, the *in vitro* determination of PCT can be used to monitor the course and prognosis of clinically relevant sepsis and to guide and control therapeutic interventions.⁴⁻⁵

The IB10 sphingotest® PCT is intended for the rapid *in vitro* quantitative measurement of PCT in the concentration range of 0.3 µg/L to 10.00 µg/L in EDTA whole blood and plasma using immunochromatographic technology in conjunction with the Nexus IB10 Analyzer.

PRINCIPLE

The Nexus IB10 immunochemistry system combines chemistry with microfluidics and centrifugal flow to rapidly prepare a cell free plasma from whole blood that can then be moved through a channel to rehydrate, solubilize and mix with freeze dried immunoconjugates. Using a combination of active flow and capillary action, the test sample is quantitatively measured in 20 minutes with an optical signal level proportional to the analyte(s) concentration.

After addition of the patient sample, the entire test is performed within the Nexus IB10 Analyzer, which provides control of the temperature of the disc, as well as the reaction sequence, centrifugal flow, mixing, incubation time, final signal measurement, quantitation and reporting of results. The Test disc includes a positive internal control to ensure that the Test disc has operated properly. Each lot is calibrated to provide accurate analyte concentration. Lot specific calibration along with additional information such as the lot expiration date is contained on a QR code label affixed to each disc. It is recommended that external controls be tested at appropriate time intervals to confirm that the system and test lot are performing within acceptable limits.

REAGENTS

The IB10 sphingotest® PCT contains all required reagents to evaluate the level of PCT including dye conjugated monoclonal anti-PCT antibody, biotin conjugated monoclonal anti-PCT antibody and streptavidin immobilized at the detection area on the disc.

MATERIALS PROVIDED

Each box contains the following:

- 10 IB10 sphingotest® PCT discs, each individually sealed in a foil pouch with a desiccant.
- 1 Instructions for Use (IFU)

MATERIALS/EQUIPMENT REQUIRED BUT NOT PROVIDED

1. Nexus IB10 Analyzer - Model #BCA-IB10
2. Commercially available PCT Controls for external Quality Control (QC). **Contact the Distributor in your area** for recommended external QC materials or related technical assistance.
3. Calibrated reusable fixed or variable volume pipette gun with high precision and accuracy capable of delivering 500 µL of whole blood or plasma.
4. Disposable pipette tips capable of accepting and delivering 500 µL of whole blood or plasma.

PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Carefully follow the Instructions for Use.
- Prior to testing controls or patient samples, ensure that the Analyzer software is updated to the latest version. (Refer to the Nexus IB10 manual for specific instructions).
- Wear disposable gloves while handling samples.
- Handle samples with care. Samples and used Test discs should be treated as potentially infectious and should be discarded as biohazardous material according to local regulations.
- Thoroughly wash hands following handling.
- The result obtained from the IB10 sphingotest® PCT does not provide a definitive diagnosis and should be interpreted by a physician in conjunction with other laboratory test results and according to current medical guidelines and patient clinical findings.
- Keep the Test disc in the sealed pouch until ready for use.
- Do not use the Test disc if the pouch is damaged or the seal is broken.
- Do not use the Test disc after the expiration date printed on the pouch.
- Prior to use, place the unopened pouch at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
- Always pay attention to cleanliness when handling the test disc. Avoid any contamination from fingerprints or foreign substances. Do not contaminate the sample channel inlet.
- Do not drop or damage the test disc.
- The test disc should be inserted with the label side up, into the Nexus IB10 Analyzer tray immediately after injecting the sample into the disc.
- Do not turn the disc upside down.
- This is a quantitative test; therefore, no visual interpretation of the results should be made.

STORAGE AND STABILITY

- Store the IB10 sphingotest® PCT between 2 and 8 °C (35 to 46 °F) until the expiration date printed on the pouch is reached.
- The IB10 sphingotest® PCT in its sealed pouch is stable at 18 to 30 °C/64 to 86 °F for 30 days, provided the expiration date printed on the pouch is not exceeded.

SAMPLE COLLECTION AND PREPARATION

- The IB10 sphingotest® PCT is to be run using EDTA whole blood or plasma samples.
- It is recommended that samples be tested as soon as possible after collection.
- Whole blood must be tested within 24 hours of collection.
- Plasma samples should be kept frozen at -20 °C (-4 °F) or lower if longer storage is required.
- Allow samples to equilibrate to room temperature (19 to 25 °C/66 to 77 °F) prior to testing.

PROCEDURE

Nexus IB10 Analyzer

Consult the Nexus IB10 Analyzer User Manual

For Analyzer installation, start up and complete instructions for use refer to the **Nexus IB10 User Manual**. Operator must consult the User Manual prior to use to become familiar with the proper operation and quality control procedures.

PERFORMING SYSTEM CHECK AND DISC CALIBRATION

Each time the Nexus IB10 Analyzer is turned on, a Self-Check is automatically performed. The QR code on each Test disc contains information for disc calibration, which the Analyzer automatically reads when running a test.

RUNNING QC WITH EXTERNAL CONTROLS

The manufacturer recommends the use of commercially available PCT Controls (please refer to section **Materials/Equipment Required But Not Provided**).

Ensure that the external controls are handled and prepared according to the corresponding Instructions for Use (IFU).

1. Remove an unopened Test pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes prior to testing.
2. Open the pouch and remove the Test disc.
3. Place the Test disc on a level surface.
4. On the Nexus IB10 Analyzer press **New Analysis**.
5. The Analyzer will perform a general system check.
6. Enter the external control ID manually (up to 20 characters can be used for the ID) or enter the external control ID by using the barcode scanner.
7. Mix the external quality control vial by gently inverting the vial several times before sampling.
8. **Testing External Control Samples on the IB10 sphingotest® PCT**
 - Using a precision pipette (fixed or adjusted to 500 µL) slowly draw well-mixed external quality control sample into the pipette tip.
 - Positioning the tapered pipette tip at a 45° angle, pierce the **X** on the red dot to expose the sample channel inlet.
 - Slowly express the external control sample into the inlet applying minimum, but continuous force, on the pipette plunger.
 - Express sample to the first stop on the pipette, at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
 - Press **QC** on the Nexus IB10 Analyzer display.
 - When the tray opens, insert the filled Test disc into the tray and press **Run**.
 - The tray will close and perform a disc validity check.
 - A screen appears for selecting quality control materials (please refer to the **Quality Control Settings** section of the **Nexus IB10 Analyzer User Manual** for how to update quality control material [external controls]).
 - Select the quality control material to be tested.
 - Press **OK** on the Nexus IB10 Analyzer display.
 - In 20 minutes, the Nexus IB10 Analyzer will display the results on the screen.
 - Results will print automatically (if selected during Set Up) or press **Print**.
 - When the test is complete, analyze and compare the result with the expected value reported in the external control IFU for the external control level as measured using the IB10 sphingotest® PCT.
 - Remove the Test disc and discard in appropriate receptacle.
 - If the external control result is outside the expected range, refer to the Quality Control section below.

Note: If the test run is cancelled before a Test result is displayed, the Test disc cannot be reused and should be disposed of appropriately.

TESTING PATIENT SAMPLES ON THE NEXUS IB10 ANALYZER

1. Remove an unopened Test pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
2. Open the pouch and remove the Test disc.
3. Place the Test disc on a level surface.
4. On the Nexus IB10 Analyzer press **New Analysis**.
5. The Analyzer will perform a general system check.
6. Enter the Patient ID manually (up to 20 characters can be used for the ID) or enter the Patient ID by using the barcode scanner.
7. Mix the whole blood patient sample tube by gently inverting the tube several times before testing.

8. Testing Patient Samples on the IB10 sphingotest® PCT

- Using a precision pipette (fixed or adjusted to 500 µL) slowly draw well-mixed patient sample into the pipette tip.
- Positioning the tapered pipette tip at a 45° angle, pierce the **X** on the red dot to expose the sample channel inlet.
- Slowly express the patient sample into the inlet applying minimum, but continuous force, on the pipette plunger.
- Express sample to the first stop on the pipette, at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
- Press **OK** on the Nexus IB10 Analyzer display.
- When the tray opens, insert the filled Test disc into the tray and press **Run**.
- In 20 minutes, the Nexus IB10 Analyzer will display the results on the screen.
- Results will print automatically (if selected during Set Up) or press **Print**.
- Remove the Test disc and discard in appropriate receptacle.

Note: If the test run is cancelled before a test result is displayed, the Test disc cannot be reused and should be disposed of appropriately.

TRACEABILITY OF CALIBRATION

Calibration of the IB10 sphingotest® PCT is traceable to in-house reference calibrators.

INTERPRETATION OF RESULTS

The range of PCT concentrations reported by the Nexus IB10 Analyzer is 0.3 µg/L to 10.00 µg/L. Results below or above this range will be shown as “< 0.3 µg/L” or “> 10.00 µg/L,” respectively.

Quality Control

EXTERNAL CONTROLS

Good laboratory practice includes the use of external controls to ensure proper test performance. It is recommended that prior to using a new lot of IB10 sphingotest® PCT, the performance of the lot should be confirmed by testing with external controls (see section **Materials/Equipment Required But Not Provided**) to ensure the Test will deliver the correct result. The frequency of quality control testing should be determined according to individual laboratory standard quality control procedures. Upon confirmation of the expected results, the Test discs are ready for use with patient samples. Controls should also be used any time the validity of test results is questionable. If external controls do not perform as expected, do not use the IB10 sphingotest® PCT and contact the Distributor in your area for Technical Assistance.

INTERNAL CONTROL

Each IB10 sphingotest® PCT has a built in positive procedural control. The Nexus IB10 Analyzer automatically determines the presence of this control thereby confirming that the test run has delivered a valid result. If the control does not form or if it is not recognized by the Analyzer, the test result is considered “invalid” and the test must be repeated.

LIMITATIONS

- For diagnostic purposes, the test result should always be evaluated and interpreted in conjunction with the patient’s medical history, symptoms and other clinical information.
- As with any assay employing mouse antibodies, the possibility exists for interference by heterophilic antibodies in the sample.^{6,7} Specimens from individuals who have been regularly exposed to animals or who have been treated with animal products may contain these antibodies.

Other substances and/or factors not listed, e.g. technical or procedural error, may interfere with the test and cause inaccurate results. Highly icteric samples may lead to invalid control line or results.

Nexus Dx, Inc. offers products for their intended use. Refer to the specific product literature for the intended use statements for each product. Product claims are subject to change. Nexus Dx, Inc.'s expressed and implied warranties (inclusive of implied warranties of merchantability and fitness) are conditional upon adherence to, or observance of Nexus Dx, Inc. published directions with regard to the use of Nexus Dx, Inc. products. Under no circumstances will Nexus Dx, Inc. be liable for any indirect or consequential damages.

For Technical Assistance please contact your local Distributor.

EXPECTED VALUES

It is recommended that each laboratory establish its own expected values that represent the population that is to be evaluated at their facility.

All measurements with EDTA whole blood samples collected from 100 apparently healthy individuals were found as < 0.3 µg/L for the IB10 sphingotest® PCT.

DIAGNOSIS OF SEPSIS⁵

SIRS, Sepsis, Severe Sepsis, and Septic Shock were categorized according to the criteria of the consensus conference of the American College of Chest Physicians/Society of Critical Care Medicine.

Table 1

µg/L PCT	Analysis
< 0.5	Local bacterial infection is possible. Systemic infection (sepsis) is not likely. Low risk for progression to severe systemic infection (severe sepsis).
	 <ul style="list-style-type: none"> • PCT levels below 0.5 µg/L do not exclude an infection, because localized infections (without systemic signs) may be associated with such low levels. • Also if the PCT measurement is done very early after following bacterial challenge (usually < 6 hours), these values may still be low. In this case, PCT should be re-assessed 6–24 hours later.
> 0.5 and < 2	Systemic infection (sepsis) is possible, but various conditions are known to induce PCT as well (see below).
	Moderate risk for progression to severe systemic infection (severe sepsis). The patient should be closely monitored both clinically and by re-assessing PCT within 6–24 hours.
> 2 and < 10	Systemic infection (sepsis) is likely, unless other causes are known.
	High risk for progression to severe systemic infection (severe sepsis).
> 10	Important systemic inflammatory response, almost exclusively due to severe bacterial sepsis or septic shock.
	High likelihood of severe sepsis or septic shock.

Note:

Increased PCT levels may not always be related to infection.

There are a few situations described where PCT can be elevated by non-infectious causes.⁸⁻¹⁶

The assay is **not intended** to be used to diagnose or determine any of these causes. It is **only intended** for the diagnosis of sepsis, the assessment of the degree of severity, and monitoring of the change of sepsis severity.

Low PCT levels do not automatically exclude the presence of bacterial infection.

Such low levels may be obtained, e.g., during the early course of infections, in localized infections and in subacute endocarditis. Therefore, follow-up and re-evaluation of PCT in clinical suspicion of infection is pivotal. The PCT measuring technique should be chosen dependent on intended clinical use.

REFERENCES

1. Russell JA. Management of sepsis. N. Engl J Med 2006 Oct 19;355(16):1699-713.
2. Christ-Crain M, Müller B. Procalcitonin in bacterial infections – hype, hope or more or less? Swiss Med Wkly 2005;135(31-32):451-60.
3. Oczenski W, Fitzgerald RD, Schwarz S. Procalcitonin: a new parameter for the diagnosis of bacterial infection in the peri-operative period.

Eur J Anaesthesiol 1998;15(2):202-209.

4. Harbarth S, Holeckova K, Froidevaux C, et al. Diagnostic value of procalcitonin, interleukin-6, and interleukin-8 in critically ill patients admitted with suspected sepsis. *Am J Respir Crit Care Med* 2001;164(3):396-402.
5. Meisner M. Procalcitonin (PCT) – Biochemistry and clinical diagnosis – First Edition. Uni-Med 2010;ISBN: 978-3-8374-1241-3.
6. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45:942-6.
7. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15:108-15.
8. Assumma M, Signore F, Pacifico L et al. Serum procalcitonin concentrations in term delivering mothers and their healthy offspring: a longitudinal study. *Clin Chem*. 2000 Oct;46(10):1583-7.
9. Sponholz C, Sakr Y, Reinhart K et al. Diagnostic value and prognostic implications of serum procalcitonin after cardiac surgery: a systematic review of the literature. *Crit Care*. 2006;10(5):R145.
10. Dehne MG, Sablotzki A, Hoffmann A et al. Alterations of acute phase reaction and cytokine production in patients following severe burn injury. *Burns*. 2002 Sep;28(6):535-42.
11. Sabat R, Höflich C, Döcke WD et al. Massive elevation of procalcitonin plasma levels in the absence of infection in kidney transplant patients treated with pan-T-cell antibodies. *Intensive Care Med*. 2001 Jun;27(6):987-91.
12. Dornbusch HJ, Strenger V, Kerbl R et al. Procalcitonin--a marker of invasive fungal infection? *Support Care Cancer*. 2005 May;13(5):343-6.
13. Uzzan B, Izri A, Durand R et al. Serum procalcitonin in uncomplicated falciparum malaria: a preliminary study. *Travel Med Infect Dis*. 2006 Mar;4(2):77-80.
14. Picariello C, Lazzeri C, Chiostrri M et al. Procalcitonin in patients with acute coronary syndromes and cardiogenic shock submitted to percutaneous coronary intervention. *Intern Emerg Med*. 2009 Oct;4(5):403-8.
15. Avrillon V, Locatelli-Sanchez M, Folliet L et al. Lung cancer may increase serum procalcitonin level. *Infect Disord Drug Targets*. 2015;15(1):57-63.
16. Bihan H, Becker KL, Snider RH et al. Calcitonin precursor levels in human medullary thyroid carcinoma. *Thyroid*. 2003 Aug;13(8):819-22.

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REFERENZEN

1. Russell JA. Management of sepsis. *N. Engl J Med* 2006 Oct 19;355(16):1699-713.
2. Christ-Crain M, Müller B. Procalcitonin in bacterial infections – hype, hope or more or less? *Swiss Med Wkly* 2005;135(31-32):451-60.
3. Oczenski W, Fitzgerald RD, Schwarz S. Procalcitonin: a new parameter for the diagnosis of bacterial infection in the peri-operative period. *Eur J Anaesthesiol* 1998;15(2):202-209.
4. Harbarth S, Holeckova K, Froidevaux C, et al. Diagnostic value of procalcitonin, interleukin-6, and interleukin-8 in critically ill patients admitted with suspected sepsis. *Am J Respir Crit Care Med* 2001;164(3):396-402.
5. Meisner M. Procalcitonin (PCT) – Biochemistry and clinical diagnosis – First Edition. Uni-Med 2010;ISBN: 978-3-8374-1241-3.
6. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45:942-6.
7. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15:108-15.
8. Assumma M, Signore F, Pacifico L et al. Serum procalcitonin concentrations in term delivering mothers and their healthy offspring: a longitudinal study. *Clin Chem*. 2000 Oct;46(10):1583-7.
9. Sponholz C, Sakr Y, Reinhart K et al. Diagnostic value and prognostic implications of serum procalcitonin after cardiac surgery: a systematic review of the literature. *Crit Care*. 2006;10(5):R145.
10. Dehne MG, Sablotzki A, Hoffmann A et al. Alterations of acute phase reaction and cytokine production in patients following severe burn injury. *Burns*. 2002 Sep;28(6):535-42.
11. Sabat R, Höflich C, Döcke WD et al. Massive elevation of procalcitonin plasma levels in the absence of infection in kidney transplant patients treated with pan-T-cell antibodies. *Intensive Care Med*. 2001 Jun;27(6):987-91.
12. Dornbusch HJ, Strenger V, Kerbl R et al. Procalcitonin--a marker of invasive fungal infection? *Support Care Cancer*. 2005 May;13(5):343-6.
13. Uzzan B, Izri A, Durand R et al. Serum procalcitonin in uncomplicated falciparum malaria: a preliminary study. *Travel Med Infect Dis*. 2006 Mar;4(2):77-80.
14. Picariello C, Lazzeri C, Chiostrì M et al. Procalcitonin in patients with acute coronary syndromes and cardiogenic shock submitted to percutaneous coronary intervention. *Intern Emerg Med*. 2009 Oct;4(5):403-8.
15. Avrillon V, Locatelli-Sanchez M, Folliet L et al. Lung cancer may increase serum procalcitonin level. *Infect Disord Drug Targets*. 2015;15(1):57-63.
16. Bihan H, Becker KL, Snider RH et al. Calcitonin precursor levels in human medullary thyroid carcinoma. *Thyroid*. 2003 Aug;13(8):819-22.

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REFERÊNCIAS

1. Russell JA. Management of sepsis. *N. Engl J Med* 2006 Oct 19;355(16):1699-713.
2. Christ-Crain M, Müller B. Procalcitonin in bacterial infections – hype, hope or more or less? *Swiss Med Wkly* 2005;135(31-32):451-60.
3. Oczenski W, Fitzgerald RD, Schwarz S. Procalcitonin: a new parameter for the diagnosis of bacterial infection in the peri-operative period. *Eur J Anaesthesiol* 1998;15(2):202-209.
4. Harbarth S, Holeckova K, Froidevaux C, et al. Diagnostic value of procalcitonin, interleukin-6, and interleukin-8 in critically ill patients admitted with suspected sepsis. *Am J Respir Crit Care Med* 2001;164(3):396-402.
5. Meisner M. Procalcitonin (PCT) – Biochemistry and clinical diagnosis – First Edition. Uni-Med 2010;ISBN: 978-3-8374-1241-3.
6. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45:942-6.
7. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15:108-15.
8. Assumma M, Signore F, Pacifico L et al. Serum procalcitonin concentrations in term delivering mothers and their healthy offspring: a longitudinal study. *Clin Chem.* 2000 Oct;46(10):1583-7.
9. Sponholz C, Sakr Y, Reinhart K et al. Diagnostic value and prognostic implications of serum procalcitonin after cardiac surgery: a systematic review of the literature. *Crit Care.* 2006;10(5):R145.
10. Dehne MG, Sablotzki A, Hoffmann A et al. Alterations of acute phase reaction and cytokine production in patients following severe burn injury. *Burns.* 2002 Sep;28(6):535-42.
11. Sabat R, Höflich C, Döcke WD et al. Massive elevation of procalcitonin plasma levels in the absence of infection in kidney transplant patients treated with pan-T-cell antibodies. *Intensive Care Med.* 2001 Jun;27(6):987-91.
12. Dornbusch HJ, Strenger V, Kerbl R et al. Procalcitonin--a marker of invasive fungal infection? *Support Care Cancer.* 2005 May;13(5):343-6.
13. Uzzan B, Izri A, Durand R et al. Serum procalcitonin in uncomplicated falciparum malaria: a preliminary study. *Travel Med Infect Dis.* 2006 Mar;4(2):77-80.
14. Picariello C, Lazzeri C, Chiostrì M et al. Procalcitonin in patients with acute coronary syndromes and cardiogenic shock submitted to percutaneous coronary intervention. *Intern Emerg Med.* 2009 Oct;4(5):403-8.
15. Avrillon V, Locatelli-Sanchez M, Folliet L et al. Lung cancer may increase serum procalcitonin level. *Infect Disord Drug Targets.* 2015;15(1):57-63.
16. Bihan H, Becker KL, Snider RH et al. Calcitonin precursor levels in human medullary thyroid carcinoma. *Thyroid.* 2003 Aug;13(8):819-22.

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- Russell JA. Management of sepsis. N. Engl J Med 2006 Oct 19;355(16):1699-713.
- Christ-Crain M, Müller B. Procalcitonin in bacterial infections – hype, hope or more or less? Swiss Med Wkly 2005;135(31-32):451-60.
- Oczenski W, Fitzgerald RD, Schwarz S. Procalcitonin: a new parameter for the diagnosis of bacterial infection in the peri-operative period. Eur J Anaesthesiol 1998;15(2):202-209.
- Harbarth S, Holeckova K, Froidevaux C, et al. Diagnostic value of procalcitonin, interleukin-6, and interleukin-8 in critically ill patients admitted with suspected sepsis. Am J Respir Crit Care Med 2001;164(3):396-402.
- Meisner M. Procalcitonin (PCT) – Biochemistry and clinical diagnosis – First Edition. Uni-Med 2010; ISBN: 978-3-8374-1241-3.
- Kricka LJ. Human anti-animal antibody interference in immunological assays. Clin Chem 1999;45:942-6.
- Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. J Clin Immunoassay 1992;15:108-15.
- Assumma M, Signore F, Pacifico L et al. Serum procalcitonin concentrations in term delivering mothers and their healthy offspring: a longitudinal study. Clin Chem. 2000 Oct;46(10):1583-7.
- Sponholz C, Sakr Y, Reinhart K et al. Diagnostic value and prognostic implications of serum procalcitonin after cardiac surgery: a systematic review of the literature. Crit Care. 2006;10(5):R145.
- Dehne MG, Sablotzki A, Hoffmann A et al. Alterations of acute phase reaction and cytokine production in patients following severe burn injury. Burns. 2002 Sep;28(6):535-42.
- Sabat R, Höflich C, Döcke WD et al. Massive elevation of procalcitonin plasma levels in the absence of infection in kidney transplant patients

treated with pan-T-cell antibodies. *Intensive Care Med.* 2001 Jun;27(6):987-91.

12. Dornbusch HJ, Strenger V, Kerbl R et al. Procalcitonin--a marker of invasive fungal infection? *Support Care Cancer.* 2005 May;13(5):343-6.
13. Uzzan B, Izri A, Durand R et al. Serum procalcitonin in uncomplicated falciparum malaria: a preliminary study. *Travel Med Infect Dis.* 2006 Mar;4(2):77-80.
14. Picariello C, Lazzeri C, Chiostrì M et al. Procalcitonin in patients with acute coronary syndromes and cardiogenic shock submitted to percutaneous coronary intervention. *Intern Emerg Med.* 2009 Oct;4(5):403-8.
15. Avrillon V, Locatelli-Sanchez M, Folliet L et al. Lung cancer may increase serum procalcitonin level. *Infect Disord Drug Targets.* 2015;15(1):57-63.
16. Bihan H, Becker KL, Snider RH et al. Calcitonin precursor levels in human medullary thyroid carcinoma. *Thyroid.* 2003 Aug;13(8):819-22.

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REFERENSER

- Russell JA. Management of sepsis. N. Engl J Med 2006 Oct 19;355(16):1699-713.
- Christ-Crain M, Müller B. Procalcitonin in bacterial infections – hype, hope or more or less? Swiss Med Wkly 2005;135(31-32):451-60.
- Oczenski W, Fitzgerald RD, Schwarz S. Procalcitonin: a new parameter for the diagnosis of bacterial infection in the peri-operative period. Eur J Anaesthesiol 1998;15(2):202-209.
- Harbarth S, Holeckova K, Froidevaux C, et al. Diagnostic value of procalcitonin, interleukin-6, and interleukin-8 in critically ill patients admitted with suspected sepsis. Am J Respir Crit Care Med 2001;164(3):396-402.
- Meisner M. Procalcitonin (PCT) – Biochemistry and clinical diagnosis – First Edition. Uni-Med 2010; ISBN: 978-3-8374-1241-3.
- Kricka LJ. Human anti-animal antibody interference in immunological assays. Clin Chem 1999;45:942-6.
- Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. J Clin Immunoassay 1992;15:108-15.
- Assumma M, Signore F, Pacifico L et al. Serum procalcitonin concentrations in term delivering mothers and their healthy offspring: a longitudinal study. Clin Chem. 2000 Oct;46(10):1583-7.

9. Sponholz C, Sakr Y, Reinhart K et al. Diagnostic value and prognostic implications of serum procalcitonin after cardiac surgery: a systematic review of the literature. *Crit Care*. 2006;10(5):R145.
10. Dehne MG, Sablotzki A, Hoffmann A et al. Alterations of acute phase reaction and cytokine production in patients following severe burn injury. *Burns*. 2002 Sep;28(6):535-42.
11. Sabat R, Höflich C, Döcke WD et al. Massive elevation of procalcitonin plasma levels in the absence of infection in kidney transplant patients treated with pan-T-cell antibodies. *Intensive Care Med*. 2001 Jun;27(6):987-91.
12. Dornbusch HJ, Strenger V, Kerbl R et al. Procalcitonin--a marker of invasive fungal infection? *Support Care Cancer*. 2005 May;13(5):343-6.
13. Uzzan B, Izri A, Durand R et al. Serum procalcitonin in uncomplicated falciparum malaria: a preliminary study. *Travel Med Infect Dis*. 2006 Mar;4(2):77-80.
14. Picariello C, Lazzeri C, Chiostrì M et al. Procalcitonin in patients with acute coronary syndromes and cardiogenic shock submitted to percutaneous coronary intervention. *Intern Emerg Med*. 2009 Oct;4(5):403-8.
15. Avrillon V, Locatelli-Sanchez M, Folliet L et al. Lung cancer may increase serum procalcitonin level. *Infect Disord Drug Targets*. 2015;15(1):57-63.
16. Bihan H, Becker KL, Snider RH et al. Calcitonin precursor levels in human medullary thyroid carcinoma. *Thyroid*. 2003 Aug;13(8):819-22.