## EN

### Intended Use

Nexus IB10 EQC is for testing performance check of the Nexus IB10 system.

### Summary

Nexus IB10 EQC consists of one disc with 3 different type of strips. Performance checks of the system should ideally be carried out on a daily basis. The Nexus IB10 analyzer should be checked using the Nexus IB10 EQC according to the applicable local or national regulations.

### Test Procedure

[New Analysis] → [QC] → Insert EQC Disc → [Run] → The results will be displayed within two minutes. "Pass" means that the system of the Nexus IB10 analyzer is functioning properly. "Fail" means that the value is outside the confidence interval. If "Fail" appears continuously, please contact the dealer or the service center.

### Precaution

- Bring the Nexus IB10 EQC to room temperature before starting the measurement.
- For information on the optimum ambient temperature for the Nexus IB10 analyzer, please read the notes in the operator's manual.
- Do not insert any sample material to the Nexus IB10 EQC.
- Do not expose the Nexus IB10 EQC to sunlight.
- Protect the Nexus IB10 EQC from dust and moisture.
- Close the EQC container immediately after performance checks to protect from the UV light, dust and moisture for later tests.

### Storage and Stability

Store at 2°C – 30°C

- The Nexus IB10 EQC is stable up to the expiry date printed on the label.
- Do not use the Nexus IB10 EQC after the expiry date. Nexus Dx does not guarantee the test results when the date is passed.

### Order Information

Nexus IB10 EQC, **REF** IVR-IB60

### Materials are not included but required

Nexus IB10 Analyzer (BCA-IB10)

Please refer to the end of the document for an explanation of the symbols used in the labels.

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Nexus Dx, Inc.
Theragenesis GmbH
SphingoTec GmbH