
English

IB10 DPP3


For the Quantitative Determination of
Dipeptidyl Peptidase 3 in EDTA Whole Blood and Plasma

Explanation of Symbols

	CE Mark of Conformity
	Manufacturer
	Catalog Number
	Expiry date/Use by
	Lot number
	<i>In Vitro</i> Diagnostic medical device
	Consult instructions for use
	Store between 2 °C and 8 °C
	European Authorized Representative
	Contains sufficient for <n> test
	Do not reuse
	Serial Number
	Insert disc label side up

Nexus Dx, Inc.,
6759 Mesa Ridge Road
San Diego, California 92121 USA
Telephone: 1 (858) 410 4600
Fax: 1 (858) 410 4700

TheraGenesis GmbH Bahnhofstrasse 5
55276 Oppenheim, Germany
Tel: +49 (0) 151 506 403 14

022-00072 REV F

IB10 DPP3

For *In Vitro* Diagnostic Use

INTENDED USE

The IB10 DPP3 is a rapid point-of-care (POC) immunoassay for the *in vitro* quantitative determination of human dipeptidyl peptidase 3 (DPP3) in EDTA whole blood and plasma. DPP3 measurement in conjunction with clinical assessments and other laboratory findings is used as an aid in the prediction of short-term survival in patients with sepsis or septic shock. This Test is intended for use in conjunction with the Nexus IB10 Analyzer and provides quantitative results in 22 minutes.

The IB10 DPP3 is designed for professional use only and may be used in hospital central laboratories and in alternate care settings such as emergency departments, critical care units, and other POC sites where near patient testing is practiced.

SUMMARY AND EXPLANATION OF THE TEST

The IB10 DPP3 quantitatively measures the level of DPP3 circulating in the bloodstream.¹ DPP3 is a cytosolic enzyme which is released into the bloodstream upon cell death accruing in various acute care settings such as cardiogenic shock, acute myocardial infarction, sepsis, trauma and burns. Degradation of its substrates, including angiotensins and enkephalins, interrupts peptide hormone signaling pathways causing organ dysfunction.²⁻⁴ In cardiogenic shock patients, elevated circulating DPP3 is a myocardial depressant factor associated to organ dysfunction.⁵ Circulating DPP3 may early predict alteration of hemodynamics in cardiogenic shock patients and rapid normalization of circulating DPP3 levels on the first day indicate low rate of refractory shock and high survival.⁶

DIPEPIDYL PEPTIDASE III

Human dipeptidyl peptidase 3 (DPP3), also known as red cell angiotensinase, is a member of the M49 family of zinc-dependent metallopeptidases.⁷ DPP3 has a molecular mass of about 83 kDa and is ubiquitously expressed and highly conserved in prokaryotes and eukaryotes.⁸ DPP3 is detected in organ homogenates and several bodily fluids, such as blood plasma, retroplacental serum, seminal plasma, and cerebrospinal fluid.^{3,9-11} The enzyme cleaves various bioactive peptides, such as angiotensins, enkephalins, and endomorphins.⁸ The most prominent substrate of DPP3 is angiotensin II (Ang II), the main effector of the renin–angiotensin system (RAS), implicating extracellular DPP3 in the regulation of the RAS system.⁴ The RAS is activated in cardiovascular diseases, sepsis and septic shock and Ang II, in particular, has been shown to modulate many cardiovascular functions including the control of blood pressure and cardiac remodeling.¹²⁻¹⁵

PRINCIPLE

The Nexus IB10 immunochemistry system combines chemistry with microfluidics and centrifugal flow to rapidly prepare a cell free plasma from whole blood that can then be moved through a channel to rehydrate, solubilize and mix with freeze dried immunoconjugates. Using a combination of active flow and capillary action, the test sample is quantitatively measured in 22 minutes with an optical signal level proportional to the analyte(s) concentration.

After addition of the patient sample, the entire test is performed within the Nexus IB10 Analyzer, which provides control of the temperature of the disc, as well as the reaction sequence, centrifugal flow, mixing, incubation time, final signal measurement, quantitation and reporting of results. The Test disc

includes a positive internal control to ensure that the Test has operated properly. Each lot is calibrated to provide accurate analyte concentration. Lot specific calibration along with additional information such as the lot expiration date is contained on a QR code label affixed to each disc. It is recommended that external controls be tested at appropriate time intervals to confirm that the system and test lot are performing within acceptable limits.

REAGENTS

The IB10 DPP3 contains all required reagents to evaluate the level of DPP3 including dye-conjugated monoclonal anti-DPP3 antibody, biotin-conjugated monoclonal anti-DPP3 antibody and streptavidin immobilized at the detection area on the disc.

MATERIALS PROVIDED

Each box contains the following:

- IB10 DPP3 discs, each individually sealed in a foil pouch with a desiccant.
- Instructions for Use (IFU).

MATERIALS/EQUIPMENT REQUIRED BUT NOT PROVIDED

1. Nexus IB10 Analyzer - Model #BCA-IB10.
2. Commercially available DPP3 Controls for external Quality Control (QC). **Contact the Distributor in your area** for recommended external QC materials or related technical assistance.
3. Calibrated reusable fixed or variable volume pipette gun with high precision and accuracy capable of delivering 500 µL of whole blood or plasma.
4. Disposable pipette tips capable of accepting and delivering 500 µL of whole blood or plasma.

PRECAUTIONS AND WARNINGS

- For *in vitro* diagnostic use only.
- Carefully follow the Instructions for Use.
- **Check your Nexus IB10 Analyzer for the latest analyzer software before testing IB10 discs (Refer to the Nexus IB10 manual for specific instructions).**
- **Not updating the software can result in error message. Login or Register on the Nexus-DX website to access Customer Downloads (all 5 files are required for the software update).**
- Wear disposable gloves while handling samples.
- Handle samples with care. Samples and used test discs should be treated as potentially infectious and should be discarded as biohazardous material according to local regulations.
- Thoroughly wash hands following handling.
- The result obtained from the IB10 DPP3 does not provide a definitive diagnosis and should be interpreted by a physician in conjunction with other laboratory test results and according to current medical guidelines and patient clinical findings.
- Keep the test disc in the sealed pouch until ready for use.
- Do not use the test disc if the pouch is damaged or the seal is broken.
- Do not use the test disc after the expiration date printed on the pouch.
- Prior to use, place the unopened pouch at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
- Always pay attention to cleanliness when handling the test disc. Avoid any contamination from fingerprints or foreign substances. Do not contaminate the sample channel inlet.

- Do not drop or damage the test disc.
- The test disc should be inserted with the label side up, into the Nexus IB10 Analyzer tray immediately after injecting the sample into the disc.
- Do not turn the disc upside down.
- This is a quantitative test; therefore no visual interpretation of the results should be made.

STORAGE AND STABILITY

- Store the IB10 DPP3 disc between 2 and 8 °C (35 and 46 °F) until the expiration date printed on the pouch is reached.
- The IB10 DPP3 in its sealed pouch is stable between 18 and 30 °C (64 and 86 °F) for 30 days, provided the expiration date printed on the pouch is not exceeded.

SAMPLE COLLECTION AND PREPARATION

- The IB10 DPP3 is to be run using EDTA whole blood or plasma samples.
- **It is recommended that EDTA whole blood is tested as soon as possible after blood draw or store it up to a maximum of 2 hours at room temperature (19 to 25 °C/66 to 77 °F) or 2-8°C (35-45 °F) until testing. If testing cannot be completed within 2 hours, the plasma should be separated and stored at -20 °C (- 4 °F) or lower until it can be tested.**
- Allow samples to equilibrate to room temperature (19 to 25 °C/66 to 77 °F) and mix gently prior to testing.

PROCEDURE

Nexus IB10 Analyzer

Consult the Nexus IB10 Analyzer User Manual

For Analyzer installation, start up and complete instructions for use refer to the **Nexus IB10 User Manual**. Operator must consult the User Manual prior to use to become familiar with the proper operation and quality control procedures.

PERFORMING SYSTEM CHECK AND DISC CALIBRATION

Each time the Nexus IB10 Analyzer is turned on, a Self Check is automatically performed. The QR code on each test disc contains information for disc calibration which the Analyzer automatically reads when running a test.

RUNNING QC WITH EXTERNAL CONTROLS

The manufacturer recommends the use of commercially available DPP3 Controls (please refer to section **Materials/Equipment Required But Not Provided**). Ensure that the DPP3 controls are handled and prepared according to the corresponding Instructions for Use (IFU).

1. Remove an unopened test pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes prior to testing.
2. Open the pouch and remove the test disc.
3. Place the test disc on a level surface.
4. On the Nexus IB10 Analyzer press **New Analysis**.
5. The Analyzer will perform a general system check.
6. Enter the external control ID manually (up to 20 characters can be used for the ID) or enter the external control ID by using the barcode scanner.

7. Mix the external quality control vial by gently inverting the vial several times before sampling.

8. Testing External Control Samples on the IB10 DPP3.

- Using a precision pipette (fixed or adjusted to 500 µL), slowly draw well-mixed external quality control sample into the pipette tip.
- Positioning the tapered pipette tip at a 45° angle, pierce the **X** on the red dot to expose the sample channel inlet.
- Slowly express the external control sample into the inlet applying minimum, but continuous force, on the pipette plunger.
- Express sample to the first stop on the pipette, at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
- Press **QC** on the Nexus IB10 Analyzer display.
- When the tray opens, insert the filled Test disc into the tray and press **Run**.
- The tray will close and perform a disc validity check.
- A screen appears for selecting quality control materials (please refer to the **Quality Control Settings** section of the **Nexus IB10 Analyzer User Manual** for how to update quality control material [external controls]).
- Select the quality control material to be tested.
- Press **OK** on the Nexus IB10 Analyzer display.
- In 22 minutes, the Nexus IB10 Analyzer will display the results on the screen.
- Results will print automatically (if selected during Set Up) or press **Print**.
- When the test is complete, analyze and compare the result with the Expected Value reported in the external control IFU for the external control level as measured using the IB10 DPP3.
- Remove the test disc and discard in appropriate receptacle.
- If the external control result is outside the expected range, refer to the Quality Control section below.

Note: If the test run is cancelled before a test result is displayed, the test disc cannot be reused and should be disposed of appropriately.

TESTING PATIENT SAMPLES ON THE NEXUS IB10 ANALYZER

1. Remove an unopened Test pouch from refrigeration and place it at room temperature (19 to 25 °C/66 to 77 °F) for at least 15 minutes.
2. Open the pouch and remove the Test disc.
3. Place the Test disc on a level surface.
4. On the Nexus IB10 Analyzer press **New Analysis**.
5. The Analyzer will perform a general system check.
6. Enter the Patient ID manually (up to 20 characters can be used for the ID) or enter the Patient ID by using the barcode scanner.
7. Mix the whole blood patient sample tube by gently inverting the tube several times before testing.

8. Testing Patient Samples on the IB10 DPP3.

- Using a precision pipette (fixed or adjusted to 500 µL), slowly draw well-mixed patient sample into the pipette tip.
- Positioning the tapered pipette tip at a 45° angle, pierce the **X** on the red dot to expose the sample channel inlet.
- Slowly express the patient sample into the inlet applying minimum, but continuous force, on

the pipette plunger.

- Express sample to the first stop on the pipette, at a rate that allows the fluid to fill the channel completely and eliminates any back pressure that could result in sample splatter or the introduction of air bubbles.
- Press **OK** on the Nexus IB10 Analyzer display.
- When the tray opens, insert the filled Test disc into the tray and press **Run**.
- In 22 minutes, the Nexus IB10 Analyzer will display the results on the screen.
- Results will print automatically (if selected during Set Up) or press **Print**.
- Remove the test disc and discard in appropriate receptacle.

Note: If the test run is cancelled before a test result is displayed, the test disc cannot be reused and should be disposed of appropriately.

INTERPRETATION OF RESULTS

The range of DPP3 concentrations reported by the Nexus IB10 is 5 ng/mL to 150 ng/mL. Results below or above this range will be shown as "<5 ng/mL" or ">150 ng/mL", respectively.

Quality Control

EXTERNAL CONTROLS

Good laboratory practice includes the use of external controls to ensure proper test performance. It is recommended that prior to using a new lot or shipment of IB10 DPP3, the performance of the lot should be confirmed by testing with external controls (see section **Materials/ Equipment Required But Not Provided**) to ensure the test will deliver the correct result. The frequency of quality control testing should be determined according to individual laboratory standard quality control procedures. Upon confirmation of the expected results, the test discs are ready for use with patient samples. Controls should also be used any time the validity of test results is questionable. If external controls do not perform as expected, do not use the IB10 DPP3 and contact the Distributor in your area for Technical Assistance.

INTERNAL CONTROL

The IB10 DPP3 has a built in positive procedural control. The Nexus IB10 Analyzer automatically determines the presence of this control thereby confirming that the test run has delivered a valid result. If the control does not form or if it is not recognized by the Analyzer, the test result is considered "invalid" and the test must be repeated.

LIMITATIONS

The test result should be evaluated and interpreted in conjunction with the patient's medical history, symptoms and other clinical information.

Hemolysis of whole blood from specimen collection, mishandling of blood, reaction to transfusion, certain disease conditions, or other variables may cause the destruction of red blood cells, releasing DPP3 into the sample and impacting test results. This could possibly generate false positive results.

As with any assay employing mouse antibodies, the possibility exists for interference by heterophilic antibodies in the sample.^{16,17} Specimens from individuals who have been regularly exposed to animals or who have been treated with animal products may contain these antibodies.

Other substances and/or factors not listed, e.g. technical or procedural error, may interfere with the test

and cause inaccurate results.

Nexus Dx, Inc. offers products for their intended use. Refer to the specific product literature for the intended use statements for each product. Product claims are subject to change. Nexus Dx, Inc.'s expressed and implied warranties (inclusive of implied warranties of merchantability and fitness) are conditional upon adherence to, or observance of Nexus Dx, Inc. published directions with regard to the use of Nexus Dx, Inc. products. Under no circumstances will Nexus Dx, Inc. be liable for any indirect or consequential damages.

For Technical Assistance please contact your local Distributor.

Performance Characteristics

MEASURING RANGE

The IB10 DPP3 has been demonstrated to provide measurable results at DPP3 levels from 5 ng/mL to 150 ng/mL.

ANALYTICAL SENSITIVITY

The LoD (Limit of Detection) of the IB10 DPP3 is 5 ng/mL, determined according to Clinical and Laboratory Standards Institute (CLSI) Approved Guideline EP17-A2.¹⁸ The proportion of false positives (α) and false negatives (β) is less than 5% (limit of blank = 0.58 ng/mL). LoQ (Limit of Quantitation) is the lowest DPP3 concentration that can be reproducibly measured with a total coefficient of variation of at most 20%. It was determined to be 5 ng/mL.

INTERFERING SUBSTANCES

The following substances were tested for potential interference with DPP3 measurement in the IB10 DPP3 based on CLSI Approved Guideline EP7-A2 (Table 1).¹⁹ No significant interference was observed at the maximum concentration for substances listed in the table below.

Table 1.

Potential Interferents	Maximum Concentration
Bilirubin (Conjugated)	40 mg/dL
Bilirubin (Free)	60 mg/dL
Biotin	800 µg/L
Cholesterol	400 mg/dL
Creatinine	2 mg/dL
Glucose	1000 mg/dL
Hemoglobin	0.5 g/dL
Heparin	3000 U/L
Rheumatoid factor (RF)	1500 IU/mL
Total Protein	6 g/dL
Triglycerides	3000 mg/dL
Urea	250 mg/dL

HOOK EFFECT

No high dose hook effect is observed for DPP3 concentrations up to 6000 ng/mL.

PRECISION

The precision of DPP3 measurement of the IB10 DPP3 was determined using samples where DPP3 was added to normal human plasma at two concentrations (Table 2). The within run and total precision was determined over a 20 day period consisting of 2 runs per day with 2 replicates per run, totaling to 80 replications for each concentration level. The within-run and total precision were calculated according to the CLSI Approved Guideline EP05-A3.²⁰

Table 2.

Sample	Test lot	Mean (ng/mL)	Within-Run Precision		Total Precision	
			Std. dev. (ng/mL)	CV (%)	Std. dev. (ng/mL)	CV (%)
1	A	44.82	2.96	6.6%	3.14	7.0%
	B	47.28	2.73	5.8%	2.85	6.0%
2	A	108.15	10.19	9.4%	10.91	10.1%
	B	110.28	9.31	8.4%	9.59	8.7%

WHOLE BLOOD Vs. PLASMA CORRELATION

A comparison study was performed using matching EDTA whole blood and EDTA plasma samples. When performing a Bland-Altman analysis comparing the whole blood concentrations versus the corresponding plasma concentrations from the same subject samples (N=58), a mean difference of 9.6% (95% C.I. = [6.3%-13.0%]), and a correlation coefficient of $r = 0.98$ were obtained.

REFERENCES

1. Rehfeld L, Funk E, Jha S, Macheroux P, Melander O, Bergmann A. Novel Methods for the Quantification of Dipeptidyl Peptidase 3 (DPP3) Concentration and Activity in Human Blood Samples. J Appl Lab Med. 2018. doi:10.1373/jalm.2018.027995.
2. Wattiaux R, Wattiaux-de Coninck S, Thirion J, Gasingirwa MC, Jadot M. Lysosomes and Fas-mediated liver cell death. Biochem J. 2007 Apr 1;403(1):89-95. PubMed PMID: 17129211; PubMed Central PMCID: PMC1828900.
3. Pang X, Shimizu A, Kurita S, Zankov DP, Takeuchi K, Yasuda-Yamahara M, et al. Novel therapeutic role for dipeptidyl peptidase III in the treatment of hypertension. Hypertension 2016;68:630–41.
4. Kumar P, Reithofer V, Reisinger M, Wallner S, Pavkov-Keller T, Macheroux P, Gruber K. Substrate complexes of human dipeptidyl peptidase III reveal the mechanism of enzyme inhibition. Sci Rep. 2016 Mar 30;6:23787. doi: 10.1038/srep23787. PubMed PMID: 27025154; PubMed Central PMCID: PMC4824452.
5. Deniau B, Rehfeld L, Santos K, Dienelt A, Azibani F, Sadoune M, Kounde PR, Samuel LS, Tolpannen H, Lassus J, Harjola VP, Vodovar N, Bergmann A, Hartmann O, Mebazaa A, Blet A. Circulating dipeptidyl peptidase-3 is a myocardial depressant factor: DPP3 inhibition rapidly and sustainably improves hemodynamics, European Journal of Heart Failure 2019, accepted.
6. Takagi K, Blet A, Levy B, Deniau B, Azibani F, Feliot E, Bergmann A, Santos K, Hartmann O, Gayat E, Mebazaa A, Kimmoun A. Circulating dipeptidyl-peptidase 3 and alteration in hemodynamics in cardiogenic shock: Results from the OptimaCC Trial, European Journal of Heart Failure 2019, accepted.
7. Information on EC 3.4.14.4 - dipeptidyl-peptidase III, BRENDA. (n.d.). <https://www.brenda->

enzymes.org/enzyme.php?ecno=3.4.14.4

8. Prajapati SC, Chauhan SS. Dipeptidyl peptidase III: a multifaceted oligopeptide N-end cutter. *FEBS J* 2011;278: 3256–76.
9. Shimamori Y, Watanabe Y, Fujimoto Y. Purification and characterization of dipeptidyl aminopeptidase III from human placenta. *Chem Pharm Bull* 1986;34:3333–40.
10. Vanha-Perttula T. Dipeptidyl peptidase III and alanyl aminopeptidase in the human seminal plasma: origin and biochemical properties. *Clin Chim Acta* 1988;177:179–95.
11. Aoyagi T, Wada T, Kojima F, Nagai M, Harada S, Takeuchi T, et al. Enzymatic changes in cerebrospinal fluid of patients with Alzheimer-type dementia. *J Clin Biochem Nutr* 1993;14:133–9.
12. Dostal DE, Hunt RA, Kule CE, Bhat GJ, Karoor V, McWhinney CD, et al. Molecular mechanisms of angiotensin II in modulating cardiac function: intracardiac effects and signal transduction pathways. *J Mol Cell Cardiol* 1997;29:2893–902.
13. Roks A, Buikema H, Pinto YM, van Gilst WH. The reninangiotensin system and vascular function: the role of angiotensin II, angiotensin- converting enzyme, and alternative conversion of angiotensin I. *Heart Vessels* 1997;Suppl 12:119–24.
14. Ferrario CM. Role of angiotensin II in cardiovascular disease—therapeutic implications of more than a century of research. *J Renin- Angiotensin-Aldosterone Syst* 2006;7:3–14.
15. Corrêa TD, Takala J, Jakob SM. Angiotensin II in septic shock. *Crit Care* 2015;19:98.
16. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45:942–6.
17. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15:108–15.
18. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
19. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
20. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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3. Pang X, Shimizu A, Kurita S, Zankov DP, Takeuchi K, Yasuda-Yamahara M, et al. Novel therapeutic role for dipeptidyl peptidase III in the treatment of hypertension. *Hypertension* 2016;68:630–41.
4. Kumar P, Reithofer V, Reisinger M, Wallner S, Pavkov-Keller T, Macheroux P, Gruber K. Substrate complexes of human dipeptidyl peptidase III reveal the mechanism of enzyme inhibition. *Sci Rep*. 2016 Mar 30;6:23787. doi: 10.1038/srep23787. PubMed PMID: 27025154; PubMed Central PMCID: PMC4824452.
5. Deniau B, Rehfeld L, Santos K, Dienelt A, Azibani F, Sadoune M, Kounde PR, Samuel LS, Tolpannen H, Lassus J, Harjola VP, Vodovar N, Bergmann A, Hartmann O, Mebazaa A, Blet A. Circulating dipeptidyl peptidase-3 is a myocardial depressant factor: DPP3 inhibition rapidly and sustainably improves hemodynamics, *European Journal of Heart Failure* 2019, accepted.
6. Takagi K, Blet A, Levy B, Deniau B, Azibani F, Feliot E, Bergmann A, Santos K, Hartmann O, Gayat E, Mebazaa A, Kimmoun A. Circulating dipeptidyl-peptidase 3 and alteration in hemodynamics in cardiogenic shock: Results from the OptimaCC Trial, *European Journal of Heart Failure* 2019, accepted.
7. Information on EC 3.4.14.4 - dipeptidyl-peptidase III, BRENDA. (n.d.). <https://www.brenda-enzymes.org/enzyme.php?ecno=3.4.14.4>
8. Prajapati SC, Chauhan SS. Dipeptidyl peptidase III: a multifaceted oligopeptide N-end cutter. *FEBS J* 2011;278: 3256–76.
9. Shimamori Y, Watanabe Y, Fujimoto Y. Purification and characterization of dipeptidyl aminopeptidase III from human placenta. *Chem Pharm Bull* 1986;34:3333–40.
10. Vanha-Perttula T. Dipeptidyl peptidase III and alanyl aminopeptidase in the human seminal plasma: origin and biochemical properties. *Clin Chim Acta* 1988;177:179–95.
11. Aoyagi T, Wada T, Kojima F, Nagai M, Harada S, Takeuchi T, et al. Enzymatic changes in cerebrospinal fluid of patients with Alzheimer-type dementia. *J Clin Biochem Nutr* 1993;14:133–9.
12. Dostal DE, Hunt RA, Kule CE, Bhat GJ, Karoor V, McWhinney CD, et al. Molecular mechanisms of angiotensin II in modulating cardiac function: intracardiac effects and signal transduction pathways. *J Mol Cell Cardiol* 1997;29:2893–902.
13. Roks A, Buikema H, Pinto YM, van Gilst WH. The reninangiotensin system and vascular function: the role of angiotensin II, angiotensin- converting enzyme, and alternative conversion of angiotensin I. *Heart Vessels* 1997;Suppl 12:119–24.
14. Ferrario CM. Role of angiotensin II in cardiovascular disease—therapeutic implications of more than a century of research. *J Renin- Angiotensin-Aldosterone Syst* 2006;7:3–14.
15. Corrêa TD, Takala J, Jakob SM. Angiotensin II in septic shock. *Crit Care* 2015;19:98.
16. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45:942-6.
17. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15:108-15.
18. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
19. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
20. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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improves hemodynamics, *European Journal of Heart Failure* 2019, accepted.

6. Takagi K, Blet A, Levy B, Deniau B, Azibani F, Feliot E, Bergmann A, Santos K, Hartmann O, Gayat E, Mebazaa A, Kimmoun A. Circulating dipeptidyl-peptidase 3 and alteration in hemodynamics in cardiogenic shock: Results from the OptimaCC Trial, *European Journal of Heart Failure* 2019, accepted.
7. Information on EC 3.4.14.4 - dipeptidyl-peptidase III, BRENDA. (n.d.). <https://www.brenda-enzymes.org/enzyme.php?ecno=3.4.14.4>
8. Prajapati SC, Chauhan SS. Dipeptidyl peptidase III: a multifaceted oligopeptide N-end cutter. *FEBS J* 2011;278: 3256–76.
9. Shimamori Y, Watanabe Y, Fujimoto Y. Purification and characterization of dipeptidyl aminopeptidase III from human placenta. *Chem Pharm Bull* 1986;34:3333–40.
10. Vanha-Perttula T. Dipeptidyl peptidase III and alanyl aminopeptidase in the human seminal plasma: origin and biochemical properties. *Clin Chim Acta* 1988;177:179–95.
11. Aoyagi T, Wada T, Kojima F, Nagai M, Harada S, Takeuchi T, et al. Enzymatic changes in cerebrospinal fluid of patients with Alzheimer-type dementia. *J Clin Biochem Nutr* 1993;14:133–9.
12. Dostal DE, Hunt RA, Kule CE, Bhat GJ, Karoor V, McWhinney CD, et al. Molecular mechanisms of angiotensin II in modulating cardiac function: intracardiac effects and signal transduction pathways. *J Mol Cell Cardiol* 1997;29:2893–902.
13. Roks A, Buikema H, Pinto YM, van Gilst WH. The reninangiotensin system and vascular function: the role of angiotensin II, angiotensin- converting enzyme, and alternative conversion of angiotensin I. *Heart Vessels* 1997;Suppl 12:119–24.
14. Ferrario CM. Role of angiotensin II in cardiovascular disease—therapeutic implications of more than a century of research. *J Renin- Angiotensin-Aldosterone Syst* 2006;7:3–14.
15. Corrêa TD, Takala J, Jakob SM. Angiotensin II in septic shock. *Crit Care* 2015;19:98.
16. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45:942-6.
17. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15:108-15.
18. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
19. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
20. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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BIBLIOGRAFIA

1. Rehfeld L, Funk E, Jha S, Macheroux P, Melander O, Bergmann A. Novel Methods for the Quantification of Dipeptidyl Peptidase 3 (DPP3) Concentration and Activity in Human Blood Samples. J Appl Lab Med. 2018. doi:10.1373/jalm.2018.027995.
2. Wattiaux R, Wattiaux-de Coninck S, Thirion J, Gasingirwa MC, Jadot M. Lysosomes and Fas-mediated liver cell death. Biochem J. 2007 Apr 1;403(1):89-95. PubMed PMID: 17129211; PubMed Central PMCID: PMC1828900.
3. Pang X, Shimizu A, Kurita S, Zankov DP, Takeuchi K, Yasuda-Yamahara M, et al. Novel therapeutic role for dipeptidyl peptidase III in the treatment of hypertension. Hypertension 2016;68:630–41.
4. Kumar P, Reithofer V, Reisinger M, Wallner S, Pavkov-Keller T, Macheroux P, Gruber K. Substrate complexes of human dipeptidyl peptidase III reveal the mechanism of enzyme inhibition. Sci Rep. 2016 Mar 30;6:23787. doi: 10.1038/srep23787. PubMed PMID: 27025154; PubMed Central PMCID: PMC4824452.
5. Deniau B, Rehfeld L, Santos K, Dienelt A, Azibani F, Sadoune M, Kounde PR, Samuel LS, Tolpannen H, Lassus J, Harjola VP, Vodovar N, Bergmann A, Hartmann O, Mebazaa A, Blet A. Circulating dipeptidyl peptidase-3 is a myocardial depressant factor: DPP3 inhibition rapidly and sustainably improves hemodynamics, European Journal of Heart Failure 2019, accepted.
6. Takagi K, Blet A, Levy B, Deniau B, Azibani F, Feliot E, Bergmann A, Santos K, Hartmann O, Gayat E, Mebazaa A, Kimmoun A. Circulating dipeptidyl-peptidase 3 and alteration in hemodynamics in cardiogenic shock: Results from the OptimaCC Trial, European Journal of Heart Failure 2019, accepted.
7. Information on EC 3.4.14.4 - dipeptidyl-peptidase III, BRENDA. (n.d.). <https://www.brenda-enzymes.org/enzyme.php?ecno=3.4.14.4>
8. Prajapati SC, Chauhan SS. Dipeptidyl peptidase III: a multifaceted oligopeptide N-end cutter. FEBS J 2011;278: 3256–76.
9. Shimamori Y, Watanabe Y, Fujimoto Y. Purification and characterization of dipeptidyl aminopeptidase III from human placenta. Chem Pharm Bull 1986;34:3333–40.
10. Vanha-Perttula T. Dipeptidyl peptidase III and alanyl aminopeptidase in the human seminal plasma: origin and biochemical properties. Clin Chim Acta 1988;177:179–95.
11. Aoyagi T, Wada T, Kojima F, Nagai M, Harada S, Takeuchi T, et al. Enzymatic changes in

- cerebrospinal fluid of patients with Alzheimer-type dementia. *J Clin Biochem Nutr* 1993;14:133–9.
12. Dostal DE, Hunt RA, Kule CE, Bhat GJ, Karoor V, McWhinney CD, et al. Molecular mechanisms of angiotensin II in modulating cardiac function: intracardiac effects and signal transduction pathways. *J Mol Cell Cardiol* 1997;29:2893–902.
 13. Roks A, Buikema H, Pinto YM, van Gilst WH. The reninangiotensin system and vascular function: the role of angiotensin II, angiotensin- converting enzyme, and alternative conversion of angiotensin I. *Heart Vessels* 1997;Suppl 12:119–24.
 14. Ferrario CM. Role of angiotensin II in cardiovascular disease—therapeutic implications of more than a century of research. *J Renin- Angiotensin-Aldosterone Syst* 2006;7:3–14.
 15. Corrêa TD, Takala J, Jakob SM. Angiotensin II in septic shock. *Crit Care* 2015;19:98.
 16. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45:942-6.
 17. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15:108-15.
 18. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
 19. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
 20. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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REFERÊNCIAS

1. Rehfeld L, Funk E, Jha S, Macheroux P, Melander O, Bergmann A. Novel Methods for the Quantification of Dipeptidyl Peptidase 3 (DPP3) Concentration and Activity in Human Blood Samples. J Appl Lab Med. 2018. doi:10.1373/jalm.2018.027995.
2. Wattiaux R, Wattiaux-de Coninck S, Thirion J, Gasingirwa MC, Jadot M. Lysosomes and Fas-mediated liver cell death. Biochem J. 2007 Apr 1;403(1):89-95. PubMed PMID: 17129211; PubMed Central PMCID: PMC1828900.
3. Pang X, Shimizu A, Kurita S, Zankov DP, Takeuchi K, Yasuda-Yamahara M, et al. Novel therapeutic role for dipeptidyl peptidase III in the treatment of hypertension. Hypertension 2016;68:630–41.
4. Kumar P, Reithofer V, Reisinger M, Wallner S, Pavkov-Keller T, Macheroux P, Gruber K. Substrate complexes of human dipeptidyl peptidase III reveal the mechanism of enzyme inhibition. Sci Rep. 2016 Mar 30;6:23787. doi: 10.1038/srep23787. PubMed PMID: 27025154; PubMed Central PMCID: PMC4824452.
5. Deniau B, Rehfeld L, Santos K, Dienelt A, Azibani F, Sadoune M, Kounde PR, Samuel LS, Tolpannen H, Lassus J, Harjola VP, Vodovar N, Bergmann A, Hartmann O, Mebazaa A, Blet A. Circulating dipeptidyl peptidase-3 is a myocardial depressant factor: DPP3 inhibition rapidly and sustainably improves hemodynamics, European Journal of Heart Failure 2019, accepted.
6. Takagi K, Blet A, Levy B, Deniau B, Azibani F, Feliot E, Bergmann A, Santos K, Hartmann O, Gayat E, Mebazaa A, Kimmoun A. Circulating dipeptidyl-peptidase 3 and alteration in hemodynamics in cardiogenic shock: Results from the OptimaCC Trial, European Journal of Heart Failure 2019, accepted.
7. Information on EC 3.4.14.4 - dipeptidyl-peptidase III, BRENDA. (n.d.). <https://www.brenda-enzymes.org/enzyme.php?ecno=3.4.14.4>
8. Prajapati SC, Chauhan SS. Dipeptidyl peptidase III: a multifaceted oligopeptide N-end cutter. FEBS J 2011;278: 3256–76.
9. Shimamori Y, Watanabe Y, Fujimoto Y. Purification and characterization of dipeptidyl aminopeptidase III from human placenta. Chem Pharm Bull 1986;34:3333–40.
10. Vanha-Perttula T. Dipeptidyl peptidase III and alanyl aminopeptidase in the human seminal plasma: origin and biochemical properties. Clin Chim Acta 1988;177:179–95.

11. Aoyagi T, Wada T, Kojima F, Nagai M, Harada S, Takeuchi T, et al. Enzymatic changes in cerebrospinal fluid of patients with Alzheimer-type dementia. *J Clin Biochem Nutr* 1993;14:133–9.
12. Dostal DE, Hunt RA, Kule CE, Bhat GJ, Karoor V, McWhinney CD, et al. Molecular mechanisms of angiotensin II in modulating cardiac function: intracardiac effects and signal transduction pathways. *J Mol Cell Cardiol* 1997;29:2893–902.
13. Roks A, Buikema H, Pinto YM, van Gilst WH. The reninangiotensin system and vascular function: the role of angiotensin II, angiotensin- converting enzyme, and alternative conversion of angiotensin I. *Heart Vessels* 1997;Suppl 12:119–24.
14. Ferrario CM. Role of angiotensin II in cardiovascular disease—therapeutic implications of more than a century of research. *J Renin- Angiotensin-Aldosterone Syst* 2006;7:3–14.
15. Corrêa TD, Takala J, Jakob SM. Angiotensin II in septic shock. *Crit Care* 2015;19:98.
16. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45:942-6.
17. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15:108-15.
18. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
19. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
20. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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REFERENCIAS

1. Rehfeld L, Funk E, Jha S, Macheroux P, Melander O, Bergmann A. Novel Methods for the Quantification of Dipeptidyl Peptidase 3 (DPP3) Concentration and Activity in Human Blood Samples. *J Appl Lab Med*. 2018. doi:10.1373/jalm.2018.027995
2. Wattiaux R, Wattiaux-de Coninck S, Thirion J, Gasingirwa MC, Jadot M. Lysosomes and Fas-mediated liver cell death. *Biochem J*. 2007 Apr 1;403(1):89-95. PubMed PMID: 17129211; PubMed Central PMCID: PMC1828900.
3. Pang X, Shimizu A, Kurita S, Zankov DP, Takeuchi K, Yasuda-Yamahara M, et al. Novel therapeutic role for dipeptidyl peptidase III in the treatment of hypertension. *Hypertension* 2016;68:630–41.
4. Kumar P, Reithofer V, Reisinger M, Wallner S, Pavkov-Keller T, Macheroux P, Gruber K. Substrate complexes of human dipeptidyl peptidase III reveal the mechanism of enzyme inhibition. *Sci Rep*. 2016 Mar 30;6:23787. doi: 10.1038/srep23787. PubMed PMID: 27025154; PubMed Central PMCID: PMC4824452.
5. Deniau B, Rehfeld L, Santos K, Dienelt A, Azibani F, Sadoune M, Kounde PR, Samuel LS, Tolpannen H, Lassus J, Harjola VP, Vodovar N, Bergmann A, Hartmann O, Mebazaa A, Blet A. Circulating dipeptidyl peptidase-3 is a myocardial depressant factor: DPP3 inhibition rapidly and sustainably improves hemodynamics, *European Journal of Heart Failure* 2019, accepted.
6. Takagi K, Blet A, Levy B, Deniau B, Azibani F, Feliot E, Bergmann A, Santos K, Hartmann O, Gayat E, Mebazaa A, Kimmoun A. Circulating dipeptidyl-peptidase 3 and alteration in hemodynamics in cardiogenic shock: Results from the OptimaCC Trial, *European Journal of Heart Failure* 2019, accepted.
7. Information on EC 3.4.14.4 - dipeptidyl-peptidase III, BRENDA. (n.d.). <https://www.brenda-enzymes.org/enzyme.php?ecno=3.4.14.4>
8. Prajapati SC, Chauhan SS. Dipeptidyl peptidase III: a multifaceted oligopeptide N-end cutter. *FEBS J* 2011;278: 3256–76.
9. Shimamori Y, Watanabe Y, Fujimoto Y. Purification and characterization of dipeptidyl aminopeptidase III from human placenta. *Chem Pharm Bull* 1986;34:3333–40.
10. Vanha-Perttula T. Dipeptidyl peptidase III and alanyl aminopeptidase in the human seminal

plasma: origin and biochemical properties. Clin Chim Acta 1988;177:179–95.

11. Aoyagi T, Wada T, Kojima F, Nagai M, Harada S, Takeuchi T, et al. Enzymatic changes in cerebrospinal fluid of patients with Alzheimer-type dementia. J Clin Biochem Nutr 1993;14:133–9.
12. Dostal DE, Hunt RA, Kule CE, Bhat GJ, Karoor V, McWhinney CD, et al. Molecular mechanisms of angiotensin II in modulating cardiac function: intracardiac effects and signal transduction pathways. J Mol Cell Cardiol 1997;29:2893–902.
13. Roks A, Buikema H, Pinto YM, van Gilst WH. The reninangiotensin system and vascular function: the role of angiotensin II, angiotensin- converting enzyme, and alternative conversion of angiotensin I. Heart Vessels 1997;Suppl 12:119–24.
14. Ferrario CM. Role of angiotensin II in cardiovascular disease—therapeutic implications of more than a century of research. J Renin- Angiotensin-Aldosterone Syst 2006;7:3–14.
15. Corrêa TD, Takala J, Jakob SM. Angiotensin II in septic shock. Crit Care 2015;19:98.
16. Kricka LJ. Human anti-animal antibody interference in immunological assays. Clin Chem 1999;45:942–6.
17. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. J Clin Immunoassay 1992;15:108–15.
18. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
19. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
20. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.

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5. Deniau B, Rehfeld L, Santos K, Dienelt A, Azibani F, Sadoune M, Kounde PR, Samuel LS, Tolpannen H, Lassus J, Harjola VP, Vodovar N, Bergmann A, Hartmann O, Mebazaa A, Blet A. Circulating dipeptidyl peptidase-3 is a myocardial depressant factor: DPP3 inhibition rapidly and sustainably improves hemodynamics, *European Journal of Heart Failure* 2019, accepted.
6. Takagi K, Blet A, Levy B, Deniau B, Azibani F, Feliot E, Bergmann A, Santos K, Hartmann O, Gayat E, Mebazaa A, Kimmoun A. Circulating dipeptidyl-peptidase 3 and alteration in hemodynamics in cardiogenic shock: Results from the OptimaCC Trial, *European Journal of Heart Failure* 2019, accepted.
7. Information on EC 3.4.14.4 - dipeptidyl-peptidase III, BRENDA. (n.d.). <https://www.brenda-enzymes.org/enzyme.php?ecno=3.4.14.4>
8. Prajapati SC, Chauhan SS. Dipeptidyl peptidase III: a multifaceted oligopeptide N-end cutter. *FEBS J* 2011;278: 3256-76.
9. Shimamori Y, Watanabe Y, Fujimoto Y. Purification and characterization of dipeptidyl aminopeptidase III from human placenta. *Chem Pharm Bull* 1986;34:3333-40.
10. Vanha-Perttula T. Dipeptidyl peptidase III and alanyl aminopeptidase in the human seminal plasma: origin and biochemical properties. *Clin Chim Acta* 1988;177:179-95.
11. Aoyagi T, Wada T, Kojima F, Nagai M, Harada S, Takeuchi T, et al. Enzymatic changes in cerebrospinal fluid of patients with Alzheimer-type dementia. *J Clin Biochem Nutr* 1993;14:133-9.
12. Dostal DE, Hunt RA, Kule CE, Bhat GJ, Karoor V, McWhinney CD, et al. Molecular mechanisms of angiotensin II in modulating cardiac function: intracardiac effects and signal transduction pathways. *J Mol Cell Cardiol* 1997;29:2893-902.
13. Roks A, Buikema H, Pinto YM, van Gilst WH. The reninangiotensin system and vascular function: the role of angiotensin II, angiotensin- converting enzyme, and alternative conversion of angiotensin I. *Heart Vessels* 1997;Suppl 12:119-24.
14. Ferrario CM. Role of angiotensin II in cardiovascular disease—therapeutic implications of more than a century of research. *J Renin- Angiotensin-Aldosterone Syst* 2006;7:3-14.
15. Corrêa TD, Takala J, Jakob SM. Angiotensin II in septic shock. *Crit Care* 2015;19:98.
16. Kricka LJ. Human anti-animal antibody interference in immunological assays. *Clin Chem* 1999;45:942-6.
17. Levinson SS. The nature of heterophilic antibodies and their role in immunoassay Interference. *J Clin Immunoassay* 1992;15:108-15.
18. Pierson-Perry J, Vaks JE, Durham AP, et al. Evaluation of detection capability for clinical laboratory measurement procedures; Approved Guideline – Second Edition. CLSI Document EP17-A2, Volume 32 Number 8, 2012.
19. McEnroe RJ, Burritt MF, Powers DM, et al. Interference testing in clinical chemistry; Approved Guideline – Second Edition. CLSI Document EP7-A2, Volume 25 Number 27, 2005.
20. McEnroe RJ, Magari R, Durham PA, et al. Evaluation of precision of quantitative measurement procedures; Approved Guideline – Third Edition. CLSI Document EP05-A3, Volume 34 Number 13, 2014.