

	DOCUMENT#	FRM-IT-008	REV	B	STATUS	RELEASED
	CHANGE REQUEST#	2021-358	EFFECTIVE DATE		2021-11-02	
	FORM					Page 1 of 1
TITLE	IB10 SOFTWARE RELEASE NOTES					

IB10 Software Release Notes		
Date	Version	Description
11/01/2018	2.00.00	Nexus DX Initial Release. The following panels added: -R70 (IB10 sphingotest® Troponin-99) -R71 (IB10 sphingotest® 3-in-1 Cardiac) -R72 (IB10 sphingotest® D-Dimer) -R73 (IB10 sphingotest® SOB) -R74 (IB10 sphingotest® beta-hCG) -R75 (IB10 sphingotest® TSH) -R76 (IB10 sphingotest® NT-proBNP)
04/08/2019	2.01.00	The following panels were added: -R77 (IB10 sphingotest® DPP3)
12/03/2019	2.02.00	The following panels were added: -R79 (IB10 sphingotest® penKid®)
01/21/2020	2.03.00	The following panels were added: -R80 (IB10 sphingotest® PCT)
04/17/2020	2.10.00	The following panels were added: -R78 (IB10 sphingotest® bio-ADM®)
08/18/2021	2.21.00	The following panels were modified: -R77 (IB10 sphingotest® DPP3) Assay time increased to 22 minutes

Notice to reader: The information above is believed to be accurate and represents the best information available to us. However neither the above named supplier nor any of its subsidiaries assume any liability for the accuracy or completeness of the information. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no event shall Nexus Dx be liable for any claims, losses or damages of any third party for loss of profits or any special, indirect, incidental, consequential or exemplary damages.